REALPURE DNA FFPE

Ref. **RBMEGS12 50 Test**

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1. INTRODUCTION

*This kit is optimized for a fast method to isolate **DNA from formalin-fixed, paraffin-embedded (FFPE) tissue specimen.***

*The procedure omits the use of flammable and malodorous xylene or d-limonene commonly used for deparaffinization, **a proprietary buffer formulation REAL DEPARAFFINIZATION SOLUTION** is used for the complete dissolution of the wax to release the tissue.*

Features:

- *Silica membrane technology with special MicroSpin columns.*
- *Very easy paraffin removal.*
- *Safe method avoids xylene and other toxic.*
- *Complete removal of contaminants and inhibitors for reliable downstream applications.*
- *Low elution volume: 15-30 µl.*
- *The quality of DNA is suitable for the following applications as quantitative PCR or Next generation sequencing (NGS).*

Applications:

- *Rapid isolation of DNA from formalin-fixed, paraffin-embedded samples.*
- *Isolation of DNA from fresh and archived FFPE samples*
- *Isolation of DNA from specimen of object slides Typical downstream application: PCR, pPCR, NGS, STR analysis.*

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(*) *These solutions must be prepared as indicated in the 3.1 section “Preliminary Preparations”*

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- * Etanol 100 %.
- * Microtubos de 1.5 ml
- * Microcentrifuga.

Equipment and additional reagents required

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3.1 Preliminary Preparations

- **ATTENTION** Both the Lysis/ Binding Buffer (REA011) and the Desinhibition Buffer (REA03) contain Guanidine hydrochloride which is an irritant agent, for this reason we recommend to use gloves and glasses for its manipulation.
- Dissolve the proteinase K in 2.75 ml of nuclease-free water and store at -20°C. It is recommended to do several aliquots to avoid many thaw/freeze cycles. At this temperature it is stable for 1 year.
- Add 10 ml of Ethanol 100 % to the Desinhibition Buffer (REA03). Keep the container closed to avoid the ethanol evaporation
- Add 40 ml of Ethanol 100 % to the Wash Buffer (REA04). Keep the container closed to avoid the ethanol evaporation.
- Pre-heat the Elution Buffer at 70°C.

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3.2 Protocol for genomic DNA extraction from FFPE samples

DEPARAFFINIZATION SAMPLE

1. Add **400 µl** of **REAL Deparaffinization solution** to one FFPE tissue slice (10µm). Make sure the sample does not comprise more than 15 mg paraffin. **Vortex for 10 seconds.**
2. Incubate **3 minutes at 60°C** to promote the melting of the paraffin.
3. Vortex the sample immediately (at 60°C) at a vigorous speed to dissolve the paraffin. After 3 minutes you will see the tissue floating in the wax remover.
4. **Centrifuge full speed for 3 minutes** to pellet the tissue.
5. **Remove the supernatant** with pipetting, avoiding the tissue.
6. **Add 1 ml ethanol 100%.** Vortex for 20 seconds.
7. **Centrifuge full speed for 3 minutes** to pellet the tissue. Remove the ethanol with pipetting, avoiding the tissue.
8. Place the tubes **at 55°C for 10 minutes** with caps open to evaporate the ethanol.

DNA ISOLATION

1. Add **180 µl** of the **Tissue Lysis Buffer + 50 µl Proteinase K.** Mix by vortexing 2-5 seconds.
2. **Incubate at 55°C** for 1 hour or until the lysis is complete.
3. **Incubate at 90°C for 1 hour.** This incubation partially reverses formaldehyde modification of nucleic acids. If using only one heating block, leave the sample at room temperature after the 55°C incubation until the heating block has reached 90°C.
4. **Briefly centrifuge** to remove drops from the inside of the lid.
5. Add **200 µl** of **Lysis/Binding Buffer.** Mix by vortexing. **Incubate at 70°C for 5 minutes.**
6. Add **200 µl** of **Ethanol (96-100%)** to the lysate. Mix by vortexing. Briefly centrifuge the 1.5 ml tube to remove drops from inside the lid.
7. **Transfer the lysate** into reservoir of a combined MicroSpin Column –collection tube.
8. **Centrifuge at 8.000 rpm for 60 seconds.** Remove the collection tube. If the sample is not drawn completely through the matrix, repeat the centrifugation step.
9. Place the MicroSpin column in a new collection tube and add **500 µl** of **Deshinhibition Buffer.**
10. **Centrifuge at 12.000-14.000 rpm for 60 seconds.** Remove the liquid.
11. Add **500 µl** of **Wash Buffer** into reservoir of MicroSpin column.
12. **Centrifuge at 12.000-14.000 rpm for 60 seconds.** Remove the liquid.
13. **Centrifuge at maximum speed for 3 minutes** to remove the residual ethanol.
14. Remove the collection tube and insert the MicroSpin column in a 1.5 ml microtube. Add **15- 30 µl** of **Elution Buffer** (preheated at 70°C) directly onto the center of the silica membrane.
15. **Incubate 1 minute.**
16. **Centrifuge at maximum speed for 60 seconds.** The microtube contains now genomic DNA.