REAL TOTAL RNA SPIN BLOOD KIT

Ref. RBMER12 50 Preps.

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1. INTRODUCTION

*This kit allows the obtaining of **DNA-free total RNA** directly from **whole blood** without previous lymphocytes isolation and using for such purpose a **Stabilizing Lysis Solution and silica membrane columns**.*

When blood is collected for a later RNA extraction, it is usually necessary to storage previous to the extraction, as such extraction must be done in a different place. If we add to this the fact that the storage of RNA from blood is difficult due to the extremely high RNA sensibility to degradation by ribonucleases, which erythrocytes are rich in, it becomes necessary to have a system to stabilize the RNA in total blood.

*For that reason we use a **Stabilizing Lysis Solution** that simultaneously lysis the cells, inactivates the RNases and stabilizes the RNA preserving its integrity. **This way, the samples can be transported and stored for 2 days at 4°C or for 15 days at -20°C.***

You can process 50 samples of 300 µl of whole blood collected according to the standard methods in tubes containing anticoagulants, the use of EDTA is recommended as anticoagulant. You can also use ACD, we don't recommend the use of heparin as it produces a low efficiency in RT-PCR.

The total RNA yield from 300 µl of total blood is usually around 0.2–0.3µg, enough quantity to do about 10 RT-PCR experiments.

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Production system certified under ISO

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3.1 Preliminary Preparations

- **Reconstitute the lyophilized RNase-free DNaseI, add 540 µl of nuclease-free water to the vial** and incubate 1 minute at room temperature, mix until the DNase I is dissolved. Dispense into aliquots and store at -20°C. Do not thaw/freeze the aliquots more than 3 times
- **Add 50 ml of 100 % Alcohol to the concentrated wash solution, mark the date.**
- **ATTENTION:** The Stabilizing Lysis Solution, the DNaseI STOP solution and the Membrane Desalting Buffer contains guanidine thiocyanate. Use gloves and glasses for protection. In case of contact wash with water.

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3.3 Protocol for RNA extraction from 300 µl of whole blood

1. Collect **300 µl of whole blood** according to the Standard methods in tubes containing anticoagulants, and mix the sample carefully. **Immediately**, put the sample into a microtube containing **500 µl Stabilizing Lysis Solution**. This way the samples can be transported and stored for 2 days at 4°C or for 15 days at – 20°C. Process the sample if desired.
2. Reduce viscosity and clean the lysate by filtration through a **filtration column**. Place the filtration column in a collection tube, apply the lysate and **centrifuge 1 minute at 11.000 x g**. Transfer the filtrate into a new centrifuge tube. **Don't shake the pellet** from cell debris, which can be observed at the bottom of the collection tube.. If this happens centrifuge once again this tube and collect the filtrated.
3. Add **700 µl of 70% Ethanol** to the filtrated obtained in step 2 and mix with vortex or pipette. After adding the ethanol you can see a precipitate that won't affect the process. Be sure to load all the precipitate on the column.
4. Take a **RNA Binding Column** plus collection tube and add **750 µl of the lysate**. Centrifuge at **8.000 x g (10.000 r.p.m.) for 30 seconds**. Repeat the process exactly the same with the **750 µl of the lysate remaining**. Place the column in a new collection tube.
5. Add **350 µl of the Membrane Desalting Buffer** and centrifuge at **maximum speed (11.000 x g) for 1 minute** to dry the membrane. The salt removal will produce a **DNase I digestion** much more effective. If the column outlet has come into contact with the flow-through for any reason, discard the flow-through and centrifuge again for 30 seconds at 11.000 x g.
6. Prepare the **DNase I Reaction mixture** in a sterile microtube, for each extraction, add **10 µl of reconstituted DNase I** to **90 µl of DNase I Reaction Buffer**. Mix.
7. Apply **95 µl of DNase Reaction mixture** directly onto the centre of the silica membrane of the column. Incubate at **room temperature for 15 minutes**.
8. Add **200 µl of STOP/DNase I Solution** to the RNA binding column. Centrifuge at **8.000 x g for 30 seconds**. Place the tube in a new collection tube. This buffer will inactivate the DNase I.
9. Add **600 µl of Wash Solution** to the column. Centrifuge at **8.000 x g for 30 seconds**. Remove the liquid and place the column back in the collection tube.
10. Add **250 µl of Washing Solution** to the column. Centrifuge for **2 minutes at maximum speed (11.000 x g)** to completely dry the membrane. Place the column in a 1.5 ml (supplied) microtube. If for any reason, the liquid touches the column after the centrifugation, repeat the process.
11. Elute the RNA in **30 µl of Nuclease-free Water** (supplied) be sure that all the membrane is well soaked. **Incubate for 1 minute at room temperature**. Centrifuge at **maximum speed (11.000 x g) for 1 minute**.
12. Add again **30 µl of Nuclease-free Water** to the same column, with the same collection tube. Centrifuge at **maximum speed (11.000 x g) for 1 minute** Finally it will be obtained **60 µl of RNA**..

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- 2.3. *The ethanol has not been added after the lysis. The RNA will only bind to the silica membrane under the presence of alcohol*
 - 2.4. *Store the kit components as indicated in the protocol. Low temperature storages can cause salts precipitation*
 - 2.5. *Keep the containers correctly closed to prevent evaporation*
 - 2.6. *The sample hasn't been correctly preserved. Use, if possible, fresh material. If it is not possible, freeze the sample in liquid Nitrogen and store them at -70°C until you add the Lysis Solution*
 - 2.7. *The initial material hasn't been enough homogenized. Always spray the samples in presence of liquid Nitrogen. Use the filtration columns.*
- 3. Overloaded membrane**
- 3.1. *Too much initial material has been used. An overloaded membrane can lead to a low yield. Reduce the sample quantity or use a bigger quantity of Lysis Solution.*
- 4. Contamination with Genomic DNA**
- 4.1. *Not active DNase I. Reconstitute and store the DNase I according to the supplied instructions.*
 - 4.2. *Pipette the DNase I solution directly on the centre of the silica membrane.*
 - 4.3. *Reduce the sample starting quantity.*
- 5. Low RNA yield in following applications or experiments**
- 5.1. *The Wash solution that contains ethanol has not been completely removed. Be sure that the end of the column does not touch the liquid after 2nd wash.*
 - 5.2. *Check that the Wash solution has been balanced at room temperature before use. Washing at low temperatures produces a decrease on the salt removal efficiency.*