REAL SPIN DTR

Ref. RBMS01 50 Preps.

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1. INTRODUCTION

Real Spin DTR are ready to use chromatographic microcolumns to obtain a fast purification and clean up of nucleic acids and proteins with the use of a microcentrifuge.

For the purification it uses the filtration in gel method, in which the molecules in solution are separate according to their different sizes when they pass through a microcolumn that has been packed with a chromatographic medium, in a gel form. The gel is an heterogeneous phase in which a watery liquid stage is inside the pores of a solid stage. These pores have a controlled range of size according to the matrix gel used.

Small molecules diffuse through the gel and are delayed on their way down the column. On the opposite, big molecules can not diffuse through the gel and leave quickly the gel followed by smaller molecules in order of size.

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Empresa Registrada

ER-0989/2/98

Production system certified under ISO

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* Microcentrifuga.

Equipment and additional reagents required

* *Microcentrifuge.*

3. PROTOCOLO GENERAL GENERAL PROCEDURE

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3.1. Specifications

- Applications:** Removal of “dideoxy terminators” and “dye-labelled” nucleotides not incorporated in sequencing reactions before the analysis in automatic sequencers
- Gel Matrix:** Special polyacrylamide grade
- Removed Dyes:** 98 % of retention
- Big dye terminators, Rodamine and d-rodamine, Well RED dye.**
- Capacity:** 10-75 µl
- Time:** 4 minutes
- Recovered DNA:** 95 % DNA > 22 pb
- Buffer:** Distilled water

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$$\text{rpm} = 1000 \times (750 / 1.12 r)^{1/2}$$

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3.2 Working Protocol

1. Dissolve the DTR Real Spin resin by turning the tube upside down several times or shake carefully with vortex.
2. Remove the lower and upper caps and place in a 2.0 ml washing tube.
3. Centrifuge for 2 minutes in a microcentrifuge at 750 xg to remove the equilibration buffer.
4. Place the DTR Real Spin in a 1.5-ml sample collecting tube.
5. Carefully, apply the sample (10-75 µl) directly on the centre of the column without damaging the gel.
6. Centrifuge for 2 minutes at 750 xg. The purified sample is collected in the 1.5 ml collecting microtube. **Molecules smaller than**

the exclusion size will be retained in the DTR Real Spin.

*** All the centrifugations are done at 750 xg in a conventional microcentrifuge. The correct speed for each microcentrifuge can be calculated with the following formula:**

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Tal y como se recomienda en el protocolo, la velocidad adecuada son 750 xg; con estas condiciones no debe encontrarse problemas en este sentido
1. **Poor results**
 - 1.1. **High centrifugation speed.** It is very important to set the correct centrifugation speed. If a higher speed is applied, the matrix gel will not be formed in the correct way.
 - 1.2. **The sample has been placed incorrectly.** You must avoid touching the gel with the micropipette end.
 - 1.3. **The DTR Real Spin has been placed incorrectly.** The gel matrix is formed in the first centrifugation, with its upper part leaned. The DTR Real Spin will be placed on the rotor with this leaned part touching the outside of the rotor.
2. **Sample dilution**

Due to the application of a lower speed than the one indicated, the equilibration buffer will be not totally removed and it will be collected in the second centrifugation step, producing the sample dilution.

To avoid mistakes when setting de work speed it is recommended to do a short centrifugation, for 30-60 seconds after step 3 of the protocol, the fist time you work with the DTR Real Spin DTR. If there is a higher buffer quantity, it will be necessary to increase the speed, or to increase the centrifugation time in step 3. As it is recommended in the protocol, the correct speed is 750 xg; you will not have problems if you follow these conditions.