REAL DNA REMOVE

Ref. **RD055 (2 x 100 ml)**

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1. INTRODUCTION:

Nucleic acid amplification techniques have been a great advance inside the Molecular Biology field: by using such simple methods as the Polymerase Chain Reaction (PCR), a single DNA or RNA can be easily detected starting from a wide variety of samples. However, in some cases and due to the high sensibility of this technique, exogen DNAs, plasmids and most of all previous amplifications products, can be target DNA in later amplifications, giving false positives. To avoid contaminations there are control measures based on the use of enzymes such as Uracile DNA Glicosilase, which jeans modifications on the standard PCR method and a more expensive technique. There have been suggested other methods to minimize contamination, such as using filter tips, work in laminar flow cabinets, divide the lab in pre-amplification and post-amplification side, autoclaved work material, etc. Some of these measures are not effective and others are not suitable in many laboratories.

DNA Remove is a mixture of detergents specially designs to remove in a very simple and effective way the DNA or RNA from work surfaces and lab material.

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2. HOW TO USE:

1. *Spray the work surfaces and lab material with **DNA Remove** (you can also use a cloth wet in **DNA Remove**).*
2. *Wait for at least 3 minutes and wash with distilled water or ethanol.*
3. *To decontaminate the pipettes, submerge the autoclavable part in **DNA Remove** and wash with distilled water.*

3. CAUTION:

- *DNA Remove is compatible with materials such as glass, PVC, PE and INOX. Avoid long time contact with surfaces like iron, nickel or copper.*
- *Use glasses and gloves to prevent splashes.*
- *In case of contact in eyes or skin wash with water.*
- *Store at room temperature. Avoid exposure to high temperature and intense light*