REALCLEAN & CONCENTRATOR MICROSPIN KIT

Ref. **RBMCS04 50 Test**
Ref. **RBMCS05 250 Test**

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1. INTRODUCTION

*This kit provides a rapid method for purification and concentration of **high-quality DNA from PCR or enzymatic reactions with an extremely small elution volume of only 10 µl** using special MicroSpin columns.*

Applications:

- *PCR products clean-up: efficient desalting of DNA with the removal of DNA polymerases, primers and free dNTPs.*
- *DNA clean-up from Enzymatic Reactions, including: Desphosphorylation, Restriction enzyme digestion, Ligation, Primed synthesis, Endlabeling and Nick translation.*
- *Isotope and Dye Removal, efficiently removes unincorporated fluorescent (i.e., AMCA, FITC, BIO, DIG, Cy3, Cy5, FAM, etc) and radiolabeled dNTP derivats from DNA following in vitro labelling reactions.*

Features:

- *The microspin columns are designed to allow elution in very small volumes (as little as 10 µl) delivering highly concentrated DNA in high yields.*
- *DNA Size Limits: From 70 pb to 23 Kb.*
- *DNA Recovery: Typically, up to 5 µg total DNA per column can be eluted into as little as 10 µl.*
- *The protocol is done in 2 minutes.*
- *Fast procedure and easy handling.*
- *Eluted DNA is well suited for use in DNA ligation, sequencing, labelling, PCR, etc.*

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Equipment and additional reagents required

- * 100 % Ethanol
- * Isopropanol
- * Water Bath
- * 1.5 ml microtubes
- * Microcentrifuge.

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3.1 Preliminary Preparations

- **Add 40 ml** (ref. RBMCS04, 50 extractions) **or 200 ml** (ref. RBMCS05, 250 extractions) of 100% ethanol to the wash solution (EP08). Label the container and keep it closed to avoid ethanol evaporation.
- **Add 10 ml** (ref. RBMCS04, 50 extractions) **or 50 ml** (ref. RBMCS05, 250 extractions) of 100% Isopropanol to the binding solution 2 (CM02). Label the container and keep it closed to avoid isopropanol evaporation.
- Pre-heat the elution buffer at 70°C.

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8. **Incubar durante 2 minutos y centrifugar durante 1 minuto a 14.000 rpm.**

3.2 Working Protocol

1. **Add 5 volumes Binding Solution 2** (ref. CM02) **with isopropanol to 1 volume of PCR** (50 -100 µl) or enzymatic reaction. Mix well.
2. **Transfer the sample to a MicroSpin column.** Put the spin column in collecting tube.
3. **Centrifuge for 1 minute at 10.000-12.000 rpm.**
4. **Remove the filtrate and add 600 µl of the Washing solution** (ref. EP08). Centrifuge for 1 minute at 14.000 rpm.
5. **Remove the filtrate and add 200 µl of washing solution.** Centrifuge for 1 minute at 14.000 rpm.
6. **Remove the residual ethanol** by centrifugation for 3 minutes at 14.000 rpm.
7. Transfer the spin column into a new receive microtube and add **12 µl of pre-warmed**

Elution Buffer (10 mM Tris.Cl, pH 8.5) or water at 70°C.

Ensure that the elution buffer is dispensed directly onto the centre of the membrane for complete elution of bounded DNA. The average elute volume is 10 µl from 12 µl elution buffer volume.

Elution efficiency is dependent on pH. When using, water, make sure that the pH value is within 7.0 and 8.5 range, and store DNA at -20°C as DNA may degrade in the absence of a buffering agent.

8. **Incubate for 2 minutes and centrifuge for 1 minute at 14.000 rpm.**

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1. We recommend contact with **REAL** laboratory technical service at Durviz s.l., for any additional question regarding the work protocols or any problems you may have during work