# REAL GELES AGAROSA 1% 7x10 CM 10P REALSAFE C/6 U.

Ref. R0504RS:

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## 1.INTRODUCTION

This kit has all the necessary components for the visualisation of DNA in processes of amplification, cloning, restriction enzyme digestion, extraction, purification, etc., by means of electrophoretic separation in agarose gel and staining with RealSafe.

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## 3. GENERAL PROCEDURE

- Place 3 or 4 µl of loading buffer on wax paper (parafilm) and mix with 20 µl of amplification product, or add the loading buffer directly to the amplified product tube in a ratio of approximately 1:6 (loading buffer:sample). Load the mixture into the wells of the gel. The gel must be completely submerged in the 0.5x electrophoresis buffer contained in the electrophoresis cuvette (this buffer is originally at 5x concentration). It is preferable to use the central wells if you have few samples. Also load a positive visualisation control, which will indicate the size of the amplified sample.
- Subject to electrophoresis for the time necessary for the visible blue front to have migrated a minimum of 3 or 4 cm from the well.
- View the agarose gel in the U.V. transilluminator at 312 nm. Samples yielding an amplified fragment of the same size as the positive amplification control are considered positive. The size of the band can be confirmed by comparison with a molecular weight marker.

**NOTE:** If the staining with the dye incorporated in the gel has lost intensity, it can be re-stained by placing the gel, after electrophoresis, in a cuvette containing REALSAFE solution.

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## 4. STORAGE AND STABILITY

- Agarose gels should be stored at 4°C until use. Keep the box containing each gel closed to avoid dehydration.

B65699985