REAL GEL/PCR KIT


Ref. **RBMCS01:** (para **50 EXTRACCIONES**)
 Ref. **RBMCS02:** (para **250 EXTRACCIONES**)
 Ref. **RBMCS08:** (para **1000 EXTRACCIONES**)

(for **50 EXTRATION**)
 (for **250 EXTRATION**)
 (for **1000 EXTRATION**)

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
1.INTRODUCTION

The **REAL GEL/PCR SPIN Kit** was designed to recover or concentrate **DNA fragments from agarose gel, PCR or other enzymatic reactions.**

This kit includes a pH indicator which is premixed with the Union Solution to ensure optimal pH, facilitate DNA binding and allow for easy observation of undissolved agarose gel. If pH exceeds the optimal level (>7.5), the color of the solution will appear purple instead of yellow. **3M Sodium Acetate** (pH5.0), which is included with the kit, can then be added to the solution to adjust pH and return the color to yellow.

Features:

- **High Efficiency:** up to 90% recovery from agarose gels, up to 95% recovery from PCR or other enzymatic reactions.
- **Sample Size:** up to 300 mg of agarose gel slice or up to 100 µl of PCR products.
- **Broad Fragment Size Range: 70 bp–20 kb.**
- Primers and primer/dimer removal.
- The protocol is done in 10 minutes

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Equipment and reagents needed and not provided

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3. GENERAL PROCEDURE

3.1 Preliminary Preparations

- Add **40 ml** (50 test) or **200 ml** (250 test) of **100% Ethanol to the Washing Solution**. Label the container and keep it closed to avoid ethanol evaporation.
- Pre-heat the *Hydration Solution* at 70°C.

3.2 Protocol 1:

- ✓ for routine purifications of **70 bp - 20 Kb dsDNA PCR fragments** from 10-40 mer primers, dNTPs, enzymes and salts mixtures.

1. **Add 5 volumes of Solubilization and Binding Solution to 1 volume of PCR (50 -100 ml).** Mix well. If the mixture has turned from yellow to purple, add **10 µl 3M Sodium Acetate** (pH 5.0) and mix thoroughly. This will adjust pH and the color return to yellow.
2. **Transfer the sample to a MicroSpin Columns.** Put the *MicroSpin Columns* in *Collection Tubes*.
3. **Centrifuge for 1 minute at 10.000-12.000 rpm.**
4. **Remove the filtrate and use 600 µl of the Washing Solution.** Centrifuge for 1 minute at 14.000 rpm.
5. **Remove the filtrate and apply 200 µl of Washing Solution.** Centrifuge for 1 minute at 14.000 rpm.
6. **Remove the residual ethanol** by centrifugation for 3 minutes at 14.000 rpm.
7. Transfer the *MicroSpin Columns* into a new *Collection Tubes* and add **at least 25 µl of pre-warmed Hydration Solution at 70°C.**
8. **Incubate for 2 minutes and centrifuge for 1 minute at 14.000 rpm.**

3.3 Protocol 2:

- ✓ For solubilization of the agarose gel slice for extractions of DNA fragments

1. To purify agarose gels, cut the band and **try to minimize the fragment size** removing the excess of agarose. **Transfer up to 300 mg of gel slice** to a 1.5 ml microcentrifuge tube.

Note: If using less than 300 mg of gel slice, the buffer does not to be scaled. If using more than 300 mg of gel slice, separate it into multiple 1.5 ml microcentrifuge tube.

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- Add 500 µl of Solubilization and Binding Solution.**
Note: If the mixture has turned from yellow to purple, add **10 µl 3M Sodium Acetate** (pH 5.0) and mix thoroughly. This will adjust pH and the color return to yellow.
- Incubate for 10 minutes at 55°C** to dissolve the agarose.
Shake with a vortex every 2-3 minutes.
- Transfer the solution to MicroSpin Columns.** Put the *MicroSpin Columns* in *Collection Tubes*.
- Centrifuge for 1 minute at 10.000-12.000 rpm.**
- Remove the filtrate and use 600 µl of the Washing solution** (ref. EP08). Centrifuge for 1 minute at 14.000 rpm.
- Remove the filtrate and apply 400 µl of Washing solution.** Centrifuge for 1 minute at 14.000 rpm.
- Remove the residual ethanol** by centrifugation for 3 minutes at 14.000 rpm.
- Transfer the *MicroSpin Columns* into a new *Collection Tubes* and add **at least 25 µl of pre-warmed Hydration Solution at 70°C.**
- Incubate for 2 minutes and centrifuge for 1 minute at 14.000 rpm.**

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- All components are stable for 12 months from date of purchase when stored as directed.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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