# REAL CLEAN PCR KIT

Ref. **RBMCS06:** (para **50 EXTRACCIONES**)  
 Ref. **RBMCS03:** (para **250 EXTRACCIONES**)  
 Ref. **RBMCS07:** (para **1000 EXTRACCIONES**)

(for **50 EXTRACTION**)  
 (for **250 EXTRACTION**)  
 (for **1000 EXTRACTION**)

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## 1.INTRODUCTION

The Real Clean PCR Kit is a versatile and effective method for fast and efficient purification and concentration of DNA from solutions, or from TAE or TBE agarose gels. **It is specially designed for rapid clean-up of PCR products.** This method allows purification of up to 15 µg of DNA.

The purification is based on the selective adsorption of nucleic acids on silica membranes in mini-spin columns in the presence of chaotropic salts.

The purified DNA is eluted in a small volume of buffer (5 mM Tris HCl pH 8.5) and can be used in enzymatic manipulations including sequencing, cloning, restriction analysis, ligation, tagging and in vitro transcription.

DanaGene Clean PCR Kit allows DNA recovery from 0.1- 10 Kb with an efficiency ranging from 75-90% depending on fragment size.

### Applications:

- Especially suitable for rapid "clean-up" of PCR products, achieving complete removal of salts, primers, nucleotides, mineral oil and enzymes.
- Purification and desalting of fragments from solutions, such as ligation, restriction or sequencing reactions.

## 2.COMPONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envases (50 Preps)	Envases (250 Preps)	Envases (1000 Preps)	T°
Disinhibition Solution (GHCL 8M)	CM02	15 ml	60 ml	240 ml	RT
Washing solution	EP08	10 ml	50 ml	200 ml	RT
Hydration Solution	REA05	2 ml	10 ml	40 ml	RT
MicroSpin columns	RSC05	50 unid.	250 unid.	1000 unid.	RT
Collection Tubes	R30	50 unid.	250 unid.	1000 unid.	RT

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### Equipment and reagents needed and not provided

- \* Isopropanol 100%.
- \* 100 % Ethanol.
- \* 1.5 ml. microtubes
- \* Microcentrifuge.

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### 3. GENERAL PROCEDURE

#### 3.1 Preliminary Preparations

- Add **40 ml** (50 test), **200 ml** (250 test) or **800 ml** (1000 test) of **100% Ethanol** to the **Washing Solution**. Label the container and keep it closed to avoid ethanol evaporation.
- Add **10 ml** (50 test), **40 ml** (250 test) or **160 ml** (1000 test) of **100% Isopropanol** to the **Disinhibition Solution (GHCL 8M)**.
- Pre-heat the **Hydration Solution** at 70°C.

#### 3.2 Protocol 1:

1. **Add 4 volumes of Disinhibition Solution (GHCL 8M) with Isopropanol to 1 volume** of PCR mix (50-100 µl). For example, 50 µl PCR mix + 200 µl **Disinhibition Solution (GHCL 8M)** with Isopropanol. Mix well by pipetting.
2. **Transfer the sample to Microspin columns.** Place the **Microspin columns** in a **Collection tube**.
3. **Centrifuge for 1 minute at 10,000-12,000 r.p.m.**
4. **Remove filtrate and add 600 µl of Washing Solution.** Centrifuge for 1 minute at 14,000 r.p.m.
5. **Remove the filtrate and add 200 µl of Washing Solution.** Centrifuge for 1 minute at 14,000 rpm.
6. **Remove the residual ethanol** by centrifugation for 3 minutes at 14,000 rpm.
7. Place the **Microspin columns** in a new **Collection tube** and **add at least 25 µl of Hydration Solution** pre-heated to 70°C. A second elution step will increase the yield.
8. **Incubate for 2 minutes and centrifuge for**

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### 4. STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

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