# REAL PLANT DNA KIT

Ref. **RBMEG04:** (para **50 EXTRACCIONES**)  
Ref. **RBMEG05:** (para **200 EXTRACCIONES**)

(for **50 EXTRACTION**)  
(for **200 EXTRACTION**)

## 1. DESCRIPCIÓN

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## 1. DESCRIPTION

### 1.1 Product description

This kit provides a method for efficient and rapid extraction of **genomic DNA from plant and fungal tissues and cells.**

It is known that plants contain quantities of different substances (polysaccharides, polyphenols, etc.) and that plants of the same genus or related genera can present enormous variability in their biochemical compositions, thus making it difficult to have a single DNA extraction method for all plants.

To solve this problem and to be able to cover as many plants as possible REAL uses a **PVP Solution** that is able to bind polysaccharides and polyphenols that are released by cell lysis and that have the ability to form complexes with nucleic acids and degrade them or cause them to precipitate with them.

The procedure includes a homogenization of the sample in an **Plants and Fungi Extraction Solution** and **PVP Solution**. Lysis is completed by incubation in **Lysis Solution** and **RNase** at 37°C for 30 minutes.

Cellular proteins and cellular debris are removed using a Protein Removal Buffer, leaving the genomic DNA in solution. Finally, the genomic DNA is isolated by isopropanol precipitation.

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### Equipment and additional reagents required

- \* Isopropanol
- \* 70 % Ethanol
- \* 1.5, or 2.0 ml microtubes.
- \* 15 or 30 ml centrifuge tubes that can withstand high speeds (for extractions greater than 100 mg of tissue).
- \* Microcentrifuge or clinic centrifuge.
- \* Vortex
- \* Water Bath
- \* Liquid nitrogen and mortar for tissue spraying.
- \* Electric homogenizer

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### 3.1 Preliminary preparations

- If the **Lysis Solution** contains a precipitate due to low temperatures. Incubate at 37°C and mix to dissolve the precipitate.
- Store the **RNase** at 4°C. If the period of use of the kit is going to be long, it is recommended to make aliquots and store at -20°C.

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### 3.2. Preliminary considerations

- For effective DNA collection from **plants**, it is recommended to pulverize the sample in a porcelain mortar with liquid nitrogen. For soft and non-fibrous tissues, such as young leaves, flowers, etc., they can be homogenized in an electric homogenizer, but always in the case of fresh samples, dissecting the sample into small pieces, adding the **Plants and Fungi Extraction Solution + PVP Solution** and homogenizing. For hard or fibrous tissues, such as stems, seeds, etc. the use of liquid Ni is recommended.
- **For fungi** collect the mycelium from the culture by filtration, wash 3 times with sterile deionized water or PBS to remove the culture medium. Pulverize the tissue in a porcelain mortar with liquid Ni, store at -80°C or process the sample.
- **NOTE:** It may be necessary to vary the amount of starting material depending on species, condition, tissue preparation or genome size.
- For processing 100 mg samples it is necessary to use 15 or 30 ml tubes that can withstand high centrifugation speeds.

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### 3.3 Protocol for extraction of phenolic DNA from 10-20 mg of plant tissue or 10-20 mg of fungal tissue.

#### Cell lysis

1. Weigh the sample **10-20 mg of plant tissue or 10-20 mg of mushroom tissue** and add 400 µl of **Plants and Fungi Extraction Solution + 200 µl of PVP Solution**. Homogenize with an electric homogenizer for 20seconds.
2. Add **60 µl Lysis Solution + 3 µl RNase**. Vortex vigorously and incubate the sample at 37°C for 30 minutes. For fungi to increase lysis efficiency add 3 µl Proteinase K (20 mg/ml) and incubate at 55°C for 1 hour or overnight. If possible, invert or vortex periodically during incubation.

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### Protein precipitation

1. Cool the sample to room temperature.
2. Add **360 µl of Protein Precipitation Solution**.
3. Vortex vigorously for 20-30 seconds.
4. Incubate at **-20°C for 10 minutes**.
5. Centrifuge at **14,000 rpm for 5 minutes**. The protein precipitate should be observed to form a pellet. If particles pass from the pellet, centrifuge again as in the previous section.

### Precipitation of DNA

1. Transfer the supernatant containing the DNA to a microtube containing **600 µl isopropanol**. Mix by inversion several times.
2. Centrifuge at **14,000 rpm for 3 minutes**.
3. Remove the supernatant. Add **600 µl of 70% ethanol** and invert several times to wash the pellet.
4. Centrifuge at **14,000 rpm for 2 minutes**. Carefully remove all ethanol. Be careful not to lose the DNA pellet. Normally this pellet has a green-brown color and not the characteristic white color.
5. Invert the tube and allow to dry on absorbent paper for 15 minutes.

### DNA hydration

1. Add **100 µl of Hydration Solution**. Depending on the species, the pellet size may be larger, in such cases use a larger amount of Hydration Solution so that the solution does not become excessively viscous.
2. Incubate at 65°C for 1 hour with periodic shaking to aid DNA dispersion.
3. **Centrifuge at 14,000 rpm for 1 minute as particles may be present in the rehydrated DNA. Transfer the supernatant containing the DNA to a new microtube.**
4. Store at 2-8°C. For long term storage store at -20°C or -80°C.

### Subsequent removal of polysaccharides

In those species with high amounts of polysaccharides, a subsequent purification can be carried out, which usually produces a reduction in the DNA obtained, but of higher quality. For this purpose, assess whether subsequent DNA manipulations are inhibited.

1. Adjust the hydrated DNA solution to 0.5 M potassium acetate and 30 % ethanol. Incubate at -20°C for 10 minutes.
2. Centrifuge at 14,000 rpm for 3 minutes.
3. Collect the supernatant with a pipette and add 1 volume of isopropanol.
4. Centrifuge at 14,000 rpm for 3 minutes.
5. To wash with ethanol 70 %.
6. Centrifuge at 14,000 rpm for 2 minutes.
7. Invert the tube and allow to dry on absorbent paper for 15 minutes.
8. Resuspend the pellet in Rehydration Buffer.

**NOTE:** To process larger samples, see the attached table with the reagents scaled according to the sample size and proceed as described in the protocol applying the corresponding amounts

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Given the wide variety of samples that can be treated to extract genomic DNA with this kit, it is difficult to generalize the possible problems and solutions. For this reason, we recommend contacting DURVIZ S.L. technical service for any additional questions regarding working protocols or problems that may arise during the work.

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All components are stable for 12 months from date of purchase when stored and used as directed.

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