REAL SALIVA DNA KIT

Ref. **RBMEG06:** (para **50 EXTRACCIONES**)
 Ref. **RBMEG07:** (para **160 EXTRACCIONES**)

(for **50 EXTRACTION**)
 (for **160 EXTRACTION**)

1. DESCRIPCIÓN

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1. DESCRIPTION

REAL SALIVA Kit provides a method for the extraction of high-quality genomic DNA from **direct saliva samples, buccal swab samples with foam brushes, saliva samples collected with the REALSALIVA Sample Collection Kit and saliva samples collected with the ORAGENE self-collection kits (DNAGenotek).**

It is a fast, safe and economical method. It uses a deproteinization step with a novel saline buffer avoiding the use of toxic organic solvents such as phenol or chloroform.

In the case of saliva, obtaining DNA in the same individual presents an intrasubject variability, i.e., there is no approximate amount of DNA / ml as it may vary depending on the time of sample collection, and it may happen that a large amount of DNA is not obtained and subsequent applications that require a large amount of DNA cannot be carried out.

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Equipment and additional reagents required

- * Isopropanol
- * 70 % Ethanol
- * 1.5 ml of microtubes.
- * 15- or 30-ml centrifuge tubes
- * Microcentrifuge or clinic centrifuge.
- * Vortex
- * Water Bath

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3.1.Preliminary preparations

- If the lysis solution contains a precipitate due to low temperatures. Incubate at 37°C and mix to dissolve the precipitate.
- Saliva samples should be collected before meals, rinse the mouth with water and wait about 30 minutes to collect the sample. Direct saliva samples should be processed, or kept at 4°C if they will be processed in less than 2 hours.

3.2. Extraction from 600 µl of direct saliva samples

Cell Lysis

1. Centrifuge **600 µl of saliva** for 90 seconds at 13,000-16,000 x g. Remove the supernatant by pipette without damaging the visible cell pellet.
2. Add **600 µl Lysis solution + 5 µl Proteinase K (20 mg/ml) + 5 µl RNase** and resuspend by pipette to lyse the cells.
3. **Incubate at 37°C for 30-60 minutes.** If possible, vortex periodically during the incubation period.

Protein precipitation.

1. Allow samples to cool for 3 min at -20°C.
2. Add **250 µl of Protein Precipitation Solution** to the cell lysate.
3. Vortex vigorously at maximum speed for 20-30 seconds.
4. Centrifuge at 13,000-16,000 x g for 5 minutes. A precipitate will form. If floating particles are observed, centrifuge again after incubating for 5 minutes on ice.

DNA precipitation.

1. Transfer the supernatant containing the DNA to a new microtube with 600µl of Isopropanol.
2. Mix by inversion about 50 times.
3. Centrifuge at 13,000-16,000 x g for 2 minutes. The DNA will be visible as a white pellet.
4. Remove the supernatant and dry the tube briefly on absorbent paper. Add 600 µl of 70% Ethanol to wash the DNA.
5. Centrifuge at 13,000-16,000 x g for 1 minute. Carefully remove the supernatant without touching the DNA pellet. Centrifuge again briefly to collect the last drops of residual ethanol with a micropipette.
6. Invert the microtube on absorbent paper and allow to dry for about 5 minutes.

DNA hydration.

1. Add 100-750 µl of **Hydration Solution** depending on the size of the DNA pellet and resuspend with micropipette. Large pellets may require incubation at 55°C for 1 hour with periodic shaking to fully resuspend or incubate overnight at room temperature with gentle agitation. For very small pellets 25-50 µl of **Hydration Solution** can be used.
2. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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3.3. Extraction from buccal swab samples with REAL SALIVA SWABS

Sample collection

1. It is recommended that the individual from whom the sample is to be taken refrain from drinking coffee and eating any food for at least 30 minutes prior to collection. If this is not possible, a gentle rinse of the mouth with water only is recommended.
2. Collect the buccal cell sample with the swab. Rub the swab on the inside of the cheek (buccal wall) and gums with firm pressure about 20 times on each side of the face and each side of the swab.
3. Use immediately for extraction. Then insert the brush into the receptacle provided for shipping. The sample can remain in this container tube for 1 week at 22-37°C before extraction. For long term storage keep the sample in the container at -20°C for up to 6 months.
4. The use of foam brushes such as those provided are suitable as they allow recovery of almost all of the lysis buffer after the incubation period.

Cell lysis

1. Add **600 µl** of **Lysis Solution** into a 1.5 ml microtube. Cut the brush head with a bit of handle and insert into the microtube. Vortex vigorously to release the cells from the brush.
2. Add **5 µl Proteinase K (20 mg/ml)** + **5 µl RNase** and incubate at 37°C for 1 hour. Vortex several times during the incubation period.
3. Remove the brush head from the lysis solution, rubbing it against the walls to collect the maximum amount of liquid.

Protein precipitation

1. Let the samples cool for 3 min at -20°C.
2. Add **250 µl** of **Protein Precipitation Solution**.
3. Vortex vigorously for 30 sec.
4. Centrifuge at 13,000-16,000 x g for 5 min. The protein precipitate will be observed to form a pellet.

DNA Precipitation

1. Transfer the supernatant containing the DNA to a 1.5 ml tube containing **600 µl** of **Isopropanol** Mix by inversion several times.
2. Centrifuge at 13,000-16,000 x g for 3 min.
3. Remove the supernatant. Add **600 µl** of **70% ethanol** and invert several times to wash the DNA pellet.
4. Centrifuge at 13,000-16,000 x g for 1 minute. Carefully remove the supernatant without touching the DNA pellet. Centrifuge again briefly to collect the last drops of residual ethanol with micropipette.
5. Invert the tube and allow to dry on absorbent paper for 5 minutes.

DNA hydration

1. Add **30-50 µl** of **DNA Hydration Solution** and resuspend with micropipette.
2. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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3.4 Extraction from saliva samples preserved in the REALSALIVA RNA Sample Collection Kit.

2 OPTIONS: *Depend on the amount of DNA needed for downstream applications or the preference of working with a microcentrifuge with 1.5 ml tubes (option A) or a clinical centrifuge and 15 ml tubes (option B).*

A) Process 1.2 ml (saliva + saliva preservative solution) which will represent processing approximately 600 l of direct saliva.

Cell Lysis

1. In the REALSALIVA Sample Collection Kit a white pellet containing the cells from which DNA is to be extracted will be observed. Shake well the tube containing the 2 ml of collected saliva. **It is important** that a homogeneous solution is observed.
2. Centrifuge **1.2 ml (saliva + saliva preservative solution)** for 90 seconds at 13,000-16,000 x g. Remove the supernatant by pipette without damaging the visible cell pellet.
3. Add **600 µl of Lysis Solution + 5 µl de Proteinase K (20 mg/ml) + 5 µl de RNase** and resuspend by pipette to lyse the cells.
4. **Incubate at 37°C for 30-60 minutes.** If possible vortex periodically during the incubation period.

Protein precipitation.

1. Allow samples to cool for 3 min at -20°C.
2. Add **250 µl of Protein Precipitation Solution** to the cell lysate.
3. Vortex vigorously at maximum speed for 20-30 sec.
4. Centrifuge at 13,000-16,000 x g for 5 minutes. A precipitate will form. If floating particles are observed, centrifuge again after incubating for 5 minutes on ice.

DNA precipitation.

1. Transfer the supernatant containing the DNA to a new microtube containing **600 µl of Isopropanol**.
2. Mix by inversion about 50 times.
3. Centrifuge at 13,000-16,000 x g for 2 minutes. The DNA will be visible as a white pellet.
4. Remove the supernatant and dry the tube briefly on absorbent paper. Add **600 µl of 70% Ethanol** to wash the DNA.
5. Centrifuge at 13,000-16,000 x g for 1 minute. Carefully remove the supernatant without touching the DNA pellet. Centrifuge again briefly to collect the last drops of residual ethanol.
6. Invert the microtube on absorbent paper and allow to dry for about 5-10 minutes.

DNA hydration.

1. Add **100-750 µl of Hydration Solution** depending on the size of the DNA pellet and resuspend with micropipette. Large pellets may require incubation at 55°C for 1 hour with periodic shaking to fully resuspend or incubate overnight at room temperature with gentle agitation. For very small pellets **25-50 µl of Hydration Solution** can be used.
2. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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B) Process the entire contents of the REAL SALIVA RNA Sample Colletion Kit (saliva + saliva preservative solution).

Cell Lysis

1. In the REAL SALIVA Sample Collection Kit a white pellet containing the cells from which the DNA is going to be extracted will be observed. Shake well the tube containing the 2 ml of saliva collected. It is important that a homogeneous solution is visible.
2. Transfer the entire contents **of the tube (4.50 ml)** from the REAL SALIVA Sample Collection Kit **into a 15 ml centrifuge tube**. Make sure that the entire contents including the white pellet have been transferred.
3. Centrifuge at 4000 rpm for 2 minutes. A white pellet should form. Remove the supernatant **without losing any cell pellet**.
4. Add **2 ml of Lysis Solution + 15 µl of Proteinase K (20 mg/ml) + 15 µl de RNase** and resuspend by pipette to lyse the cells. Incubate at 37°C for 60 minutes. **Very important that the entire cell pellet is resuspended**. If possible vortex periodically during the incubation period.

Protein precipitation.

1. Allow the samples to cool for 3 min at -20°C.
2. Add **750 µl of Protein Precipitation Solution** to the cell lysate.
3. Vortex vigorously at maximum speed for 20-30 sec.
4. Centrifuge at 4000 rpm for 5 minutes. A precipitate will form. If floating particles are observed, centrifuge again after incubating for 5 minutes on ice.

Precipitation of DNA.

1. Transfer the supernatant containing the DNA to **a fresh 15 ml centrifuge tube** containing **2 ml Isopropanol**.
2. Mix by inversion about 25 times.
3. Centrifuge at 4000 rpm for 3 minutes. The DNA will be visible as a white pellet on the bottom or wall of the tube.
4. Remove the supernatant and dry the tube briefly on absorbent paper. Add **2 ml of 70% ethanol** to wash the DNA.
5. Centrifuge at 4000 rpm for 1 minute. Carefully remove the supernatant without touching the DNA pellet. Centrifuge again briefly to collect the last few drops of residual ethanol.
6. Invert the tube on absorbent paper and allow to dry for 5-10 minutes.

Hydration of DNA

1. Add 750-1,000 µl of sterile water or **Hydration Solution** and resuspend by pipetting.
2. Incubate at 50°C for 1 hour with periodic shaking to aid DNA rehydration, or incubate "overnight" at room temperature with gentle agitation.
3. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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3.5 Extraction from saliva samples stored in the ORAGENE self collection kits (DNA Genotek). **OG-500 (2ml saliva) / OG-501 (1ml saliva)**

Cell Lysis

1. Incubate the ORAGENE/saliva samples at 50°C in a water bath for 1 hour or in an air incubator for a minimum of 2 hours.
2. Transfer the entire contents **of the tube (4 ml / 2 ml) to a 15 ml centrifuge tube.**
3. Add **(1 ml / 0,50 ml) of Lysis Solution + (15 µl / 7,50 µl) of RNase.** Vortex to mix well. Incubate at room temperature for 15 minutes.

Protein Precipitation.

1. Add **(1.60 ml / 0,80 ml)** of **Protein Precipitation Solution** to the cell lysate.
2. Vortex vigorously at maximum speed for 20-30 seconds.
3. Centrifuge at 4000 rpm for 5 minutes. A precipitate will form. If floating particles are observed, centrifuge again after incubating for 5 minutes on ice.

Precipitation of the DNA.

1. Transfer the supernatant containing the DNA **to a fresh 15 ml centrifuge tube** containing **(5 ml / 2,50 ml) Isopropanol.**
2. Mix by inversion about 25 times.
3. Centrifuge at 4000 rpm for 3 minutes. The DNA will be visible as a white pellet.
4. Remove the supernatant and dry the tube briefly on absorbent paper. Add **(5 ml / 2,50 ml) of 70% Ethanol** to wash the DNA.
5. Centrifuge at 4000 rpm for 1 minute. Carefully remove the supernatant without touching the DNA pellet. Centrifuge again briefly to collect the last drops of residual ethanol.
6. Invert the tube on absorbent paper and allow to dry for about 5-10 minutes.

DNA hydration.

1. Add **(750-1.000 µl / 350-500 µl)** of sterile water or **Hydration Solution** and resuspend by pipette.
2. Incubate at 50°C for 1 hour with periodic shaking to aid DNA rehydration, or incubate "overnight" at room temperature with gentle agitation.
3. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service.

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All components are stable for 12 months from date of purchase when stored and used as directed.

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