REAL SWABS DNA KIT

Ref. RBMEG20:	(para 100 EXTRACCIONES)	(for 100 EXTRACTION)
Ref. RBMEG21:	(para 500 EXTRACCIONES)	(for 500 EXTRACTION).
Ref. RBMEG22:	(para 1000 EXTRACCIONES)	(for 1000 EXTRACTION).

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REAL SWABS DNA Kit provides a method for the extraction

1. SAMPLE TAKING

of high quality genomic DNA from **buccal swab samples with foam brushes and saliva samples collected with our REALSWABS Sample Collection Kit system.**

It is a fast, safe and economical method. It uses a deproteinization with a novel saline buffer avoiding the use of toxic organic solvents such as phenol or chloroform.

In the case of saliva, obtaining DNA in the same individual presents an intrasubject variability, that is, there is not an approximate amount of DNA / ml since it can vary according to the time of collection of the sample, given the case of that a large amount of DNA is not obtained using swabs and subsequent applications that require a large amount of DNA cannot be carried out.

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Equipment and reagents needed and not provided

- * Isopropanol.
- * 70% ethanol.
- * Microtubes of 1.5 ml and 2.0 ml.
- * Microcentrifuge or clinical centrifuge.
- * Vortex.
- * Water bath

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
3.1 Preliminary preparations

Sample collection

1. It is recommended that the individual from whom the sample is to be taken refrain from drinking coffee and eating any food at least 30 minutes prior to collection. If this is not possible, a gentle rinse of the mouth with water only is recommended.
2. Collect the buccal cell sample with the foam swab. Rub the swab on the inside of the cheek (buccal wall) and gums with firm pressure about 30 times on each side of the face and each side of the swab. **Important: apply firm, solid and reasonable pressure with the swab. The hand can be placed on the cheek to provide a more solid surface.**
3. Proceed with extraction or the indicated preservation method depending on the swab used.
4. If using our **REALSWABS Sample Collection Kit** proceed as follows.
5. Insert the swab into the microtube with **the Preservative Solution**. Rotate the swab in a rapid manner to release the buccal cells into the solution **with the same stirring motion of a teaspoon**. Press the swab against the wall of the microtube and rotate to ensure that most of the liquid remains in the microtube. It is also recommended to pass the foam head of the swab several times over the top of the microtube to release the last drops.

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6. The solution should be observed to become slightly turbid indicating the presence of buccal cells. **If the solution remains clear, a new sample should be collected with a new swab at a different time.**
7. The preserved sample is stable at room temperature (15-25°C) for 1 year and can be sent to the laboratory for genomic DNA purification. If stored at -20°C or -80°C the sample is stable indefinitely.

If the Lysis Solution contains a precipitate due to low temperatures, incubate at 37°C and mix to dissolve the precipitate.

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3.2 Extraction from preserved specimens with our REALSWABS Sample Collection Kit

1. Those samples that have been left in an upright position for several days will show a white pellet containing the buccal cells. Using a micropipette resuspend the pellet completely and transfer the entire solution to a new 2.0 ml microtube.
2. Add **1000 µl of Resuspension Solution** and Centrifuge at 13,000-16,000 x g for 2 minutes.
3. Remove the supernatant by decanting **taking care not to lose the cell pellet**. Centrifuge again briefly and remove all the liquid by micropipetting.
4. Add **300 µl of Lysis Solution + 5 µl of Proteinase K + 5 µl of RNase** and resuspend the pellet with micropipette. Incubate at 37°C for 30 minutes. Vortex several times during the incubation period.
5. Add **100 µl of Protein Precipitation Solution**.
6. Vortex vigorously for 30 seconds.
7. Centrifuge at 13,000-16,000 x g for 5 minutes. The protein precipitate will be observed to form a very faint pellet (if there was little cell pellet, sometimes the protein pellet is not observed).
8. Transfer the supernatant containing the DNA to a 1.5 ml tube containing **300 µl of isopropanol**. Mix by inversion several times.
9. Centrifuge at 13,000-16,000 x g for 3 min.
10. Remove the supernatant. Add **300 µl of 70% ethanol** and invert several times to wash the DNA pellet. **Watch not to lose the DNA pellet which will be very small or undetectable and will be more visible after the 70% ethanol wash.**
11. Centrifuge at 13,000-16,000 x g for 2 minutes. Carefully remove all ethanol. Be careful not to lose the DNA pellet which will be very small or undetectable.
12. Invert the tube and allow to dry on absorbent paper for 5 minutes. Alternatively, pulse centrifuge and remove the liquid with a micropipette.
13. Add **30-40 µl of DNA Hydration Solution** and resuspend with micropipette.
14. Store at 2-8°C. For long term storage store at -20°C or -80°C.

3.3 Extraction from swab samples

The use of foam brushes such as those provided are suitable as they allow recovery of almost all of the **Lysis Solution** after the incubation period.

Cell lysis

1. Add **600 µl of Lysis Solution** into a 2.0 ml microtube. Cut the brush head with a bit of handle and insert into the microtube. Vortex vigorously to release the cells from the brush.
2. Add **5 µl of Proteinase K + 5 µl of RNasa** and incubate at 37°C for 30 minutes. **Vortex several times during the incubation period.**
3. Remove the brush head from the **Lysis solution**, rubbing it against the walls to collect the maximum amount of liquid.

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Protein Precipitation

1. Add **200 µl of Protein Precipitation Solution**
2. Vortex vigorously for 30 seconds.
3. Centrifuge at 13,000-16,000 x g for 5 minutes. The protein precipitate will be observed to form a pellet.

DNA Precipitation

1. Transfer the supernatant containing the DNA to a 1.5 ml tube containing **600 µl of isopropanol**. Mix by inversion several times.
2. Centrifuge at 13,000-16,000 x g for 3 minutes.
3. Remove the supernatant. Add **600 µl of 70% ethanol** and invert several times to wash the DNA pellet. **Watch not to lose the DNA pellet which will be very small or undetectable and will be more visible after the 70% ethanol wash.**
4. Centrifuge at 13,000-16,000 x g for 2 minutes. Carefully remove all ethanol. Be careful not to lose the DNA pellet.
5. Invert the tube and allow to dry on absorbent paper for 5 minutes.

DNA hydration

1. Add **30-40 µl of DNA Hydration Solution** and resuspend with micropipette.
2. Store at 2-8°C. For long term storage store at -20°C or -80°C.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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