# REAL SPIN FOOD-STOOL KIT

Ref. **RBMEGS05:** (para **50 EXTRACCIONES**)  
 Ref. **RBMEGS06:** (para **250 EXTRACCIONES**)

(for **50 EXTRACTION**)  
 (for **250 EXTRACTION**).

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## SAMPLE TAKING

This kit has been optimized for efficient and rapid **DNA purification from fresh or frozen stool and various food samples (raw materials and processed foods)**.

After homogenizing the samples, the DNA can be extracted with the Extraction Buffer, the lysis mixture is centrifuged to remove contaminants and residual debris. The supernatant is mixed with Binding Buffer + Proteinase K, incubated at 70°C for 10 minutes and then isopropanol to create the appropriate conditions for optimal DNA binding to the silica membrane of the column. Then 2 washes with different buffers for efficient removal of potential PCR inhibitors. Finally, the DNA is eluted and ready to be used in subsequent reactions.

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### Equipment and reagents needed and not provided

- \* 100% ethanol.
- \* 1.5 ml microtubes.
- \* Microcentrifuge or clinical centrifuge.
- \* Hand-held electric homogenizer and/or liquid nitrogen to pre-process the samples.
- \* Vortex.
- \* Bath for incubation at 80°C.

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**2.1 Preliminary preparations**

- Dissolve Proteinase K in **1.3 ml** (50 extractions) or in **2 x 3.35 ml** (250 extractions) in nuclease-free water and store at  $-20^{\circ}\text{C}$ . **It is recommended to make several aliquots** to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Verify that the Lysis/Binding Buffer does not have precipitates due to the low temperatures. If necessary, dissolve by heating to  $37^{\circ}\text{C}$ .
- **Add 100 % Ethanol to the Disinhibition Buffer** indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- **Add the 100% Ethanol to the Wash Buffer** indicated on the label, about **40 ml** (kit 50 extractions) and about **200 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- **Pre-heating the Elution Buffer to  $70^{\circ}\text{C}$  may increase the yield of DNA obtained. For some subsequent application it may be necessary for the DNA to be concentrated, elution in volumes smaller than 200  $\mu\text{l}$  will increase the final concentration of DNA in the eluate but will reduce the overall yield of DNA obtained. For samples containing  $< 3 \mu\text{g}$  DNA, elution in 100  $\mu\text{l}$  is recommended. For samples containing  $< 1 \mu\text{g}$  DNA, a 35-50  $\mu\text{l}$  elution is recommended.**

**2.2 Protocol for DNA extraction from 200 mg of fresh or dried feces**

**NOTE:** It is recommended to start with large amounts of stool when DNA is not homogeneously distributed or is found in small amounts in the sample. Smaller sample sizes than indicated can be processed especially when maximum removal of potential PCR inhibitors is required.

In general, for optimal PCR results use the minimum amount of eluted DNA possible, the volume of eluate should never exceed 10 % of the final volume of the PCR mixture. It is recommended to add BSA to a final concentration of  $0.1 \mu\text{g} / \mu\text{l}$  of the PCR mixture and use HOT Star polymerase.

1. Weigh **180-200 mg of feces** into a 2.0 ml microtube and add **1.20 ml of CTAB Extraction Buffer**. Vortex vigorously for 1 minute. If a homogeneous solution is not observed it will be necessary to use a hand held electric homogenizer.  
*If the sample is liquid pipette 200  $\mu\text{l}$  into a microtube.*
2. **Incubate at  $70^{\circ}\text{C}$  for 30-60 minutes.** Repeat vortexing several times during incubation.  
*For detection of human DNA it is sufficient to incubate at  $70^{\circ}\text{C}$ , for detection of bacteria or parasites incubate at  $80^{\circ}\text{C}$  and For detection of difficult to lyse cells, such as some bacteria or parasites, the incubation temperature can be increased to  $95^{\circ}\text{C}$  if necessary.*
3. **Centrifuge at 14,000 rpm for 10 minutes.** A pellet will appear and on the surface a layer of fat, introduce the pipette tip through this surface layer, trying to collect only **500  $\mu\text{l}$  of supernatant** which is the transparent liquid with color (avoid catching pellet and surface layer) and place in a 2.0 ml microtube.

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4. Add **25 µl Proteínasa K** to the 500 µl supernatant. Mix well. **Incubate at 70°C for 10 minutes.**
5. Add **250 µl of Binding Buffer** to the 500 µl of supernatant I. Mix well by vortexing.
6. Add to the Spin microcolumn reservoir with its collection tube. **Centrifuge at 10,000 rpm for 60 seconds.** Discard the collection tube.
7. Place the MicroSpin column in a new collection tube and add **500 µl of Deinhhibition Buffer to the reservoir.** **Centrifuge at 12,000 rpm for 60 seconds.** Discard the liquid.
8. **Add 700 µl of Wash Buffer** to the Spin microcolumn reservoir. **Centrifuge at 14,000 rpm for 60 seconds.** Discard the liquid.
9. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol.**
10. Remove the collection tube and insert the spin microcolumn into a 1.5 ml microtube. **Add 50-200 µl of Elution Buffer** (preheated to 70°C) to the spin microcolumn reservoir. Incubate for 2 minutes.
11. **Centrifuge at maximum speed for 60 seconds.** The microtube now contains the bacterial DNA.

### **2.3 Protocol for DNA extraction from food samples.**

Given the wide variety of samples that include food, plant origin, animal origin, processed food, raw materials, etc., it is difficult to present a universal protocol for all samples. This is why DURVIZ Technical Department can study and fine-tune your specific protocol for your particular type of sample.

1. Weigh **100-200 mg** of the sample in a 2.0 ml microtube and add **1.2 ml of CTAB Buffer + 25 µl Proteínasa K** Vortex vigorously.

**The main and most important step to obtain good yields is a good sample breakage and homogenization which will be specific for each type of sample. In all cases and for greater effectiveness, liquid nitrogen should be used to pulverize the sample.**

As a general rule for solid samples (sausages, sausages, etc.), prepare several fragments and homogenize with a hand-held electric homogenizer; For solid powder samples (flours, etc.) homogenize with a hand-held electric homogenizer; For large solid samples (corn flakes, chocolate, cookies, etc.) use a coffee grinder to pulverize a large sample and then weigh the required amount of powder; For liquid samples use 200 µl directly. You can also use the different rotors + particles available on the market.

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1. **Incubate at 65°C for 30 minutes.** Repeat vortexing several times during incubation.
2. **Centrifuge at 14,000 rpm for 5-10 minutes.** A pellet will appear and a fatty layer will appear on the surface, introduce the pipette tip through this surface layer, trying to collect only **500 µl of supernatant** which is the clear colored liquid (avoid catching pellet and surface layer) and place in a new 1.5 ml microtube.
3. Add **250 µl of Binding Buffer** to the 500 µl of supernatant. Mix well by vortexing.
4. Add to the Spin microcolumn reservoir with its collection tube. **Centrifuge at 10,000 rpm for 60 seconds.** Discard the collection tube.
5. Place the MicroSpin column in a new collection tube and add **500 µl of Deinhibition Buffer to the reservoir. Centrifuge at 12,000 rpm for 60 seconds.** Discard the liquid.
6. **Add 700 µl of Wash Buffer** to the Spin microcolumn reservoir. **Centrifuge at 14,000 rpm for 60 seconds.** Discard the liquid.
7. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol.**
8. Remove the collection tube and insert the spin microcolumn into a 1.5 ml microtube. **Add 50-200 µl of Elution Buffer** (preheated to 70°C) to the spin microcolumn reservoir. Incubate for 2 minutes.
9. **Centrifuge at maximum speed for 60 seconds.** The microtube now contains the DNA.

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#### STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

**CAUTION:** Lysis/Union Buffer and Deinhibition Buffer contain guanidine hydrochloride which is irritating, use gloves and goggles.

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