REAL SPIN VIRAL DNA/RNA

Ref. **RBMEGS07:** (para **100 EXTRACCIONES**)(for **100 EXTRACTION**)

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1.INTRODUCTION

REALPURE Spin Viral DNA/RNA Kit is designed for rapid simultaneous purification of viral DNA and RNA from cell-free samples such as serum, plasma and cerebrospinal fluid.

Viruses when lysed with a detergent and proteinase K, release viral nucleic acids. Then, in the presence of chaotropic salts, the viral nucleic acids selectively bind to a silica membrane on a spin column. The nucleic acids remain bound while a series of washes and centrifugation steps will remove contaminating cellular components. Finally, with an elution buffer the viral nucleic acids will be eluted from the membrane. This process does not require alcohol precipitation, organic solvent extractions, or extensive nucleic acid handling.

The REALPURE Spin Viral DNA/RNA Kit can be used for the extraction of viral DNA and RNA from a wide range of DNA and RNA viruses. However, success cannot be guaranteed for all virus species and must be validated by the user.

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Equipment and reagents needed and not provided

- * 100% ethanol.
- * Microtubes of 1.5 ml.
- * Microcentrifuge or clinical centrifuge

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3.1 Preliminary considerations

Viral Lysis Buffer (REA01V) and Viral Wash Buffer 1 (REA03V) contain guanidine hydrochloride, which can form reactive components when combined with bleach. Both buffers are irritants, for this reason we recommend the use of gloves and goggles for handling. In case of contact with skin or eyes, wash with plenty of water.

3.2 Preliminary preparations

- Dissolve Proteinase K in **5.2 ml** of nuclease-free water and store at -20°C . It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Add **20 ml** 100% ethanol to Viral Wash Buffer 1 (REA03V). Keep the container tightly closed to prevent evaporation of the ethanol.
- Add **80 ml** 100 % Ethanol to Viral Wash Buffer 2 (REA04V). Keep the container tightly closed to avoid evaporation of ethanol.
- **Carrier RNA: The recommended purification protocol uses 5.6 µg carrier RNA per sample.** If you wish to use less carrier RNA per sample, it is necessary to validate the amount of carrier RNA needed for each sample type and subsequent applications.
- The eluates obtained with this kit contain both viral nucleic acids and carrier RNA, **and the amounts of carrier RNA may exceed the amounts of viral nucleic acids.** Calculations of how much eluate should be added to subsequent applications should be based on the amounts of carrier RNA added. To obtain high levels of sensitivity in amplification reactions, it may be necessary to adjust the amount of carrier RNA added to the Viral Lysis Buffer.

➤ **Carrier RNA preparation protocol:**

- 1 Add 500 µl of Elution Buffer to the microtube containing the carrier RNA to obtain a stock solution
- 2 Mix and aliquot in 25-50 µ, store at -20°C . Avoid repeated heating-cooling cycles.
- 3 Calculate the volume of Viral Lysis Buffer/carrier RNA mixture required to process all samples to be analyzed simultaneously. For example, add 50µL of carrier to 2.5 ml of Viral Lysis Buffer.
- 4 In a sterile tube add the volume of stock solution to the volume of Viral Lysis Buffer required. Mix gently with a pipette.
- 5 Store at 4°C until use. **Use this buffer thus prepared within 1 hour.**

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3.3 Protocol for the extraction of viral DNA and RNA from serum, plasma and biological fluids

- 1 Pipette **50 µl Proteinase K** into a sterile centrifuge microtube.
- 2 **Add 200 µL of plasma or serum** into the microtube. If you process samples smaller than 200 µl, adjust the final volume to 200 µl using PBS or 0.9% NaCl.
- 3 **Add 200 µl of Viral Lysis Buffer (containing the RNA carrier)**. Close the microtube and vortex vigorously for 20 seconds.
- 4 **Incubate at 70°C for 10 minutes.**
- 5 **Briefly centrifuge the microtube to collect the droplets from the cap.**
- 6 **Add 250 µl 100% ethanol to the sample.** Vortex the sample for 15 seconds.
- 7 **Add the lysate to the Viral Spin Column** with a collection tube.
- 8 **Centrifuge the column at 8,000-10,000 rpm for 1 min.** Discard the collection tube. Place the spin column in a new collection tube.
- 9 **Wash the column with 500 µl of Viral Wash Buffer 1 with ethanol.** Centrifuge the column at 12,000 rpm for 1 min.
- 10 **Wash the column with 500 µl Viral Wash Buffer 2 with ethanol.** Centrifuge the column at 12,000 rpm for 1 min.
- 11 **Wash the column again with 500 µl of Viral Wash Buffer 2 with ethanol.** Centrifuge the column at 14,000rpm for 1 min.
- 12 **Centrifuge the spin column at maximum speed for 2 min to remove residual ethanol.**
- 13 **Elute with 30-50 µl of Elution Buffer by adding the nuclease-free water to the center of the column membrane.**
- 14 **Incubate at room temperature for 2 minutes.**
- 15 **Centrifuge the spin column at 10,000rpm for 60 seconds.** Collect the 30-50 µl and redeposit in the center of the membrane. This increases the yield.
- 16 Incubate for 2 minutes and centrifuge at maximum speed. If problems are observed in subsequent detections, the volume of eluate added to the PCR or RT-PCR can be changed.

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: Lysis/Union Buffer and Disinhibition Buffer contain guanidine hydrochloride which is irritating, use gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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