REAL SPIN BLOOD DNA KIT

Ref. **RBMEGS08:** (para **50 EXTRACCIONES**)
Ref. **RBMEGS09:** (para **250 EXTRACCIONES**)

(for **50 EXTRACTION**)
(for **250 EXTRACTION**)

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1.INTRODUCTION

This kit is designed for rapid extraction and purification of high quality **genomic DNA from whole blood, serum, plasma and biological fluids** using Spin columns with silica membrane that selectively binds DNA.

This Kit utilizes a new Lysis / Binding Buffer specifically formulated for DNA extraction from blood samples for high throughput.


Features:

- For rapid obtaining of high quality, ready-to-use DNA from blood.
- Sample size: 300 µl of whole blood, plasma, serum and biological fluids.
- No organic extractions or alcohol precipitations are used.
- Complete removal of inhibitors or contaminants.
- Typical yield: 6-9 µg genomic DNA.
- Elution volume: 50-200 µl.
- High quality DNA is obtained that can be used directly in PCR, Southern, cloning and any enzymatic reaction.

Applications:

- Genomic, viral and bacterial DNA extractions.
- DNA from whole blood (human or animal blood, fresh or frozen).
- DNA from blood treated with citrate, EDTA, heparin.
- DNA from serum, plasma, platelets, buffy coat, biological fluids and dried blood spots.

Procedure: In the case of blood samples and if genomic DNA is required, the best option is a previous and selective lysis of erythrocytes with RBC Lysis Buffer to process only lymphocytes which will provide better results in terms of quality and yield. If we are looking for other types of DNA (bacterial, viral) lysis is achieved by incubating the whole blood in a chaotropic solution in the presence of proteinase K at 70°C. The appropriate conditions for DNA to bind to the silica fiber membrane will be created by adding ethanol to the lysate. Contaminants are removed by 2 different washes and finally the DNA is eluted with elution buffer.

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Equipment and reagents needed and not provided

- * 100% ethanol.
- * Microtubes of 1.5 ml.
- * Microcentrifuge or clinical centrifuge

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3.1 Specimen collection and storage

Whole blood samples should be stored at 4°C immediately after collection. They are stable for weeks at 4°C. They can also be shipped in refrigerated containers.

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Plasma and serum samples should be cooled and centrifuged immediately within one hour of collection. For blood serum preparation centrifuge at 3000 rpm for 10 minutes and for plasma samples centrifuge at 3500 rpm for 15 minutes.

The plasma samples should be separated from the cells and transferred to a new 1.5 ml microtube, the intermediate layer that includes the white cells and platelets should not be transferred with the plasma.

Serum samples are usually obtained with plain glass tubes without anticoagulants that allow blood clots to form after the serum is collected and transferred to a new microtube for transport.

If collection is not possible within three days of collection, the plasma and serum should preferably be frozen at -80°C or at least -20°C.

Buffy coat is a leukocyte-enriched fraction from whole blood. Preparing this fraction from whole blood is simple and has a yield of approximately 5-10 times more DNA than an equivalent volume of whole blood. It is prepared by centrifugation of the blood at 2500 x g for 10 minutes at room temperature. After centrifugation, 3 fractions are observed. The top fraction is the plasma; the middle fraction is the byffy coat containing the concentrated leukocytes; and the bottom fraction contains the erythrocytes.

3.2 Preliminary preparations

- Dissolve Proteinase K in **1.3 ml** (50 extractions) or in **2 x 3.35 ml** (250 extractions) in nuclease-free water and store at -20°C. Several aliquots are recommended to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Verify that the Lysis/Binding Buffer does not have precipitates due to the low temperatures. If necessary, dissolve by heating to 37°C.
- Add 100 % Ethanol to the Disinhibition Buffer indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Add the 100% Ethanol to the Wash Buffer indicated on the label, about **40 ml** (kit 50 extractions) and about **200 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Pre-heat the Elution Buffer to 70°C.

3.3 Protocol for genomic DNA extraction from blood (leukocytes)

This protocol is for the purification of genomic DNA from leukocytes by selective lysis of erythrocytes with RBC Lysis Buffer. If only genomic or mitochondrial DNA is required this is the ideal method as it produces better results in terms of quality and yield.

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1. Pipette **300 µl of blood** into a 1.5 ml microtube. Add **900 µl of RBC Lysis Buffer**. Vortex and incubate at room temperature for 10 minutes.
2. Centrifuge at maximum speed for 1 minute. Remove the supernatant by decanting better than with micropipette which can aspirate the small cell pellet not visible and leave 10-20 µl of residual liquid. Vortex the microtube to resuspend the pellet.
3. Add **180 µl of Tissue Lysis Buffer + 25 µl Proteinase K**. Resuspend the cell pellet with micropipette.
4. **Incubate at 56°C** for 15 minutes.
5. Add **200 µl of Lysis / Binding Buffer**. Mix by vortexing. **Incubate at 70°C for 10 minutes**.
6. Add **200 µl of Ethanol (96-100%) to the lysate**. Mix by vortex.
7. **Transfer the sample to a MicroSpin column** with its collection tube.
8. **Centrifuge at 8,000 rpm for 60 seconds**. Discard the collection tube. If the sample has not passed completely, repeat the centrifugation step.
9. Place the MicroSpin column in a new collection tube and add **500 µl of Deinhibition Buffer**. **Centrifuge at 12,000-14,000 rpm for 60 seconds**. Discard the liquid.
10. **Add 700 µl of Wash Buffer** to the Spin column reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Discard the liquid.
11. **Centrifuge at maximum speed for 2 minutes to remove residual ethanol**.
12. Remove the collection tube and insert the spin column into a 1.5 ml microtube. **Add 50-200 µl of Elution Buffer** (pre-warmed to 70°C) to the spin microcolumn reservoir. **Incubate for 1 minute**.
13. **Centrifuge at maximum speed for 60 seconds**. The microtube now contains the genomic DNA.

3.4 Protocol for DNA extraction from whole blood, buffy coat and biological fluids

This protocol is for the purification of total DNA (genomic, mitochondrial, bacterial and viral) from whole blood, buffy coat and biological fluids.

1. **Pipette 25 µl proteinase K** into a 1.5 ml microtube.
2. **Add 300 µl of sample** to the microtube. Use up to 300 µl of whole blood, plasma, serum, buffy coat or biological fluid. For samples smaller than 300 µl, add PBS to adjust the volume to 300 µl. If purifying viral DNA, we recommend starting with 200 µl of serum or plasma.
3. **Add 300 µl of Lysis/Binding Buffer**. Vortex the mixture vigorously (10-20 s). **Incubate at 70°C for 15 minutes**. The lysate will turn brown during incubation. Increase the incubation period with proteinase K to 30 minutes and vortex vigorously more times during the incubation process if processing old or clotted blood.
4. **Add 300 µL ethanol (96-100 %)** to each sample and vortex.
5. Pipette half of the lysate into the Spin column reservoir with its collection tube. **Centrifuge at 10,000 rpm for 60 seconds**. Discard the collection tube.
6. Repeat step 5 with the other half of the lysate.

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7. Place the Spin column in a new collection tube and add **500 µl of Deinhhibition Buffer** to the reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Discard the liquid.
8. **Add 700 µl of Wash Buffer** to the Spin column reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Remove the liquid.
9. **Centrifuge at maximum speed for 2 minutes to remove residual ethanol**.
10. Remove the collection tube and insert the spin column into a 1.5 ml microtube. **Add 50-200 µl of Elution Buffer** (pre-warmed to 70°C) to the reservoir of the spin microcolumn. **Incubate for 1 minute**.
11. **Centrifuge at maximum speed for 60 seconds**. The microtube now contains the genomic DNA.

3.5 Protocol for extraction of genomic DNA from dried blood spots

Step Sample collection

Place 3 circles obtained from puncturing a dried blood spot in a 1.5 ml microtube.

Pre-treatment step

Add 200 µl of PBS and vortex vigorously. Incubate at 85°C for 10 minutes. Briefly centrifuge to remove droplets from the stopper.

1. **Pipette 25 µl proteinase K** into the sample.
2. **Add 200 µl of Lysis/Binding Buffer**. Vortex the mixture vigorously (10-20 s).
3. **Incubate at 70°C for 1 hr**.
4. **Add 200 µl of ethanol (96-100 %)** to each sample and vortex. Continue at point 5 of protocol 3.4.

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: Lysis/Union Buffer and Disinhibition Buffer contain guanidine hydrochloride which is irritating, use gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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