REAL SPIN BLOOD DNA KIT

Ref. **RBMEGS10:** (para **50 EXTRACCIONES**)
 Ref. **RBMEGS11:** (para **250 EXTRACCIONES**)

(for **50 EXTRACTION**)
 (for **250 EXTRACTION**)

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1.INTRODUCTION

This kit allows efficient extraction of **genomic and mitochondrial DNA from small sample sizes such as different cell types and tissues, laser microdissected samples**, small amounts of blood using a special column design,

The special column design is connected with a reduced dead volume and a small diameter of the silica membrane which **allows an elution in only 10 µl.**

SAMPLE	SIZE
Tissues	< 10 mg
Cells in culture	< 10 ⁵
Blood	< 100 µl
Laser microdissected tissues	A sample
"Buccal Swabs"	1

Features:

- Silica membrane technology but with a special design of MicroSpin columns.
- Rapid purification of high quality DNA from small sample sizes.
- No organic extractions or alcohol precipitations are used.
- Complete removal of inhibitors or contaminants.
- Elution volume: 10-30 µl.
- High quality DNA is obtained that can be used directly in PCR, Southern, cloning and any enzymatic reaction.

Applications:

- DNA extraction from tissues (human or mouse tissue, laser microdissected tissues).
- DNA extraction from cells (cells in culture).
- DNA extraction from clinical samples (blood samples, biopsies).
- DNA extraction from forensic samples ("buccal swabs").

2.COMONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase. (50 Preps)	Envase. (250 Preps)	T°
RBC Lysis Solution	E20T	8 ml	40 ml	RT
Lysis Solution	E21	10 ml	50 ml	RT
Washing solution *	EP08	6 ml	30 ml	RT
Lysis/Binding Solution	REA01	10 ml	50 ml	RT
Proteinase K *	REA02	22 mg	105 mg	-20°C
Disinhibition Solution (GHCL 8M) *	REA03	18 ml	85 ml	RT
Hydration Solution	REA05	2 ml	10 ml	RT
microSpin columns	RSC06	50 unid.	250 unid.	RT
Collection Tubes	R30	100 unid.	500 unid.	RT

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* See the section on preliminary preparations on how to prepare these solutions.

Equipment and reagents needed and not provided

- * 100% ethanol.
- * Microtubes of 1.5 ml.
- * Microcentrifuge or clinical centrifuge

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3.1 Preliminary preparations

- Dissolve **Proteinase K** in **1.1 ml** of nuclease-free water (kit 50 extractions) and in **5.2 ml** (kit 250 extractions) and store at -20°C. It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Verify that the **Lysis Solution**, **Lysis/Binding Solution** do not have precipitates due to low temperatures. If necessary, dissolve by heating to 37°C.
- Add 100% Ethanol to the **Disinhibition Solution (GHCL 8M)** indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Add 100% Ethanol to the **Washing Solution** indicated on the label, about **24 ml** (kit 50 extractions) and about **120 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Pre-heat the **Hydration Solution** to 70°C.

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3.2 Protocol for extraction of genomic DNA from blood

1. Pipette **30 or 50 µl of blood** into a 1.5 ml microtube. Add **150 µl of RBC Lysis Solution**. Vortex and incubate at room temperature for 10 minutes.
2. Centrifuge at maximum speed for 1 minute. Remove the supernatant by decanting better than with micropipette which can aspirate the small non-visible cell pellet and leave 10-20 µl of residual liquid. Vortex the microtube to resuspend the pellet.
3. Add **180 µl of Lysis Solution + 20 µl Proteinase K**. Mix by vortexing for 2-5 seconds.
4. **Incubate at 56°C** for 10 minutes.
5. Add **200 µl of Lysis/Binding Solution**. Mix by vortexing. Incubate at 70°C for 10 minutes.
6. Add **200 µl of Ethanol (96-100%) to the lysate**. Mix by vortex.
7. **Transfer the sample to a MicroSpin Columns** with its *Collection Tubes*.
8. **Centrifuge at 8,000 rpm for 60 seconds**. Discard the *Collection Tubes*. If the sample has not passed completely, repeat the centrifugation step.
9. Place the *MicroSpin columns* in a new *Collection Tubes* and add **500 µl of Desinhibition Solution (GHCL 8M)**. **Centrifuge at 12,000-14,000 rpm for 60 seconds**. Discard the liquid.
10. **Add 500 µl of Whasing Solution**. **Centrifuge at 14,000 rpm for 60 seconds**. Remove the liquid.
11. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol**.
12. Place the *MicroSpin column* in a 1.5 ml microtube and add **10-50 µl of Hydration Solution (10 mM Tris.HCl, pH 8.5)**. Ensure that the *Hydration Solution* is dispensed directly into the center of the membrane for complete elution of bound DNA. The average eluate volume is 10 µl from 12 µl of *Hydration Solution* used.
13. **Incubate for 1 minute. Centrifuge at maximum speed for 1 minute**. The microtube now contains the genomic DNA.

3.3 Protocol for extraction of genomic DNA from cultured cells

1. Resuspend up to **10⁵ cells** in a final volume of **80 µl of Lysis Solution + 10 µl Proteinase K**. Vortex mix for 2-5 seconds.
2. **Incubate at 56°C for 10 minutes**.
3. Add **80 µl Lysis/Binding Solution**. Vortex mix. **Incubate at 70°C for 5 minutes**. Allow the lysate to cool to room temperature.
4. Add **80 µl of Ethanol (96-100%) to the lysate**. Vortex mix. Briefly centrifuge to collect any droplets present in the microtube.
5. **Transfer the sample to a MicroSpin columns** with its *Collection Tubes*.
6. **Centrifuge at 8000 rpm for 60 seconds**. Remove the *Collection Tubes*. If the sample has not passed completely, repeat the centrifugation step.
7. Place the *MicroSpin columns* in a new collection tube and add 80 µl of *Disinhibition Solution (GHCL 8M)*.

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8. **Centrifuge at 12,000-14,000 rpm for 60 seconds.** Eliminate the liquid.
9. **Add 80 µl of *Washing Solution*.** Centrifuge at **12,000-14,000 rpm for 60 seconds.** Remove the liquid.
10. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol.**
11. Place the ***MicroSpin column*** in a 1.5 ml microtube and add **10-50 µl of *Hydration Solution* (10 mM Tris.HCl, pH 8.5).** Ensure that the ***Hydration Solution*** is dispensed directly into the center of the membrane for complete elution of bound DNA. The average eluate volume is 10 l from 12 l of elution buffer used.
12. **Incubate for 1 minute. Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

3.4 Protocol for extraction of genomic DNA from less than 10 mg of animal tissue

1. Transfer a tissue sample of less than 10 mg to a 1.5 ml microtube.
2. Add **180 µl of *Lysis Solution* + 20 µl *Proteinase K*.** Vortex mix for 2-5 seconds.
3. **Incubate at 56°C** for 1-4 hours or until lysis is complete, samples can be incubated overnight. All samples and especially samples that are difficult to lyse can be ground with liquid Ni or can be treated directly with a Polytron type homogenizer.
4. Add **200 µl of *Lysis/Binding Solution*.** Mix by vortexing. **Incubate at 70°C for 10 minutes.** If insoluble particles are observed, centrifuge 5 minutes at maximum speed and transfer the supernatant to a new microtube.
5. Add **200 µl of Ethanol (96-100%) to the lysate.** Mix by vortexing. Briefly centrifuge to collect any droplets present in the microtube.
6. **Transfer the sample to a *MicroSpin column*** with collection tube.
7. **Centrifuge at 8,000 rpm for 60 seconds.** Discard the ***Collection Tubes***. If the sample has not passed completely, repeat the centrifugation step.
8. Place the ***MicroSpin column*** in a new collection tube and add **500 µl of *Desinhibition Solution* (GHCL 8M)..**
9. **Centrifuge at 12,000-14,000 rpm for 60 seconds.** Remove the liquid.
10. **Add 500 µl of *Washing Solution*.** Centrifuge at **12,000-14,000 rpm for 60 seconds.** Remove the liquid.
11. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol.**
12. Place the ***MicroSpin columns*** in a 1.5 ml microtube and add **10-50 µl of *Hydration Solution* (10 mM Tris.HCl, pH 8.5).** Ensure that the ***Hydration Solution*** is dispensed directly into the center of the membrane for complete elution of bound DNA. The average eluate volume is 10 µl from 12 µl of elution buffer used.
13. **Incubate for 1 minute. Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

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3.5 Protocol for extraction of genomic DNA from laser microdissected tissue

NOTE: This protocol for extraction of genomic DNA from laser microdissected tissue is quite difficult, the sample is very small and DNA quality may be adversely affected by the fixation and staining processes, and it may be necessary to either modify fixation protocols or use cryosections to minimize this problem.

1. Transfer a laser microdissected tissue sample to a microtube.
2. Add **80 µl of Lysis Solution + 10 µl Proteinase K**. Vortex mix for 2-5 seconds.
3. **Incubate at 56°C** for 1-4 hours or until lysis is complete, samples can be incubated overnight.
4. Add **80 µl of Lysis / Binding Solution**. Mix by vortexing. **Incubate at 70°C for 5 minutes**. Allow the lysate to cool to room temperature.
5. Add **80 µl of Ethanol (96-100%) to the lysate**. Vortex mix. Briefly centrifuge to collect any droplets present in the microtube.
6. **Transfer the sample to a MicroSpin columns** with its collection tube.
7. **Centrifuge at 8000 rpm for 60 seconds**. Discard the collection tube. If the sample has not passed completely, repeat the centrifugation step.
8. Place the **MicroSpin columns** in a new collection tube and add **80 µl of Disinhibition Solution (GHCL 8M)**.
9. **Centrifuge at 12,000-14,000 rpm for 60 seconds**. Remove the liquid.
10. **Add 80 µl of Washing Solution**. **Centrifuge at 12,000-14,000 rpm for 60 seconds**. Remove the liquid.
11. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol**.
12. Place the **MicroSpin columns** in a 1.5 ml microtube and add **10-30 µl of Hydration Solution (10 mM Tris.HCl, pH 8.5)**. Ensure that the **Hydration Solution** is dispensed directly into the center of the membrane for complete elution of bound DNA. The average eluate volume is 10 µl from 12 µl of **Hydration Solution** used.
13. **Incubate for 1 minute**. **Centrifuge at maximum speed for 1 minute**. The microtube now contains the genomic DNA.

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3.6 Protocol for extraction of genomic DNA from buccal swabs

Sample collection

1. It is recommended that the individual from whom the sample is to be taken refrain from drinking coffee and eating any food at least 30 minutes prior to collection. If this is not possible, a gentle rinse of the mouth with water only is recommended.
2. Collect the buccal cell sample with the swab. Rub the swab on the inside of the cheek (buccal wall) and gums with firm pressure about 20 times on each side of the face and each side of the swab.
3. Use immediately for extraction. If the specimen is to be transported, allow the brush to dry at room temperature for 30 minutes. Then insert the brush into the receptacle provided for shipping. The sample can remain in this container tube for 1 week at 22-37°C before extraction. For long term storage keep the sample in the container at -20°C for up to 6 months.
4. Add **250 µl of RBC Lysis Solution + 20 µl of Proteinase K (20 mg/ml)** into a 1.5 ml microtube. Cut the brush head with a bit of handle and insert into the microtube. Vortex vigorously to release the cells from the brush. Ensure that the sample is completely covered with buffer during incubation.
5. **Incubate at 56°C for 30-60 minutes.** Centrifuge briefly. Remove the brush head from the lysis solution, rubbing it against the walls to collect the maximum amount of liquid. **Measure the amount of lysate.**
6. Add **1 volume of Lysis/Binding Solution** equal to the amount of lysate. Mix by vortexing. **Incubate at 70°C for 10 minutes.**
7. Add **1 volume of 96-100% Ethanol** equivalent to the amount of lysate. Vortex.
8. **Transfer the sample to a MicroSpin columns** with its **Collection Tubes.**
9. **Centrifuge at 8,000 rpm for 60 seconds.** Discard the **Collection Tubes.** If the sample has not passed completely, repeat centrifugation step.
10. Place the **MicroSpin columns** in a new **Collection Tubes** and add **200 µl of Disinhibition Solution.**
11. **Centrifuge at 12,000-14,000 rpm for 60 seconds.** Discard the liquid.
12. **Add 200 µl of Washing Solution. Centrifuge at 12,000-14,000 rpm for 60 seconds.** Remove the liquid.
13. **Centrifuge at full speed for 3 minutes to remove residual ethanol.**
14. Place the **MicroSpin columns** in a 1.5 ml microtube and add **10-30 l of Hydration Solution (10 mM Tris.HCl, pH 8.5).** Ensure that the **Hydration Solution** is dispensed directly into the center of the membrane for complete elution of bound DNA. The average eluate volume is 10 l from 12 l of **Hydration Solution** used.
15. **Incubate for 1 minute. Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

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- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: Lysis/Union Buffer and Disinhibition Buffer contain guanidine hydrochloride which is irritating, use gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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