REAL SPIN PLANT DNA KIT

Ref. **RBMEGS13:** (para **50 EXTRACCIONES**)
 Ref. **RBMEGS14:** (para **250 EXTRACCIONES**)

(for **50 EXTRACTION**)
 (for **250 EXTRACTION**)

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
1.INTRODUCTION

This kit provides a method for an efficient and fast method for **genomic DNA extraction from tissues from plants and fungi using MiniSpin columns.**

The kit includes two optimized, alternative lysis solution based on the established CTAB and SDS lysis methods. As plants are very heterogenous and contain a lot of different metabolites like polyphenols, polysaccharides, or acidic components, REALPURE SPIN "Plants and Fungi" **offers two different lysis procedures for optimal processing of various samples.**

In addition we also **use a PVP solution that can bind the polysaccharides and polyphenols.**

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Plant samples are first homogenized and then lysed in a highly optimized solution system, containing chaotropic salt, denaturing agent and detergents. A choice of two lysis solutions based on the established CTAB or SDS method are provided.

Crude lysate is cleared by centrifugation and the cleared lysate is then mixed with the Lysis/Binding Solution and processed through a Columns Spin DNA containing a silica membrane to which the plant genomic DNA binds. Contaminants and impurities such as salts, metabolites and cellular components are removed by simple Washing Solution steps with two different Soluciones.

High-quality purified plant genomic DNA is then eluted in a low *Hydration Solution*. The DNA is ready-to-use for a wide variety of applications.

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Features:

- Silica Membrane Technology using columns Spin DNA.
- Choice of two optimized lysis solution and PVP solution.
- High-purity DNA: typical A_{260}/A_{280} ratio 1.6 - 1.9.
- Plant genomic DNA isolated in 30 minutes

Applications:

- Isolation of genomic DNA from: fresh/frozen/lyophilized plant tissue and fungi.
- Isolated DNA is ready for downstream applications such as PCR, real-time PCR and genotyping.

2.COMONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase (50 Preps)	Envase (250 Preps)	T°
Neutralization Solution	E08	12 ml	60 ml	4°C
Proteinase K *	E40	20 mg	100 mg	-20°C
Washing solution *	EP08	10 ml	50 ml	RT
Plants and Fungi Extraction Solution	EPH1	50 ml	250 ml	RT
Lysis Solution	EPH2	12 ml	60 ml	RT
PVP Solution	EPH3	10 ml	50 ml	4°C
Lysis/Binding Solution	REA01	35 ml	175 ml	RT
Disinhibition Solution (GHCL 8M) *	REA03	16,5 ml	82,5 ml	RT
Hydration Solution	REA05	10 ml	50 ml	RT
CTAB Extraction Solution	RES03	60 ml	300 ml	RT
Columns Spin DNA	RSC03	50 unid.	250 unid.	RT
Collection Tubes	R30	100 unid.	500 unid.	RT

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Equipment and reagents needed and not provided

- * Isopropanol
- * Ethanol 100%
- * Ethanol 70 %
- * Microcentrifuge
- * 1.5 ml and 2.0 ml microtubes
- * Electronic homogenization
- * Water bath
- * Vortex
- * Liquid nitrogen and pestle and mortar to grind tissues

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3.1 Preliminary preparations

- Dissolve *Proteinase K* in **1.0 ml** of nuclease-free water (kit 50 extractions) and in **5.0 ml** (kit 250 extractions) and store at -20°C . It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Verify that the *Lysis/Binding Solution* and *CTAB Extraction Solution* do not have precipitates due to low temperatures. If necessary, dissolve by heating to 37°C .
- Add 100% Ethanol to the *Disinhibition Solution* indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Add 100% Ethanol to the *Washing Solution* indicated on the label, about **40 ml** (kit 50 extractions) and about **200 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Pre-heat the *Hydration Solution* to 70°C .
- Both the *Lysis/Binding Solution* and the *Disinhibition Solution* contain irritating salts, for this reason we recommend the use of gloves and goggles.

3.2 Preliminary Considerations

- It can be necessary to reduce the starting quantity of sample depending on the specie, state, tissue preparation or genome size.
- The first time that one works with a sample of plant tissue, we recommend to use the 2 isoaltion methods (A and B) for establishing which is the method that works better with this certain sample.

METHOD A. Based on Lysis Solution SDS

1. Weigh **180-200 mg of plant tissue**. Add **1000 μl of *Plants and Fungi E.S* + 200 μl *PVP solution*** and place them in 2.0 ml microtube. Vortex vigorously. Homogenize with an electric homogenizer for 20 seconds.

For an effective DNA isolation from **plant tissues**, it is recommended to grind the sample a fine powder with liquid nitrogen. Soft and not fibrous tissues, like young leaves, flowers, etc. can be homogenized in an electrical homogenizer but always when you have fresh samples, dissect the sample in small pieces, add ***Plants and Fungi E.S* + *PVP Solutio*** and homogenize. For hard or fibrous tissues, like seeds, steams, etc. It's recommended to use liquid Ni.

In case of **fungi**, collect the micellium from the culture by filtration, wash 3 times with deionized sterile water or PBS to remove the culture medium. Grind the sample a fine powder with liquid Ni, store at -80°C or process the sample.

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2. Add **200 µl Lysis Solution.** **Incubate at 80°C for 30-45 minutes.** If possible, shake with vortex or invert during the incubation. Verify that the **Lysis Solution** do not have precipitates due to the low temperatures. If necessary, dissolve heating at 37°C.
3. **Centrifuge at 14.000 rpm for 5 minutes.**
4. Pour the supernatant to a new 2.0 ml microtube with **200 µl Neutralization Solution.** **Incubate at -20°C for 3 minutes.**
5. **Centrifuge at 14.000 rpm for 5 minutes.** A pellet will appear and in the surface a layer of fat, introduce the pipette tip crossing this superficial layer, only trying to pick up **600-700 µ l of supernatant** that it is the transparent liquid with color (to avoid to catch pellet and superficial layer) and to place in a 2.0 ml microtube.
6. Pass the supernatant (**between 600-700 µl**) to a new microtube.
7. Add **600 µl-700 Lysis/Binding Solution + 20 µl Proteinase K** to 600-700 µl of supernatant. Mix well. **Incubate at 70°C for 10 minutes.**
8. Add **300 µl of Isopropanol.** Mix well.
9. Pipette **750 µl** the lysate into reservoir of a combined **Column Spin DNA – Collection Tubes** assembly. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the **Collection Tubes.**
10. Repeat step N.8 with the **750 µl** remaining. If there is more volume it can be passed through the column a third time.
11. Place the **Column Spin DNA** in a new **Collection Tubes** and add **500 µl of Disinhibition Solution** to the reservoir. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
12. **Add 500 µl of Washing Solution** into reservoir of **Column Spin DNA.** **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
13. 2º Wash. **Add 500 µl of Washing Solution** into reservoir of **Column Spin DNA.** **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
14. **Centrifuge at maximum speed for 2 minutes to remove the residual ethanol.**
15. Remove the **Collection Tubes** and insert the **Column Spin DNA** in a 1.5 ml microtube. **Add 50µl Hydration Solution** (preheated at 70°C) into reservoir of **Column Spin DNA.** **Incubate 1 minute.**
16. **Centrifuge at maximum speed for 60 seconds.**
17. **Add another 50µl of Hydration Solution** (preheated at 70°C) into reservoir of **Column Spin DNA.** **Incubate 1 minute.**
18. **Centrifuge at maximum speed for 60 seconds.** The microtube contains now genomic DNA.

METHOD B. Based on Lysis Buffer CTAB

1. Weight **180-200 mg of plant tissue.** **Add 1.20ml of CTAB Extraction Solution** and place them in 2.0 ml microtube. Shake with vortex vigorously.
*For an effective DNA isolation from **plant tissues**, it is recommended to grind the sample a fine powder with liquid nitrogen. Soft and not fibrous tissues, like young leaves, flowers, etc. can be homogenized in an electrical homogenizer but always when you have fresh samples, dissect the sample in small pieces, add the **CTAB Extraction Solution** and homogenize. For hard or fibrous tissues, like seeds, steams, etc. It's recommended to use liquid Ni.*

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In case of **fungi**, collect the micellium from the culture by filtration, wash 3 times with deionized sterile water or PBS to remove the culture medium. Grind the sample a fine powder with liquid nitrogen, store at -80°C or process the sample.

2. **Incubate at 80°C for 30-45 minutes.** If possible, shake with vortex or invert during the incubation.
3. **Centrifuge at 14.000 rpm for 5 minutes.** A pellet will appear and in the surface a layer of fat, to introduce the pipette tip crossing this superficial layer of fat, only trying to pick up **700 μl of supernatant** that it is the transparent liquid with color (to avoid to catch pellet and superficial layer) and to place in a 2.0 ml microtube.
4. Add **700 Lysis/ Binding Solution + 20 μl Proteinase K** to **700 μl of supernatant**. Mix well. **Incubate at 70°C for 10 minutes.**
5. Add **300 μl of Isopropanol**. Mix well.
6. Pipette **750 μl** the lysate into reservoir of a combined **Column Spin DNA**. –**Collection tube** ssembly. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the **Collection tube**.
7. Repeat step N.6 with the **750 μl** remaining. If there were more quantity it can be passed through a column for third time
8. Place the **Column Spin DNA** in a new **Collection tube** and add **500 μl of Dishinhibition Solution** to the reservoir. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
9. **Add 500 μl of Washing Solution** into reservoir of **Column Spin DNA**. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
10. 2º Wash. **Add 500 μl of Wash Buffer** into reservoir of **Column Spin DNA**. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
11. **Centrifuge at maximum speed for 2 minutes to remove the residual ethanol.**
12. Remove the **Collection tube** and insert the **Column Spin DNA** in a 1.5 ml microtube. **Add 50 μl of Hydration Solution** (preheated at 70°C) into reservoir of MicroSpin column. **Incubate 1 minute.**
13. **Centrifuge at maximum speed for 60 seconds.**
14. **Add another 50 μl of Hydration Solution** (preheated at 70°C) into reservoir of **Column Spin DNA**. **Incubate 1 minute.**
15. **Centrifuge at maximum speed for 60 seconds.** The microtube contains now genomic DNA

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: **Lysis/Binding Solution** and **Disinhibition Solution** contain guanidine hydrochloride which is irritating, use gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ S.L. technical service at durviz@durviz.com.

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