# REAL MICROBIOME SOIL DNA KIT

Ref. **RBMEGS18:** (para **50 EXTRACCIONES**)(for **50 EXTRATION**)

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## 1.INTRODUCTION

The DANAGENE MICROBIOME Soil DNA kit is designed for rapid and efficient purification of **microbial DNA from environmental samples such as soil samples.**

In this protocol, microorganisms are efficiently lysed by a combination of heat, chemical reactions and mechanical disruption using beads. Inhibitors are removed by precipitation using a unique wash buffer. Samples are transferred to a *Columns Spin DNA* where DNA is bound, thus requiring only one wash step prior to elution.

### Features:

- Designed for simple and fast purification of microbial DNA from different types of soil samples.
- Optimized lysis procedure using a combination of heat, chemistry and *Beads*, allows purification of DNA from yeast, fungi, Gram-negative or Gram-positive bacteria.
- Removes inhibitory substances such as humic acid derivatives and other inhibitors.
- Avoids the use of phenol/chloroforms used in precipitation extractions.

### Applications:

- Microbiome analysis
- PCR Applications
- RFLP analysis
- Pathogen identification
- Mutation analysis.

## 2.COMPONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase	T°
Washing solution	EP08	10 ml	RT
Proteinase K	REA02	30 mg	-20°C
Bead Microtubes	RBD	50 unid.	RT
Hydration Solution	REA05	10 ml	RT
EC Solution	RES05	8 ml	RT
Soil lysis solution	RES06	30 ml	RT
Enhancer solution	RES07	5 ml	RT
Soil binding solution	RES08	50 ml	RT
Columns Spin DNA	RSC03	50 unid.	RT
Collection Tubes	R30	100 unid.	RT

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\* *These solutions require pretreatment as indicated in the Preliminary Preparations section of this protocol.*

**CAUTIONS:** Soil lysis solution and Disinhibition buffer contain guanidine salts and should be handled with care. Guanidine salts form highly reactive compounds if combined with chlorine

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### Equipment and reagents needed and not provided

- \* Ethanol 100%
- \* Microcentrifuge
- \* 1.5 ml and 2.0 ml microtubes
- \* Water bath, dry bath or heated block (65°C).
- \* Vortex for homogenization of the bead microtubes, we recommend the use of the Vortex Genie 2 or similar with a horizontal rack.
- \* Bead mill homogenizer (optional).

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## 3. GENERAL PROCEDURE

### 3.1 Preliminary preparations

- Dissolve *Proteinase K* in 1.3 ml of nuclease-free water and store at -20°C. It is recommended to make several aliquots to avoid many thawing/freezing cycles. At this temperature it is stable for 1 year.
- Add 40 ml of 100% Ethanol to the *Washing Solution*. Keep the container closed to avoid evaporation of the ethanol.
- Preheat the *Hydration Solution* to 70°C.

### 3.2 General considerations

#### Amount of starting material

Generally, a reduction of starting material also helps to improve lysis efficiency and increase DNA purity.

Very dry material can absorb large volumes of *Lysis Solution* and *Enhancer Solution*. In this case, reduce the amount of sample material or add more *Lysis Solution* and *Enhancer Solution*.

If possible, remove foreign material such as leaves, stones or twigs (e.g. by sieving) as well as excess water (e.g. by discarding the supernatant after centrifuging sediment samples).

#### Sample lysing

Processing time will vary depending on sample input and disruptor speed. Times can be as little as 5 minutes when using high speed cell disrupters (*FastPrep*) or 10-20 minutes when using lower speeds (*Disruptor Genie*).

The procedure is optimized by using "Beads" in a horizontally agitated vortex (*Vortex Genie 2* or similar). Make sure that the vortex adapter allows horizontal agitation; this is a cost-effective method for the recovery of high quality microbial DNA.

Adapters with vertical tube orientation do not shake properly.

You can use "Bead mill" homogenizers such as *FastPrep*, *Precellys* and others, but follow the manufacturer's instructions to optimize sample lysis. IMPORTANT: Many modern disruption devices can cause very high energy input into the bead tubes. Depending on the type of bead tube and its contents (beads, liquid volume, sample type), a particularly high shaking frequency and/or prolonged shaking duration may cause the bead tubes to rupture. It is the responsibility of the user to perform the initial stability test of the bead tubes used under the conditions employed. Perform the initial test with water instead of *Lysis Solution* and a moderate machine setting (low frequency, short duration) to avoid spillage of chaotropic lysis buffer in case of tube breakage.

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### Sample by hand

Collect samples according to your laboratory guidelines and experimental needs.

Ensure that the samples are thoroughly mixed with the *Lysis Solution* and *Enhancer Solution* to create a homogeneous sample.

One way to ensure thorough mixing is to shake the tube with the cap down.

### 3.3 Protocol for an extraction of 250 mg of DNA from environmental samples

1. Weigh **250 mg of soil sample** into a 2.0 mL *Beads Microtubes* and add **600 µL of soil lysis solution**. Vortex.  
Depending on the nature of the soil it may be necessary to increase the *Soil Lysis Solution* and *Enhancer Solution* proportionally with the goal of recovering 400 µL of supernatant at point 5. Additional buffers may be requested.
2. Add **100 µL of Enhancer Solution**. Mix well by gently shaking the microtube.
3. Add **25µ of Proteinase K**. **Incubate at 65°C for 10 minutes**.
4. **Homogenize** together with the *Beads Microtubes* for **10 minutes** at maximum speed in the **Vortex Genie 2 or similar using a horizontal adapter**.  
The lysis time should be as short as necessary to avoid DNA shearing and minimize the release of humic acids. However, depending on the sample, it may be advantageous to increase the lysis time to 10, 20 or 30 min. **Centrifuge at 14,000 rpm for 5 minutes**.
5. Transfer up to **400 µL of the supernatant** to a clean microcentrifuge tube.  
**IMPORTANT:** There may be a layer of debris on top of the bead pellet. Avoid transferring this debris with the supernatant.
6. **Add 150 µl of EC Solution**. **Vortex**. Incubate at 0-4°C for 5 minutes.
7. **Centrifuge at 14,000 rpm for 3 minutes**. Transfer up to **500 µl of supernatant** into a new 1.5 ml microtube, avoiding touching the pellet.
8. Add **900 µl of Soil Binding Solution** and vortex briefly.
9. Load **700 µl of sample** mix into the reservoir of a combined *Collection tubes* and *Columns Spin DNA* assembly. **Centrifuge at 8,000 rpm for 30 seconds**. Remove the collection tube.
10. Discard the flow-through and repeat step 9 with the remainder of the sample mixture.  
Ensure that all sample mixture has passed into the *Collection tubes* by inspecting the *Columns Spin DNA*. If sample remains in the column, centrifuge again at 14,000 rpm for 1 minute.
11. Add **700 µl Washing Solution**. **Centrifuge at 14,000 rpm for 1 minute**. Discard the flow.
12. **Dry the silica membrane**. Centrifuge at 14,000 rpm for 3 minutes.
13. Place the *Microbial Columns Spin DNA* in a 1.5 mL nuclease-free tube (not supplied) and add **50-100 µL of Hydration solution pre-warmed to 70°C**. Incubate at room temperature for **2 minutes**.
14. **Centrifuge** the column-spin tube assembly at **14,000 rpm for 1 minute** and discard the column. The purified DNA is in the tube.

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