REAL MICROBIOME VAGINAL DNA KIT

Ref. **RBMEGS20:** (para **50 EXTRACCIONES**)

(for **50 EXTRATION**)

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1. INTRODUCTION

The **vaginal microbiome** is a specific part of the human microbiome. The unique conditions of the vagina give rise to the growth of diverse species of microorganisms, usually lactobacilli.

The vaginal microbiome harbors diverse communities of microorganisms, which form the vaginal flora, which has an important impact on the health of both the woman and her future children.

The cervicovaginal ecosystem is made up of diverse microorganisms that coexist in a dynamic balance establishing complex connections both between the different colonies and with the host. In women of reproductive age, the vaginal microbiome is generally dominated by the presence of *Lactobacillus genus*, and many women also have a prevalence of *L. crispatus*, *L. iners*, *L. jensenii* and *L. gasseri*¹. Lactobacilli contribute to the maintenance of vaginal homeostasis and prevent the colonization and growth of other undesirable microorganisms, such as those that cause sexually transmitted infections (STI).

The composition of the vaginal microbiota depends on age, stage of the menstrual cycle, hormonal fluctuations, sexual behaviors, or even the use of various drugs such as probiotics or antibiotics, which may cause an imbalance in the microbiota. As part of the Human Microbiome Project, the study of the vaginal microbiome has shown the relationship between the bacteria present in the vagina and various diseases.

REAL has developed a complete system for the study of the **VAGINAL MICROBIOME**:

1. VAGINAL SELF-COLLECTION SWAB for the collection and stabilization of microbial DNA from the vagina for a microbiome analysis.

2. REAL MICROBIOME VAGINAL DNA Kit has been designed for rapid and efficient purification of microbial DNA for vaginal samples.

Features:

VAGINAL SELF-COLLECTION SWAB

The new solution for at-home sampling, easy to perform, transport and process in the laboratory.

The vaginal self-collection swab is a novel collection method that offers several advantages compared to other traditional methods:


- Minimizes the time required
- Absolute privacy
- Maximum convenience
- Reduces anxiety

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REAL MICROBIOME VAGINAL DNA Kit

- Designed for rapid purification of microbial DNA for microbiome analysis. - Silica membrane based technology with *Columns Spin DNA*
- For efficient lysis the kit includes Bead Microtubes in combination with Proteinase K liquid.
- Sampling material: vaginal swabs.
- Typical yield: depends on the patient.
- Preparation time: 35 min.
- Elution volume: 100 µl.

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
Starting with a Bead-beating protocol, cell lysis is produced by a combination of mechanical forces, heat and detergents. After agitation using a horizontal adapter for the Vortex Genie 2 Vortex (recommended) or other rupture devices, efficient cell lysis is achieved.

Due to the addition to the lysate of a series of chaotropic salts contained in the *Binding Solution*, proper DNA binding to the column matrix is achieved. Contaminants are removed by two efficient washing steps. The resulting DNA is then recovered with a DNA-free Tris buffer for further use.

Microbial composition of room temperature preserved vaginal samples from 4 different women obtained using the VAGINAL SELF-COLLECTION SWAB. Vaginal samples were collected using our system and stored at room temperature. They were processed with the REAL MICROBIOME VAGINAL DNA Kit. The extracted DNA was subjected to microbial composition profiling by targeted sequencing of the 16s rRNA gene. Fig.1

Microbial Composition(Species)

2.COMONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envases.	T°
Washing solution	EP08	10 ml	RT
Bead Microtubes	RBD1	50 unid.	RT
Proteinase K	REA02	30 mg	-20°C
Disinhibition Solution (GHCL 8M)	REA03	18 ml	RT
Hydration Solution	REA05	10 ml	RT
CTAB Extraction Solution	RES03	50 ml	RT
Bonding Solution	RES04	15 ml	RT
Columns Spin DNA	RSC03	50 unid.	RT
Collection Tubes	R30	100 unid.	RT
Swabs	REGS20	50 unid.	RT

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(*) These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

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Equipment and reagents needed and not provided

- * Ethanol 100%
- * Microcentrifuge
- * 1.5 ml and 2.0 ml microtubes
- * Water bath, dry bath or heated block (70°C).
- * Vortex for homogenization of the bead microtubes, we recommend the use of the Vortex Genie 2 or similar with a horizontal rack.
- * Bead mill homogenizer (optional).

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3.1 Preliminary preparations

- Dissolve the *Proteinase K* in **1.3 ml of nuclease-free water** and store at -20°C. It is recommended to do several aliquots to avoid many thaw/freeze cycles. At this temperature it is stable for 1 year.
- **Add 10 ml of Ethanol 100 %** to the *Disinhibition Solution (GHCL 8M)*. Keep the container closed to avoid the ethanol evaporation.
- **Add 40 ml of Ethanol 100 %** to the *Washing Solution*. Keep the container closed to avoid the ethanol evaporation.
- Pre-heat the *Hydration Solution* at 70°C.

3.2 General considerations

- INSTRUCTIONS FOR SAMPLE COLLECTION:

Open the *Swabs* only when it is to be used taking care not to touch surfaces with the tip of the *Swabs*. Do not use if pregnant or within three months of delivery. Wash your hands. Get into a comfortable position and collect the sample. The collected specimen can be stored in the tube at room temperature. Return the *Swabs* according to the instructions provided separately. NOTES: If unexpected events not described in the instructions occur during the self-test, contact the physician.

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- Vaginal swabs should be taken with. The *Swabs* have the ability to **stabilize nucleic acids at room temperature for up to 4 weeks**.
- The procedure is optimized by using Beads that should be agitated with a Vortex that allows horizontal agitation (Vortex Genie 2 or similar). It is important that this agitation is horizontal.
- Adapters with a vertical orientation may not agitate the sample adequately.
- Ball homogenizers such as FasPrep, Precellys and others can be used, but follow the manufacturer's instructions to optimize sample lysing. **IMPORTANT:** Many modern disruption devices can apply too high an energy input to Bead tubes. Depending on the type of tube and its contents (Beads, liquid volume, sample type, etc.), especially with high vibration frequencies or excessive shaking duration, the tubes may rupture **It is the user's responsibility to test the stability of the tubes under the conditions to be used.** An initial test should be performed by replacing the lysis buffer with water and adjusting the machine to operate under moderate conditions (low frequency, short time) in order to avoid spilling the lysis buffer in case of tube breakage.
- Various modifications can be applied in addition to the standard method to increase yield, concentration and at the user's convenience.
Convenient elution (standard elution): For convenience, elution can be carried out by the single addition of 100 µL of elution buffer to the column.
High yield: Two serial elutions of 100 µL each for a total elution volume of 200 µL.
High concentration: initially use 100 µL of eluate for a second elution - total eluate volume: 100 µL, 2 elutions.

3.3 Protocol for isolation of microbial DNA from vaginal samples

- Add 900 µl of **CTAB Extraction Solution** to a 2.0 ml microtube.

Cut off the *Swabs* head together with part of the handle and insert into the microtube so that the microtube can be tightly closed. Vortex vigorously to release the cells from the swab fibers.

- Add 25 µl of Proteinase K. Incubate at 70°C for 10 minutes.**
- Remove the *Swabs* fragment from the lysis solution, pressing with the walls to recover as much liquid as possible. **Transfer 900 µl of the lysate to a microtube with Beads.**
- Homogenize** for 10 minutes at maximum speed in a Vortex Genie 2 or similar using a **horizontal adapter**.
- Centrifuge at 14,000 rpm for 5 minutes.**
- Transfer 500 µL of the supernatant to a clean microcentrifuge tube.

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7. Add **250 µL** of **Binding Solution** and vortex briefly.
8. Load the sample-containing mixture onto a **previously assembled Columns Spin DNA** with a **Collection Tube** and discard the obtained flow-through. Inspect the **Columns Spin DNA** to ensure that all the sample has passed into the **Collection Tube**. If any sample remains in the **Columns Spin DNA**, centrifuge at 14,000rpm for 1 minute.
9. Assemble the **Columns Spin DNA** to a new clean **Collection Tube** and add **500 µl** of **Disinhibition Solution (GHCL 8M)**. Centrifuge at **14,000 rpm for 1 minute**. Discard the obtained flow.
10. Add **700 µl** of **Washing Solution**. Centrifuge at **14,000 rpm for 1 minute**.
11. **Dry the silica membrane**. Centrifuge at 14,000 rpm for 3 minutes.
12. Place the **Columns Spin DNA** in a 1.5 ml nuclease-free tube (not included) and add **100 µL** of **Hydration Solution** preheated to 70°C. Incubate at room temperature for 2 minutes.
13. **Centrifuge** the assembled **Columns Spin DNA** with the tube at **14,000 rpm for 1 minute**. The purified DNA is in the microtube.

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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