# REAL SPIN MINIPREP KIT

Ref. RBMEPS06: (para 250 EXTRACCIONES)  
Ref. RBMEPS07: (para 1000 EXTRACCIONES)

(for 250 EXTRATION)  
(for 1000 EXTRATION)

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
## 1.INTRODUCTION

This kit is designed for rapid extraction and purification of **plasmid DNA from *E.coli* cultures by means of MiniSpin columns.**

The method described employs a modified alkaline lysis and **columns with silica fiber membranes** that selectively bind plasmid DNA in the presence of chaotropic salts. Plasmid DNA is purified with several washing steps to remove impurities and finally eluted.

Introduces the **TrueBLUE Lysis control reagent**, a color indicator that provides visual identification of optimal buffer mixing. It prevents errors that lead to efficient cell lysis or incomplete precipitation of SDS, genomic DNA or cell debris. This is ideal for researchers who do not have much experience with plasmid preparation and for experienced researchers who want to ensure maximum throughput.

**Visualizing efficient cell lysis and SDS precipitation using the TrueBLUE Lysis control reagent**

### Features:

- **Purifies plasmid DNA in as little as 15 minutes.**
- **Convenient:** MiniSpin columns with silica membrane.
- **Plasmid size:** 1-15 kb.
- **High Yield:** Up to 20 µg.
- **Sample volume:** 1.5-3.0 ml of bacterial cultures.

## 2.COMONENTES DEL KIT / KIT COMPONENTS

Componentes REAL	Ref.	Envase (250 PREPS)	Envase (1000 PREPS)	Tª
Resuspension solution	E25	65 ml	260 ml	RT
Lysis solution	EP02	65 ml	260 ml	RT
Washing solution	EP08	2x35 ml	4x35 ml	RT
Neutralization/binding solution	EPS03	90 ml	360 ml	RT
True Blue lysis control	RBLR	750 µl	750 µl	RT
RNase	REA04	10 mg	40 mg	-20°C
Elution solution	REA05	50 ml	200 ml	RT
Columns Spin DNA	RSC04	250 unid.	1000 unid.	RT
CollectionTubes	R30	500 unid.	1000 unid.	RT

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- \* Etanol 100 %.
- \* Microcentrifuga.
- \* Microtubos de 1.5 ml.

### Equipment and reagents needed and not provided

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### 3.1 Preliminary considerations

- Several factors can influence the obtaining of plasmid DNA. These include vector copy number, DNA insert, host strain, growth conditions and medium.
- The use of host strains containing a mutation in the Endonuclease I (end A) gene are recommended, such as JM109, DH5alpha, DH10B, XL1-Blue. DNA extraction with strains containing the endonuclease I gene product, such as HB101 and MC106, can produce samples with traces of endonucleases, therefore these strains should be avoided. If it is necessary to use these strains, please contact DanaGen-BioTed technical support who will provide you with an additional protocol.

### 3.2 Preliminary Preparations

- **CAUTION:** The **Neutralization/Binding Solution** contains guanidine hydrochloride which is an irritant agent, therefore the use of gloves and safety goggles is recommended for handling.
- The **Resuspension Solution** already contains **RNase**. The kit includes an extra vial of **RNase** for those users who need a larger amount of this enzyme to obtain a plasmid free of RNA contamination. Thus, if necessary, dissolve the **RNase** powder in the **Resuspension Solution**, then store at **-20°C**. This solution thus prepared is stable for 6 months.
- Verify that the **Lysis Solution** does not have precipitated SDS due to low temperatures. If necessary, dissolve the SDS by heating to 37°C.

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- **Add 140 ml of 100% Ethanol to each of the bottles of Washing Solution.** Keep the bottle tightly closed to prevent evaporation of the ethanol.
- Elution can be carried out at room temperature, although the yield can be increased by 20 % if the **Elution Solution is heated to 55°C.**

### **3.3 Protocol for extraction of plasmid DNA from 1-3 ml cultures**

1. Centrifuge **1.5 ml** of an overnight culture of E.coli at > 12,000 x g for 60 seconds. Remove the supernatant by aspiration without touching the pellet. 1.5 ml of the culture can be collected in the same microtube for repeat collection.

**NOTE:** The use of bacterial cultures with OD600= 2-6 units is recommended. Do not use overgrown bacterial cultures (< 16 hours incubation at 37°C with shaking). Use only fresh cultures.

2. **Resuspend the pellet in 250 µl Resuspension Solution + 2.50 µl TRUEBLUE Lysis control reagent by vortexing or pipetting,** ensuring complete resuspension of the cells.

**NOTE:** Prepare sufficient EP01/ TRUEBLUE Lysis control reagent **Resuspension Solution** for the number of minipreps to be performed. A precipitate will form after addition of TRUEBLUE lysis control reagent which will not affect. Shake to dissolve.

3. **Lyse the cells with 250 µl of Lysis Solution.** Gently invert the tube 8- 10 times until the mixture appears clear and viscous. **Do not vortex.** Incubate for 3 minutes, but never more than 5 minutes.
4. **Neutralize by adding 350µl of Neutralization / Binding Solution.** Gently invert the tube 8-10 times. Incubation for 5 minutes on ice is recommended.
5. **Centrifuge at maximum speed in a microcentrifuge for 5 minutes.** If there are particles floating in the supernatant, centrifuge again.
6. **Carefully transfer the supernatant to Columns Spin DNA with its collection tube.**
7. **Centrifuge at maximum speed for 30-60 seconds.** Remove the Spin microcolumn from the **collection tube** and place in a new **collection tube**.
8. **Add 600 µl of Washing Solution to the Columns Spin DNA reservoir. Centrifuge at maximum speed for 60 seconds.** Remove the liquid.
9. **Perform a 2nd wash. Add 600 µl of Washing Solution to the Columns Spin DNA reservoir. Centrifuge at maximum speed for 60 seconds.** Eliminate the liquid.
10. **Centrifuge at maximum speed for 3 minutes to remove residual ethanol.**
11. Insert the **Columns Spin DNA** into a new 1.5 ml microtube. **Add 50- 100 µl of Elution Solution into the spin microcolumn reservoir.**
12. **Incubate for 1 minute.** It is recommended that the buffer be at 55°C.
13. Centrifuge at maximum speed for 60 seconds. The microtube now contains the plasmid DNA.

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## STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

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### 1. Low or negative DNA yield:

- 1.1. Due to insufficient number of cells:** possibly the culture is old. Prepare a new culture. Confirm cell density.
- 1.2. Due to poor plasmid replication:** Confirm that the cells were grown in appropriate medium under optimal conditions.
- 1.3. Due to insufficient antibiotic activity:** Many antibiotics are light sensitive and degrade in prolonged storage at 2-8°C.
- 1.4. Due to prolonged alkaline lysis:** Reduce the lysis time.
- 1.5. Due to incomplete precipitation of cell debris:** Reduce initial culture volume.
- 1.6. Due to incomplete lysis:** Reduce the initial culture volume or Increase the lysis time.

### 2. Contaminated or poor quality DNA:

- 2.1. Genomic DNA:** Mixing of bacterial lysates has been performed too vigorously. Do not use vortexing. Do not use cultures that have been grown for more than 24 h.
- 2.2. Degraded and genomic DNA:** Make sure that no particles pass from the lysate to the supernatant.
- 2.3. RNA:** Digestion with RNAsa was insufficient. Confirm that RNase has been added to the EP01 solution. If the EP01 solution is older than 6 months, it is recommended to add more RNAsa.

If the EP01 solution is older than 6 months, it is recommended to add more RNAsa.

- 2.4. Nucleases:** Confirm that the solutions are free of nucleases.

### 3. Low A260/280 of purified DNA:

- 3.1. Purification is incomplete due to a large amount of DNA:** Reduce the initial volume of the cell culture.

### 4. Difficulty to redissolve plasmid DNA:

- 4.1. pellet was over dried:** air dry and heat TE or sterile water to 60-70°C.
- 4.2. Low pH:** Verify that the pH of the TE is >8.0.
- 4.3. Resuspension volume too low.**

### 6. Additional plasmid bands.

### 7.

- 5.1. Denatured plasmid DNA appears as a band below the supercoiled DNA:** Do not allow the lysis reaction to exceed 3 min.

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We recommend contacting the technical service of **DURVIZ S.L** laboratory for any additional questions regarding work protocols or problems that may arise during the work.

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