REAL DNA/RNA PURIFICATION KIT

Ref. RBMER16: (para 50 EXTRACCIONES)

(for 50 EXTRATION)

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1.INTRODUCTION

1.1 Product Description

REAL DNA/RNA PURIFICATION Kit provides a rapid method for the simultaneous extraction and purification of genomic DNA and total RNA from a single sample of cultured cells and small samples of tissues, blood, bacteria, yeast, fungi and plants.

The process first involves lysis of the cells or tissue of interest in a guanidine isothiocyanate buffer that rapidly inactivates RNases and DNases to ensure intact DNA and RNA extraction. The lysate is passed through a DNA Spin column, this column allows selective binding of genomic DNA. The column is washed and the DNA is subsequently eluted.

Ethanol is added to the liquid that has passed through the DNA Spin column to create favorable conditions for RNA to bind to the RNA Spin column, and the sample is passed through the RNA Spin column, where the total RNA binds to the membrane and contaminants are removed by washing and then the total RNA is eluted with nuclease-free water. The kit allows purification of all RNA species, from ribosomal RNA and mRNA to microRNA (miRNA) and interference RNA (RNAi).

Purified RNA can be used in a variety of applications including real time PCR, RT-PCR, Northern blotting, RNase protection and primer extension and expression array assays. The genomic DNA is of high quality and can be used in PCR reactions, sequencing, Southern blotting and SNP analysis.

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(*) These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

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Equipment and reagents needed and not provided

- * Ethanol 100%
- * Ethanol 70%.
- * Microcentrifuge
- * 1.5 ml and 2.0 ml microtubes
- * General homogenizers for homogenizing animal and plant tissues.
- * Liquid nitrogen.
- * β-mercaptoethanol.

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3.1 Preliminary preparations

- Add 100 % Ethanol to the RNA Wash Buffer indicated on the label, about **50 ml** so that the final volume will be 72 ml. Keep the container tightly closed to prevent evaporation of the ethanol.
- Add **10 ml of 100% Ethanol** to the Disinhibition Buffer. Keep the container tightly closed to avoid evaporation of the ethanol.
- Add **40 ml of 100 % Ethanol** to the DNA Wash Buffer. Keep the container tightly closed to avoid evaporation of ethanol.
- Preheat the DNA Elution Buffer to 70°C.
- Ensure that all solutions are at room temperature.
- **OPTIONAL:** The use of β-mercaptoethanol in lysis is highly recommended for animal tissues, particularly those known to have a high RNase content (e.g. pancreas), as well as plant tissues. It is also recommended for users who want to isolate RNA for very sensitive applications or microRNA enrichment. Prepare an appropriate amount of RNA Lysis Buffer by adding 10 µl of β-mercaptoethanol per 1 ml of RNA Lysis Buffer. Alternatively, RNA Lysis Solution can be used as supplied.
- **DO NOT use more sample than indicated on the DNA Spin column, as this will lead to copurify DNA with RNA. Also do not overload the RNA Spin column as this will reduce the yield and purity of the total RNA obtained. Some tissues or samples may contain high amounts of DNA which will overload the DNA Spin column, for these samples, it is recommended to perform a DNAase digestion if the total RNA eluate is to be used later in applications sensitive to the presence of small amounts of DNA.**

3.2 Protocol for DNA/RNA extraction from cells in culture

Note: The maximum recommended number of cells is 1 x 10⁶. Cell pellets can be stored at -70°C for later use (determine the number of cells before freezing). Frozen pellets should not be stored for more than 2 weeks to ensure that the integrity of the RNA will not be compromised. RNA Lysis Buffer should be added directly onto the frozen pellets without allowing them to reach room temperature.

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A) Total RNA from cells growing as a monolayer

1. Aspirate the medium and wash the cells with an appropriate amount of PBS. Aspirate the PBS.
2. Directly **add 400µl of RNA Lysis Buffer** directly to the culture plate.
3. Lyse the cells, cover the plate and shake the buffer around the surface of the plate for 5 minutes.
4. Transfer the lysate to a microtube. **IMPORTANT: Note the amount of lysate.**
5. Add a volume of 100% ethanol that is equivalent to **(15 µl of 100% ethanol is added per 100 µl of lysate)**. Vortex mix for 10 seconds.
6. Assemble a **DNA Spin column** with a collection tube provided.
7. Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
8. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
9. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
10. **Add 450 µl of 95-100% Ethanol to the retained liquid containing total RNA.** Vortex mix for 10 seconds. After the addition of ethanol, a precipitate will be visible, make sure that all of the precipitate is added to the RNA Spin column.
11. Assemble an **RNA Spin column** with a collection tube provided.
12. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute.** Discard the liquid and reassemble the column with its collection tube.
13. Repeat 12 until complete capture of all total RNA.
14. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all wash buffer has passed through the column, otherwise centrifuge for an additional minute.
15. Remove the liquid and assemble the column with the collection tube. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
16. Perform a third wash. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
17. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
18. Add **50 -100 µl of RNA Elution Buffer** to the column.
19. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if the 50 µL have not been eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
20. Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

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Extraction of genomic DNA from cells growing in a monolayer.

1. Add **500 µl of Deinhibition Buffer to the DNA Spin column. Centrifuge for 1 minute.** Remove the liquid.
2. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove liquid.
3. 2nd Wash. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove the liquid.
4. **Centrifuge at maximum speed for 2 minutes to remove all residual ethanol.**
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. **Add 50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. **Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

B) Extraction of total RNA from cells growing in suspension

1. Transfer the cell suspension to an RNase-free tube (not supplied) and centrifuge at no more than **200x g (approx. 2,000 rpm)** for 10 min to pellet the cells.
2. Carefully remove the supernatant.
3. Add **400 µl of RNA Lysis Buffer** to the pellet. Lyse the cells by micropipetting up-down or vortexing. Ensure that the entire pellet is dissolved before proceeding to the next step.
4. Transfer the lysate to a microtube. **IMPORTANT: Note the amount of lysate.**
5. Add a volume of 100% ethanol that is equivalent to **(15 µl of 100% ethanol is added per 100 µl of lysate)**. Vortex mix for 10 seconds.
6. Assemble a **DNA Spin column** with a collection tube provided.
7. Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
8. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
9. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
10. **Add 450 µl of 95-100% Ethanol to the retained liquid containing total RNA.** Vortex mix for 10 seconds. After the addition of ethanol, a precipitate will be visible, make sure that all of the precipitate is added to the RNA Spin column.
11. Assemble an **RNA Spin column** with a collection tube provided.
12. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute.** Discard the liquid and reassemble the column with its collection tube.
13. Repeat 12 until complete capture of all total RNA.
14. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all the wash buffer has passed through the column, otherwise centrifuge for an additional minute.

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15. Remove the liquid and assemble the column with the collection tube. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
16. Perform a third wash. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
17. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
18. Add **50 -100 µl of RNA Elution Buffer** to the column.
19. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm), followed by 1 minute at 14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µl are eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
20. Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

Extraction of genomic DNA from cells growing in solution.

1. Add **500 µl of Deinhibition Buffer to the DNA Spin column. Centrifuge for 1 minute.** Remove the liquid.
2. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute. Remove liquid.**
3. 2nd Wash. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove the liquid.
4. **Centrifuge at maximum speed for 2 minutes to remove all residual ethanol.**
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. **Add 50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. **Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

3.3 Protocol for extraction of DNA/RNA from animal tissues

Note: RNA in animal tissues is not protected until after collection it is broken and homogenized in Lysis Solution. Fresh or frozen tissues can be used for this procedure; tissues should be frozen in liquid nitrogen and transferred immediately to -70°C. When isolating microRNA from frozen tissues ensure that the temperature does not drop too low during weighing and spraying with a pestle and mortar.

1. Determine the amount of tissue to be weighed, it is recommended not to use more than 10 mg of tissue.
2. Place the tissue in a pestle containing an appropriate amount of liquid nitrogen to cover the sample. Pulverize the tissue using a pestle and mortar. Allow the liquid nitrogen to evaporate.
3. Transfer the sample to an RNase-free microtube (not supplied) and add **400 µl of RNA Lysis Buffer**. Homogenize by passing the lysate 5 to 10 times through a syringe with a 25-gauge needle or with a hand-held electric homogenizer.

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4. Centrifuge for 2 minutes to pellet any cell debris that may remain. Transfer the supernatant to a new RNase-free microtube (not supplied). **IMPORTANT: Note the amount of supernatant/lysate.**
5. Add a volume of 100% ethanol that is equivalent to **(15 µl of 100% ethanol is added per 100 µl of lysate)**. Vortex mix for 10 seconds.
6. Assemble a **DNA Spin column** with a collection tube provided.
7. Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
8. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
9. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
10. Add **450 µl of 95-100% Ethanol to the retained liquid containing total RNA.** Vortex mix for 10 seconds. After addition of ethanol a precipitate will be visible, make sure that all precipitate is added to the RNA Spin column.
11. Assemble an **RNA Spin column** with a collection tube provided.
12. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute.** Discard the liquid and reassemble the column with its collection tube.
13. Repeat step 12 until capture of all total RNA is complete.
14. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all of the wash buffer has passed through the column, otherwise centrifuge for an additional minute.
15. Remove the liquid and assemble the column with the collection tube. Add 400 µl of RNA Wash Buffer to the column and centrifuge for 1 minute.
16. Perform a third wash. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
17. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
18. Add **50 -100 µl of RNA Elution Buffer** to the column.
19. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µL are eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
20. The purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

Extraction of genomic DNA from animal tissues

1. Add **500 µl of Deinhibition Buffer to the DNA Spin column. Centrifuge for 1 minute.** Remove the liquid.
2. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove liquid.
3. 2nd Wash. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove the liquid.

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4. **Centrifuge at maximum speed for 2 minutes to eliminate all the residual ethanol.**
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. **Add 50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. **Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

3.4 Protocol for DNA/RNA extraction from blood.

Note: It is recommended that no more than 100 µl of whole blood be used and that EDTA has been used as an anticoagulant.

1. Transfer 100 µl of unclotted blood to an RNase-free microtube (not supplied).
2. Add **300 µl of RNA Lysis Buffer** to the blood. Vortex the cells for 15 seconds. Ensure that the mixture is clear before proceeding to the next step. **IMPORTANT: Record the amount of lysate.**
3. Add a volume of 100% ethanol that is equivalent to **(15 µl of 100% ethanol is added per 100 µl of lysate)**. Vortex mix for 10 seconds.
4. Assemble a **DNA Spin column** with a collection tube provided.
5. Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
6. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
7. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
8. Add **450 µl of 95-100% Ethanol to the retained liquid containing total RNA.** Vortex mix for 10 seconds. After the addition of ethanol, a precipitate will be visible, make sure that all of the precipitate is added to the RNA Spin column.
9. Assemble an **RNA Spin column** with a collection tube provided.
10. Apply half of the liquid + ethanol to the RNA Spin **column and centrifuge for 1 minute.** Discard the liquid and reassemble the column with its collection tube.
11. Repeat step 12 until capture of all total RNA is complete.
12. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all the wash buffer has passed through the column, otherwise centrifuge for an additional minute.
13. Remove the liquid and assemble the column with the collection tube. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
14. Perform a third wash. Add **400 µl RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
15. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
16. Add **50 -100 µl of RNA Elution Buffer** to the column.

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- Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µl are eluted, centrifuge for 1 additional minute at 14,000 x g (approx. 14,000 rpm).
- Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

Extraction of genomic DNA from blood

- Add **500 µl of Disinhibition Buffer to the DNA Spin column. Centrifuge for 1 minute.** Discard the liquid.
- Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove liquid.
- 2nd Wash. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute.** Remove the liquid.
- Centrifuge at full speed for 2 minutes to remove all residual ethanol.**
- Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. **Add 50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
- Centrifuge at maximum speed for 1 minute.** The microtube now contains the genomic DNA.

3.5 Protocol for DNA/RNA extraction from bacteria

Note: It is recommended not to use more than 1 x 10⁹ cells /ml, bacterial growth can be measured by spectrophotometer, as a rule of thumb, an E.coli culture of 1 x 10⁹ cells /ml has an OD 600 of 1.0.

Prepare an appropriate amount of TE containing lysozyme, this solution should be prepared with RNase-free TE and kept on ice until use. For Gram-negative cells the lysozyme concentration will be 1 mg/ml and Gram-positive cells the lysozyme concentration will be 3 mg/ml.

- Pellet the cells by centrifugation at **14,000 x g (approx. 14,000 rpm) for 1 minute.**
- Remove the supernatant and vortex resuspend the pellet in **100 µl of TE with the appropriate amount of lysozyme.** Incubate at room temperature for 10 minutes.
- Add **200 µl of RNA Lysis Buffer** and vortex for 15 sec.
- Add 50 µl of 95-100% Ethanol to the lysate. Vortex mix for 10 seconds. Ensure that the pellet is completely dissolved before proceeding to the next step.
- Assemble a **DNA Spin column** with a collection tube provided.
- Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
- Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
- Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.

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9. Add **450 µl of 95-100% Ethanol to the retained liquid containing total RNA**. Vortex mix for 10 seconds. After addition of ethanol a precipitate will be visible, make sure that all precipitate is added to the RNA Spin column.
10. Assemble an **RNA Spin column** with a collection tube provided.
11. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute**. Discard the liquid and reassemble the column with its collection tube.
12. Repeat step 11 until capture of all total RNA is complete.
13. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all of the wash buffer has passed through the column, otherwise centrifuge for an additional minute.
14. Remove the liquid and assemble the column with the collection tube. **Add 400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
15. Perform a third wash. **Add 400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
16. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
17. Add **50 -100 µl of RNA Elution Buffer** to the column.
18. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µL are eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
19. Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

Extraction of genomic DNA from bacteria

1. Add **500 µl of Deinhhibition Buffer to the DNA Spin column. Centrifuge for 1 minute**. Remove the liquid.
2. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute**. Remove liquid.
3. 2nd Wash. **Add 500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute**. Remove the liquid.
4. **Centrifuge at full speed for 2 minutes to remove all residual ethanol**.
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. **Add 50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. **Centrifuge at maximum speed for 1 minute**. The microtube now contains the genomic DNA.

3.6 Protocol for DNA/RNA extraction from yeasts

Note: It is recommended to use no more than 10⁷ yeast cells or 1ml of culture.

Prepare an appropriate amount of Resuspension Buffer containing lyticase: 50mM tris, pH 7.5, 10 mM EDTA, 1M Sorbitol, 0.10 % β-mercaptoethanol and 1 unit / µl of lyticase. This solution should be prepared with sterile, RNase-free reagents and kept on ice until used.

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1. Pellet the yeasts by centrifugation at **14,000 x g (approx. 14,000 rpm)** for **1 minute**.
2. Remove the supernatant and vortex resuspend the pellet in **100 µl of Resuspension Buffer with the appropriate amount of litycase**. Incubate at 37°C for 10 minutes.
3. Add **300 µl of RNA Lysis Buffer** to the pellet. Vortex the cells for 15 seconds. Ensure that the entire pellet is dissolved before proceeding to the next step.
4. Transfer the lysate to a microtube. **IMPORTANT: Note the amount of lysate**.
5. Add a volume of 100% ethanol that is equivalent to **(15 µl of 100% ethanol is added per 100 µl of lysate)**. Vortex mix for 10 seconds.
6. Assemble a **DNA Spin column** with a collection tube provided.
7. Add the lysate + ethanol to the column. **Centrifuge for 1 minute**.
8. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA**. Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
9. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
10. Add **450 µl of 95-100% Ethanol to the retained liquid containing total RNA**. Vortex mix for 10 seconds. After addition of the ethanol a precipitate will be visible, make sure that all of the precipitate is added to the RNA Spin column.
11. Assemble an **RNA Spin column** with a collection tube provided.
12. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute**. Discard the liquid and reassemble the column with its collection tube.
13. Repeat 12 until complete capture of all total RNA.
14. Add 400 µl of RNA Wash Buffer to the column and centrifuge for 1 minute. Make sure that all wash buffer has passed through the column, otherwise centrifuge for an additional minute.
15. Remove the liquid and assemble the column with the collection tube. **Add 400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
16. Perform a third wash. **Add 400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
17. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
18. Add **50 -100 µl of RNA Elution Buffer** to the column.
19. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µl are eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
20. Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

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Yeast Genomic DNA Extraction

1. Add 500 µl of Deinhhibition Buffer to the DNA Spin column. Centrifuge for 1 minute. Remove the liquid.
2. Add 500 µl of DNA Wash Buffer to the reservoir of the DNA Spin column. Centrifuge for 1 minute. Remove liquid.
3. 2nd Wash. Add 500 µl of DNA Wash Buffer to the reservoir of the DNA Spin column. Centrifuge for 1 minute. Remove the liquid.
4. Centrifuge at full speed for 2 minutes to remove all residual ethanol.
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. Add 50-100 µl DNA Elution Buffer (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. Centrifuge at maximum speed for 1 minute. The microtube now contains the genomic DNA.

3.7 Protocol for DNA/RNA extraction from plants

Note: It is recommended not to use more than 10 mg of tissue or 1 x 10⁶ plant cells.

Fresh or frozen tissues can be used for this procedure. Tissues should be frozen in liquid nitrogen and transferred immediately to -70°C. When isolating microRNA from frozen tissues ensure that the temperature does not drop too low during weighing and spraying with a pestle and mortar to ensure that the integrity of the microRNA is not compromised.

1. Place the tissue or cells in a mortar containing an appropriate amount of liquid nitrogen to cover the sample. Spray the tissue using a pestle and mortar. Allow the liquid nitrogen to evaporate.
2. Transfer the sample to an RNase-free microtube (not supplied) and add **400 µl of RNA Lysis Buffer.** Homogenize with a handheld electric homogenizer.
3. Centrifuge for 2 minutes to pellet any cell debris that may remain. Centrifuge for 2 minutes to pellet any cell debris that may remain. Transfer the supernatant to a new RNase-free microtube (not supplied). **IMPORTANT: Note the amount of supernatant/lysate.**
4. Add a volume of 100% ethanol that is equivalent to **15 µl of 100% ethanol is added per 100 µl of lysate.** Vortex mix for 10 seconds.
5. Assemble a **DNA Spin column** with a collection tube provided.
6. Add the lysate + ethanol to the column. **Centrifuge for 1 minute.**
7. **Retain the liquid that has passed through the column and is in the collection tube containing total RNA.** Transfer the retained liquid to a new 1.5 ml microtube (not supplied with the kit).
8. Assemble the DNA Spin column with a collection tube and store at room temperature (15-25°C) or at 4°C for later DNA purification.
9. **Add 450 µl of 95-100% Ethanol to the retained liquid containing total RNA.** Vortex mix for 10 seconds. After addition of ethanol a precipitate will be visible, make sure that all precipitate is added to the RNA Spin column.

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6. **Centrifugar a máxima velocidad durante 1 minuto**. El microtubo contiene ahora el ADN genómico.
10. Assemble an **RNA Spin column** with a collection tube provided.
11. Apply half of the liquid + ethanol to the RNA Spin column and **centrifuge for 1 minute**. Discard the liquid and reassemble the column with its collection tube.
12. Repeat step 12 until the capture of all total RNA is complete.
13. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Make sure that all the Wash Buffer has passed through the column, otherwise centrifuge for an additional minute.
14. Remove the liquid and assemble the column with the collection tube. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute.
15. Perform a third wash. Add **400 µl of RNA Wash Buffer** to the column and centrifuge for 1 minute. Remove the liquid from the collection tube and centrifuge for an additional 2 minutes to remove all residual ethanol. Discard the collection tube.
16. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute total RNA.
17. Add **50 -100 µl of RNA Elution Buffer** to the column.
18. Centrifuge for 2 minutes at **200x g (approx. 2,000 rpm)**, followed by 1 minute at **14,000 x g (approx. 14,000 rpm)**. Check the amount of eluate on the column, if not all 50 µl are eluted, centrifuge for an additional 1 minute at 14,000 x g (approx. 14,000 rpm).
19. Purified RNA can be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for long term storage.

Extraction of genomic DNA from plants

1. Add **500 µl of Deinhhibition Buffer to the DNA Spin column**. **Centrifuge for 1 minute**. Remove the liquid.
2. Add **500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute**. Remove liquid.
3. 2nd Wash. Add **500 µl of DNA Wash Buffer** to the reservoir of the DNA Spin column. **Centrifuge for 1 minute**. Remove the liquid.
4. **Centrifuge at full speed for 2 minutes to remove all residual ethanol**.
5. Remove the collection tube and insert the DNA Spin column into a 1.5 ml microtube. Add **50-100 µl DNA Elution Buffer** (pre-warmed to 70°C) to the reservoir of the DNA Spin column. Incubate for 1 minute.
6. Centrifuge at maximum speed for 1 minute. The microtube now contains the genomic DNA.

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: RNA Lysis Buffer contains guanidine isothiocyanate and Disinhibition Buffer contains guanidine hydrochloride which are irritants, wear gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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