REAL SPIN VIRAL RNA KIT

Ref. RBMER19: (para 100 EXTRACCIONES)

(for 100 EXTRACTION)

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1.INTRODUCTION

REAL Spin Viral RNA Kit is designed for rapid purification of viral RNA from cell-free samples such as serum, plasma and other biological fluids.

The sample is first lysed in the presence of chaotropic salts that inactivate RNases ensuring extraction of intact viral RNA and allow selective binding to a glass fiber membrane located on a spin column. The viral RNAs remain bound while a series of washes and centrifugation steps will remove contaminating cellular components. Finally, with an elution buffer or nuclease-free water the viral RNAs will be eluted from the membrane. This process does not require alcohol precipitation, organic solvent extractions, or extensive nucleic acid handling.

The REAL Spin Viral RNA Kit can be used for the extraction of viral RNA from a wide range of RNA viruses. **However, success cannot be guaranteed for all virus species and must be validated by the user.**

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(*) These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

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- * Etanol 100 %.
- * Microcentrifuga.
- * Microtubos de 1.5 ml y 2.0 ml.

Equipment and reagents needed and not provided

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3.1 Preliminary Preparations

- **Lysis Buffer, microRNA Binding Buffer and Wash Buffer 1** contain **thiocyanate guanidinium** which is an irritant agent, for this reason, we recommend the use of goggles and gloves for handling.
- **Add 80 ml of 100% Ethanol Virus Wash Buffer 2***. Keep the container tightly closed to avoid evaporation of the ethanol.

3.2 General recommendations

Types of samples

- Plasma/serum: Avoid working with samples where hemolysis is observed. After obtaining plasma or serum, it is important to centrifuge the sample to obtain a cell-free starting material.
- For solid samples such as tissue (5-10 mg) it can be homogenized in 300-400 µl of PBS using a hand-held electric homogenizer or bead-based homogenization systems. Centrifuge the sample and use 200 µl of the clear, particle-free supernatant.
- For stool prepare a suspension with PBS, 10% (w/v) and use the best system to lyse all viral particles. Centrifuge the sample and use 200 µl of the clear particle-free supernatant.
- Swabs: Incubate the swab in an appropriate amount of buffer (e.g. PBS) or cell-free medium for 30 minutes with motion. Remove swab and proceed with 200 µl of the clear particle-free supernatant.
- Transport Medium/Viral Transport Medium: Vortex the tubes containing the swab at maximum speed for 1 minute. Use 200 µl as the input sample.

In tissues and feces it should be considered that other RNA can be copurified which may inhibit the following PCR assays.

3.3 Protocol for the extraction of viral RNA from cell-free biological fluids

1. **Add 200 µL of sample** to a microtube. If you process < 200 µl samples adjust the final volume to 200 µl using PBS.
2. **Add 400 µl Viral Lysis Buffer**. Close the microtube and vortex vigorously for 20 seconds.
3. **Incubate at room temperature for 10-15 minutes**.

Incubation time and temperature is critical for lysis as well as for RNA stability. Incubation at room temperature is usually sufficient without significant loss of sensitivity. It is recommended to optimize these points for the type of sample to be used. Protocols with and without the use of proteinase K as well as different incubation times and temperatures can be compared.

If the sample is very viscous (sputum) the use of proteinaseK or an incubation at 70°C is recommended.

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4. Add **30 µl of Protein Precipitation Buffer**. Vortex 5 seconds. Incubate 1 minute at room temperature.
SKIP this step and go to step 5 if processing samples from Viral Transport Media such as our REALSALIVA/SWAB VIRAL Sample Collection Kit, eNAT, DNA/RNA Shield, etc.
5. Centrifuge for **2 minutes at 11,000 x g**.
6. **Transfer the supernatant**, avoiding touching the pellet that may form, to a new microtube.
7. Add **350 µl 100% ethanol**. Mix well.
8. Transfer half of the lysate to a **Spin RNA column** with its collection tube. **Centrifuge at 8,000 rpm for 30 seconds**.
9. **Transfer the other half and centrifuge at 8,000 rpm for 30 seconds**.
10. Add **100 µl Wash Buffer 1**. Centrifuge at 11,000 x g for 1 minute.
11. Add **700 µl Wash Buffer 2**. Centrifuge at 11,000 x g for 1 minute.
12. **Centrifuge 3 minutes at maximum speed** to remove all ethanol.
13. Elution with **50 µl Nuclease-free water**. 2 minutes incubation.
It is very important to add the nuclease-free water in the center of the membrane so that it is completely wetted.
14. **Centrifuge at 10,000 rpm for 60 seconds**. Collect the 50 µl and redeposit in the center of the membrane. This increases the yield.
15. **Incubate for 2 minutes and centrifuge at maximum speed**.
If problems are observed in subsequent detections, the volume of the eluate can be changed.
the volume of eluate added to the RT-PCR.

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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