# REAL TISSUE/CELLS RNA KIT

Ref. **RBMER20:** (para **100 EXTRACCIONES**)  
 Ref. **RBMER21:** (para **150 EXTRACCIONES**)

(for **100 EXTRATION**)  
 (for **150 EXTRATION**)

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## 1. INTRODUCTION

This kit is optimized to be a fast and efficient method to isolate **total RNA from tissues and cells using MicroSpin Columns** without the use of toxic reagents such as chloroform or phenol.

The REAL TISSUE / CELLS RNA kit integrates a genomic DNA removal column. This column quickly and efficiently removes most genomic DNA without the need for DNase digestion.

In the first step, cells and tissues are lysed without the use of  $\beta$ -mercaptoethanol. The chaotropic salt included in the **RNA Lysis Solution** immediately inactivates RNases. The lysate is added to the **MicroSpin Columns** to clarify the lysate and remove contaminating gDNA. After addition of the **Protein Precipitation Solution (RNA)**, the RNA binds to the **Spin RNA Columns**. Subsequently, two wash steps remove salts, metabolites and cellular components. The high-quality RNA is eluted with RNAase-free H<sub>2</sub>O.

### Features:

- Efficient and rapid isolation of total RNA.
- Lysate filtration and genomic DNA removal on the same column.
- Complete removal of contaminants and inhibitors.
- Sample size: < 1x 10<sup>7</sup> cells in culture; 25 mg animal or human tissue.
- Elution volume: 30  $\mu$ l.
- No phenol/chloroform, CICs gradients, LiCl or ethanol precipitation.

## 2. COMPONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase (100 Preps)	Envase (150 Preps)	T°
Washing solution 2	EP09	20 ml	4x20 ml	RT
Washing solution 1	ER14	10 ml	50 ml	RT
Nuclease-Free Water	ER7	8 ml	40 ml	RT
MicroSpin columns	RSC05	100 unid.	500 unid.	RT
Spin RNA Columns	RSCV02	100 unid.	500 unid.	RT
RNA lysis solution	V01	45 ml	2x100 ml	RT
Protein precipitation solution (RNA)	V02	4 ml	20 ml	RT
Collection Tubes	R30	200 unid.	1000 unid.	RT

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(\*). These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

**RNA Lysis Solution and Washing Solution 1 contain isothiazacyanate guanidinium which is a potent irritant, wear gloves and goggles. Both buffers can form hazardous reactive components when combined with bleach.**

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### Equipment and additional reagents required

- 100 % Ethanol.
- Microcentrifuge.
- Micropipettes
- RNase-free 1.5 ml or 2.0 ml microfuge tubes.
- General homogenisers for homogenising animal tissues.
- Liquid Nitrogen , mortar and pestle.
- RNAase-free environment.

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### 3.1 Preliminary Preparations

- Add **80 ml 100 % ethanol** (kit 100 extractions) to **Washing Solution 2** and about **80 ml** to each container (kit 500 extractions) indicated on the label. Keep the container tightly closed to avoid evaporation of the ethanol.

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### 3.2 General recommendations

#### Sample collection and RNase inhibition

RNA is not protected until the sample material is flash frozen or cleaved/smoothed in the presence of RNase inhibiting or denaturing agents.

Sample collection methods:

- Use a freshly collected sample for immediate lysis and RNA purification.
- Samples can be stored in **RNA Lysis Solution** after lysis at -70°C for one year, at 4°C for up to 24 hours or up to several hours at room temperature. Samples frozen in **RNA Lysis Solution** should be thawed slowly before starting RNA isolation.
- Freeze the sample in liquid Ni immediately after collection and store at 70 °C. This can be done with a mortar and pestle. This can be done with a mortar and pestle by spraying the sample in a frozen state. Ensure that the sample does not thaw before contact with the lysis buffer.
- Samples can be immersed and stored in **RNAlater**. Before using such samples, remove excess RNAlater from the tissue before use.

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#### Sample breakage/lysis and homogenisation

Efficient sample breakage/lysis and homogenisation is the most important requirement for a successful total RNA isolation process. Both terms are two distinct steps:

Breakage/Lysis: Complete breakage of the cell walls and plasma membranes of cells and organelles is required to release all RNA present in the sample. Insufficient breakage will result in low yields.

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Homogenisation: This is necessary to reduce the viscosity of the lysate produced. Incomplete homogenisation will result in inefficient binding of RNA to the membrane and therefore low yields.

Some breakage methods simultaneously homogenise the sample as well (**manual electric homogenisers**) while others require an additional homogenisation step.

### **3.3 Protocol for RNA purification from animal or human tissue**

Process samples of up to **25 mg** of fresh or frozen tissue and spray with liquid Ni.

**IMPORTANT:** It is essential for efficient RNA preparation that all RNA contained in the sample is released from the cells by homogenisation with a mechanical homogeniser (Polytron type), take care to keep the rotor submerged to avoid excessive foaming and choose a homogeniser with a 5-7 mm rotor that can be used in microtubes.

For fresh and soft tissues use the homogeniser; for fresh hard or RNase-rich tissues spray with liquid Ni; for frozen soft tissues or small pieces use the homogeniser; for all other frozen tissues spray with liquid Ni.

Purification of RNA from skeletal muscle, heart and skin can be difficult due to the abundance of contractile proteins, connective tissue and collagen. To remove these proteins that can interfere with purification, the sample should be treated with proteinase K.

1. Add **400 µl of ARN Lysis Solution** to the tissue sprayed with liquid Ni. Homogenise with a manual electric homogeniser or pass the lysate 10 times through a syringe with a 20-G needle (0.90 mm). Incubate for 3 minutes at room temperature.
2. Add **30 µl Protein Precipitation Solution (RNA)**. Vortex and incubate for 1 minute.
3. Centrifuge **3 minutes at maximum speed. Transfer the supernatant to a MicroSpin Columns.**
4. **Centrifuge for 1 minute at 8000 rpm.** Transfer the filtrate to a new centrifuge tube. This step also removes a large part of the contaminating genomic DNA, not being total for those applications that require a total elimination, for which a DNase I treatment must be carried out on the column or once the RNA has been eluted.
5. Add **350 µl of 100% ethanol to the lysate collected in step 4.**
6. Take an **Spin RNA Columns** plus its **Collection Tubes** and add the lysate. Centrifuge at **8,000-10,000 rpm** for **60 seconds.**
7. Add **100 µl Washing Solution 1.** Centrifuge at maximum speed for 1 minute.
8. Add **700 µl Washing Solution 2.** Centrifuge at maximum speed for 1 minute.
9. **Centrifuge at maximum speed for 3 minutes** to remove all ethanol.
10. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute the total RNA.

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11. Elute the RNA in **30 µl Nuclease-Free Water**. Incubate for 2 minutes and centrifuge **at maximum speed** for **1 minute**.
12. Repeat step 11 using another **30 µl of Nuclease-Free Water**. Incubate 2 minutes and centrifuge **at maximum speed** for **1 minute**. The eluate from step 11 can also be used if a high concentration of RNA is required but in this case the yield will be 25 % lower than if a second elution is performed.

### **3.4 Protocol for RNA purification from cultured cells (<math>1 \times 10^7</math>)**

#### **Adherent cultured cells in a dish**

Remove the culture medium and wash the cells with PBS. Aspirate PBS and add 0.10-0.30 % trypsin in PBS and incubate until cells detach from the plate surface. Then add culture medium and transfer the cells to an appropriate tube and centrifuge to obtain the cell pellet. Remove supernatant and proceed to step 1.

#### **Suspension cultured cells**

Transfer cultured cells (up to  $1 \times 10^7$ ) to a sterile 1.5 ml microtube. Centrifuge at 6000 x g for 1 minute. Discard the supernatant and proceed to step 1.

1. Add **400 µl of RNA Lysis Solution** and resuspend the cells with the micropipette. Incubate for 15-20 minutes at room temperature.
2. Add **30 µl Protein Precipitation Solution (RNA)**. Vortex and incubate for 1 minute.
3. Centrifuge 3 minutes at maximum speed. Transfer the supernatant to a **MicroSpin Columns**.
4. **Centrifuge 1 minute at 8000 rpm.** Transfer the filtrate to a new centrifuge tube. This step also removes a large part of the contaminating genomic DNA, not being total for those applications that require a total elimination, for this, carry out a treatment with DNase I in the column or once the RNA has been eluted.
5. Add **350 µl of 100% ethanol to the lysate collected in step 4.**
6. Take an **Spin RNA Columns** plus its collection tube and add the lysate. Centrifuge at **8,000 - 10,000 rpm for 60 seconds.**
7. Add **100 µl Washing Solution 1.** Centrifuge at maximum speed for 1 minute.
8. Add **700 µl Washing Solution 2.** Centrifuge at maximum speed for 1 minute.
9. **Centrifuge at maximum speed for 3 minutes** to remove all ethanol.
10. Place the column in a new 1.5 ml microtube (not supplied with the kit) to elute the total RNA.

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11. Elute the RNA in **30 µl Nuclease-Free Water**. Incubate for 2 minutes and centrifuge **at maximum speed** for **1 minute**.
12. Repeat step 11 using a further **30 µl of Nuclease-Free Water**. Incubate for 2 minutes and **centrifuge at maximum speed** for **1 minute**. The eluate from step 11 can also be used if a high concentration of RNA is required but in this case the yield will be 25 % lower than if a second elution is performed.

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#### APPENDIX

##### DNase digestion in solution

Passing the lysed sample through a **MicroSpin Columns** according to the standard protocol is very efficient in binding DNA, resulting in minimal residual DNA in the purified RNA. Residual DNA will not be detectable in most downstream applications. Despite this, there are still certain applications that require even lower residual DNA contents. However, the removal of DNA to a completely undetectable level is complicated and the efficiency of the Gdna removal column is sometimes not sufficient for downstream applications that require a lower residual DNA content. This can occur especially in cases where a large amount of sample or a sample containing a lot of DNA is processed. The amount of residual DNA detected depends on the type of sample, the amount and its DNA content and the detection sensitivity of the method used to analyse the residual DNA.

Digestion of DNA in solution can efficiently destroy contaminating DNA. However, strict control of RNases and subsequent repurification of RNA (in order to remove buffer, salts, DNAase and digested DNA) is usually required.

We recommend the use of our **REAL DNA Removal Kit** which provides a method to remove contaminating genomic DNA in RNA preparations using 2 sequential filtrations with different columns. It should be noted that this method reduces the amount of RNA, which is why it is important to decide whether complete removal of possible contaminating DNA is necessary for the realisation of our subsequent application.

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#### STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed

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