REAL MICRORNA KIT

Ref. RBMER28: (para 50 EXTRACCIONES)

(for 50 EXTRATION)

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1.INTRODUCTION

REAL MicroRNA Kit provides a rapid and efficient method for the isolation and purification of small RNA molecules (<200 nt) from cultured animal cells, small tissue samples, and bacterial cells. These small RNAs include regulatory RNA molecules such as **microRNA (miRNA)** and **short interfering RNA (siRNA)**.

Most commercial RNA purification kits do not recover RNA molecules smaller than < 200 nucleotides, using an approach consisting of two sequential filtrations with different ethanol concentrations, an RNA fraction highly enriched in RNA species < 200 nucleotides can be obtained with of REAL MicroRNA Kit.

Features:

- microRNA isolation from cultured cells and tissues.
- Fast protocol: 25 minutes.
- microRNA isolated without the use of harmful chemicals as phenol or chloroform.
- Efficient isolation of small RNA species using 2 column process, resulting in minimal contamination of larger RNA and gDNA.
- Elution volume: 30 µl.
- Purified RNA can be used in a number of downstream applications including: qPCR, reverse transcription PCR, Northern blotting, RNase protection and expression array assays.

2.COMPONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase. (50 Preps)	T°
RNA Lysis Solution	V01	20 ml	RT
Washing Solution 1	ER14	6 ml	RT
Washing Solution 2 *	ER6	10 ml	RT
Nuclease-Free Water	ER7	4 ml	RT
Columns RNA Large Removal	RSC02	50 unid.	RT
Spin RNA Columns	RSCV01	50 unid.	RT
Collection Tubes	R30	150 unid.	RT

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* These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

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Equipment and reagents needed and not provided

- 100% ethanol.
- Microcentrifuge
- Micropipettes
- 1.5 ml and 2.0 ml microtubes.
- General homogenisers for homogenising animal tissues.
- Liquid nitrogen.
- RNAase-free environment.
- β -mercaptoethanol.
- PBS
- Trypsin
- Buffer TE: Tris-EDTA
- Buffer resuspension: 50mM Tris, pH 7.5, 10 mM EDTA, 1M Sorbitol, 0.10 % μ -mercaptoethanol and 1 unit / μ l of lyticase

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3.1 Preliminary preparations

- Add **40ml of Ethanol 100 %** to the **Washing Solution 2**. Keep the container closed to avoid the ethanol evaporation.
- **OPTIONAL:** The use of β -mercaptoethanol in lysis is highly recommended for most animal tissues, particularly those known to have a high RNAse content (ex.pancreas), as well as for most plant tissues. It is also recommended for users who wish to isolate RNA for sensitive downstream applications. **Add 10 μ l of β -mercaptoethanol** (provided by the user) to each 1 ml of **RNA Lysis Solution** required. Alternatively, the **RNA Lysis Solution** can be used as provided.

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3.2 General important recommendations

Sample collection and RNase inhibition

The RNA is not protected until the sample material freezes instantaneously or breaks / lyses in the presence of inhibitors or denaturing agents of the RNases.

Immediate stabilization of the DNA expression pattern is a prerequisite for accurate gene-expression analysis.

Methods for sample collection

- Use a freshly collected sample for immediate lysis and RNA purification.
- Samples can be stored in the **RNA Lysis Solution** after lysis at -80°C for one year, at 4°C for up to 24 hours or up to several hours at room temperature. Frozen samples in **RNA Lysis Solution** should be thawed slowly before starting RNA isolation.
- Freeze the sample in liquid N_2 immediately after collection and store at 80°C . This can be done with a mortar and pestle by spraying the sample in a frozen state. Make sure that the sample does not thaw before contact with the **RNA Lysis Solution**.
- Samples can be submerged and stored in or **RNAlater**.

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Disruption and homogenization of starting materials for RNA isolation

Efficient disruption and homogenization of the starting material is an absolute requirement for all total RNA isolation procedures. Disruption and homogenization are two distinct steps.

Disruption: Complete disruption of cells walls and plasma membranes of cells and organelles is absolutely required to release all the RNA contained in the sample. Different samples requires different methods to achieve complete disruption. Incomplete disruption results in significantly reduced yields.

Homogenization: Is necessary to reduce the viscosity of the cells lysates produced by disruption. Homogenization shears the high-molecular-weight genomic DNA and other high-molecular-weight cellular components to create a homogeneous lysate. Incomplete homogenization results in efficient binding of RNA and therefore reduced yields.

Some disruption methods simultaneously homogenize the sample (**rotor-stator-homogenizer**) while others require an additional homogenization step

3.3 Protocol for microRNA purification from human or animal tissues

Process samples of up to **25mg** of fresh or frozen tissue. Grind the sample to a fine powder in the presence of liquid Nitrogen with a pestle and mortar.

- Transfer the tissue into a mortar that contains an appropriate amount of liquid nitrogen to cover the sample. Grind the tissue thoroughly using a pestle. Allow the liquid nitrogen to evaporate, without allowing the tissue to thaw.
- Transfer the sample to an RNase-free tube (not provided) and add 400µl de **RNA Lysis Solution.** Homogenize by passing the lysate 5-10 times through a 20-G (0.90 mm) gauge needle attached to a syringe or an electrical homogenizer. **Incubate at room temperature for 5 minutes.**
- Centrifuge for **3 minutes a maximum speed.**
- Transfer de supernatant to another RNase-free microcentrifuge tube (not provided). **IMPORTANT: Note the volume of the supernatant/lysate.**
- Add a volume of **100% ethanol** that is equivalent to **(55µl of 100 % ethanol is added to every 100µl of supernatant/lysate).** Mix by vortex for 10 seconds.
- Assemble a **Columns RNA Large Removal** with one the provided **Collection tubes.**
- Apply the clarified lysate with ethanol into the columns. **Centrifuge at 10.000 rpm for 1 minute.**
- Retain the liquid that has passed through the column and that is in the Collection tubes**, which contains the **small RNA species.**
- Add **700µl de etanol 100% to the lysate.** Mix well.
- Pass half of the lysate to a **Spin RNA Columns** with its **Collection tubes.** **Centrifuge at 10.000 rpm for 30 seconds.**
- Place the column in a new **Collection tube.** **Pass the other half and centrifuge at 10,000 rpm for 30 seconds.**
- Add **100µl of RNA Washing Solution 1.** Centrifuge at maximum speed for 1 minute.
- Add **700µl of RNA Washing Solution 2.** Centrifuge at maximum speed for 1 minute.
- Centrifuge for 3 minutes at maximum speed (11.000 rpm) to completely dry the membrane.

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15. Place the column in a 1.5 ml (not supplied) Microtube.
16. **Add 25-30µl Nuclease-Free Water** into reservoir of Spin column. Dispense buffer directly onto the silica membrane. **Incubate at room temperature for 2 min.**
It is very important to add the elution buffer in the centre of the membrane to be completely wet.
17. **Centrifuge at 10.000 rpm for 60 seconds.** Collect 25-30µl and redeposit in the centre of the membrane. This increases yield.
18. **Incubate for 2 minutes and centrifuge at maximum speed.** The microtube now contains the MicroRNAs. Use or store at -20°C or -80°C.

3.4 Protocol for microARN purification from cultured cells (< 1 x 10⁷)

Adherent Cultured Animal Cells

Remove the culture medium and wash the cells with PBS. Aspirate cell-culture medium and add an equal amount of PBS in order to wash the cells. Aspirate PBS. Add 0.1-0.3% trypsin PBS and incubate for an appropriate time to detach the cells the dish surface. After cell detachment, add cell culture medium, transfer cells to an appropriate tube and pellet by centrifugation. Remove supernatant and continue with the point 1.

Suspension Cultured Animal Cells

Transfer cells (up to 1 x 10⁷) to a 1.5 ml microcentrifuge tube. Harvest cells by centrifugation (6000 rpm), remove the supernatant and continue with the point 1.

1. Add **400µl of RNA Lysis Solution**. Lyse the cells using up-down the micropipette to dissolve the cell pellet. Make sure that all pellets are dissolved before proceeding with the next step.
2. Add **210µl ethanol 100%** to the lysate. Mix by vortex for 10 seconds.
3. Assemble a **Columns RNA Large Removal** with one the provided **Collection tubes**.
4. Apply the clarified lysate with ethanol into the columns. **Centrifuge at 10.000 rpm for 1 minute.**
5. **Retain the liquid that has passed through the column and that is in the collection tube, which contains the small RNA species.**
6. Add **700µl de etanol 100%** to the lysate. Mix well.
7. Pass half of the lysate to a **Spin RNA Columns** with its **Collection tubes**. **Centrifuge at 10.000 rpm for 30 seconds.**
8. Place the column in a new **Collection tube**. **Pass the other half and centrifuge at 10,000 rpm for 30 seconds.**
9. Add **100µl of RNA Washing Solution 1**. Centrifuge at maximum speed for 1 minute.
10. Add **700µl of RNA Washing Solution 2**. Centrifuge at maximum speed for 1 minute.
11. **Centrifuge for 3 minutes at maximum speed (11.000 rpm) to completely dry the membrane.**
12. Place the column in a 1.5 ml (not supplied) microtube. To avoid **Spin RNA columns**.

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13. **Add 25-30µl Nuclease-Free Water** into reservoir of Spin column. Dispense buffer directly onto the silica membrane. **Incubate at room temperature for 2 min.**
It is very important to add the elution buffer in the centre of the membrane to be completely wet.
14. **Centrifuge at 10.000 rpm for 60 seconds.** Collect 25-30µl and redeposit in the centre of the membrane. This increases yield.
15. **Incubate 2 minutes and centrifuge at maximum speed for 1 minute** The microtube now contains the microRNAs. Use or store at -20°C or -80°C.

3.5 Protocol for microARN purification from bacteria cultures.

Note: It is recommended that no more than 1×10^9 bacteria cells/ml be used. Bacterial growth can be measured using a spectrophotometer, as a general rule, an E.coli culture containing 1×10^9 bacteria cells has an OD₆₀₀ of 1.0.

Prepare the appropriate lysozyme-containing TE buffer, this solution should be prepared with sterile, Rnase free TE buffer, and kept on ice until needed. For Gram-negative the concentration of lysozyme will be of 1 mg/ml and for Gram-positive the concentration of lysozyme will be of 3 mg/ml.

1. Pellet bacteria by centrifugation at (**aprox. 14.000 rpm for 1 minute.**)
2. Decant supernatant, and carefully remove any remaining media by aspiration. Resuspended the bacteria in **100µl of appropriate lysozyme-containing TE by vortexing.** **Incubate at room temperature for 10 minutes.**
3. Add **300µl of RNA Lysis Solution.** Lyse cells by vortex for 15 seconds.
4. Add **210µl of ethanol 100%** to lysate. Mix by vortex for 10 seconds. Ensure that the entire pellet is completely dissolved before proceeding to the next step.
5. Assemble a **Columns RNA Large Removal** with one the provided **collection tubes.**
6. Apply the clarified lysate with ethanol into the columns. **Centrifuge at 10.000 rpm for 1 minute.**
7. **Retain the liquid that has passed through the column and that is in the collection tubes, which contains the small RNA species.**
8. Add **700µl de ethanol 100% to the lysate.** Mix well.
9. Pass half of the lysate to a **Spin RNA Columns** with its **collection tubes.** **Centrifuge at 10.000 rpm for 30 seconds.**
10. Place the column in a new **collection tube.** **Pass the other half and centrifuge at 10,000 rpm for 30 seconds.**
11. Add **100µl of RNA Washing Solution 1.** Centrifuge at maximum speed for 1 minute.
12. Add **700µl of RNA Washing Solution 2.** Centrifuge at maximum speed for 1 minute.

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13. Centrifuge for **3 minutes at maximum speed (11.000 rpm)** to completely dry the membrane. Place the column in a 1.5 ml (not supplied) microtube.
14. Remove the collection tube and insert the **Spin RNA Columns** in a 1.5 ml microtube.
15. Add **25-30µl Nuclease-free water** into reservoir of Spin column. Dispense buffer directly onto the silica membrane. **Incubate at room temperature for 2 min**.
It is very important to add the **Nuclease-free water** in the centre of the membrane to be completely wet.
16. **Centrifuge at 10.000 rpm for 60 seconds**. Collect 25-30µl and redeposit in the centre of the membrane. This increases yield.
17. **Incubate 2 minutes and centrifuge at maximum speed for 1 minute**.

3.6 Protocol for microARN purification from yeast cultures

Note: It is recommended to use no more than 107 yeast cells or 1ml of culture.

Prepare the appropriate amount of lyticase resuspension buffer have the following composition: 50mM Tris, pH 7.5, 10 mM EDTA, 1M Sorbitol, 0.10 % μ-mercaptoethanol and 1 unit /µl of lyticase. This solution should be prepared with sterile, RNase free reagents, and kept on ice until needed.

1. Pellet yeast for centrifugation (**aprox. 14.000 rpm) for 1 minute**.
2. Decant supernatant, and carefully remove any remaining media by aspiration. Resuspended the yeast in **100 µl of Lyticase resuspension buffer by vortexing**. Incubate at 37°C for 10 minutes.
3. Add **300µl of RNA Lysis Solution**. Lyse cells by vortex for 15 seconds.
4. Add **210µl of ethanol 100%** to lysate. Mix by vortex for 10 seconds. Ensure that the entire pellet is completely dissolved before proceeding to the next step.
5. Assemble a **Columns RNA Large Removal** with one the provided **Collection tube**.
6. Apply the clarified lysate with ethanol into the columns. **Centrifuge at 10.000 rpm for 1 minute**.
7. **Retain the liquid that has passed through the column and that is in the Collection tube**, which contains the **small RNA species**.
8. Add **700 µl de ethanol 100%** to the lysate. Mix well.
9. Pass half of the lysate to a **Spin RNA Columns** with its **Collection tube**. **Centrifuge at 10.000 rpm for 30 seconds**.
10. Place the column in a new **Collection tube**. **Pass the other half and centrifuge at 10,000 rpm for 30 seconds**.
11. Add **100µl of RNA Washing Solution 1**. Centrifuge at maximum speed for 1 minute.
12. Add **700µl of RNA Washing Solution 2**. Centrifuge at maximum speed for 1 minute.

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13. Centrifuge for **3 minutes at maximum speed (11.000 rpm)** to completely dry the membrane. Place the column in a 1.5 ml (not supplied) microtube.
14. Remove the **Collection tube** and insert the **Spin RNA Columns** in a 1.5 ml microtube. Add **25-30µl Nuclease-free water** into reservoir of **Spin RNA Columns**. Dispense buffer directly onto the silica membrane. **Incubate at room temperature for 2 min.**
15. It is very important to add the **Nuclease-free water** in the centre of the membrane to be completely wet.
16. **Centrifuge at 10.000 rpm for 60 seconds.** Collect 25-30µl and redeposit in the centre of the membrane. This increases yield.
17. **Incubate 2 minutes and centrifuge at maximum speed for 1 minute.**

3.7 Protocol for microARN extraction from plant tissues.

Note: The maximum recommended input of plant tissue is 30 mg o 5 x 10⁶ plant cells.

Both fresh and frozen plant samples can be used for this protocol. Samples should be flash frozen in liquid nitrogen and transfer immediately to a -70°C freezer for long-term storage. When isolating MicroRNA from frozen tissues ensure that the temperature does not drop too low during weighing and spraying with a pestle and mortar to ensure that the integrity of the MicroRNA is not compromised.

1. Transfer the tissue into a mortar that contains an appropriate amount of liquid nitrogen to cover the sample. Grind the tissue thoroughly using a pestle. Allow the liquid nitrogen to evaporate, without allowing the tissue to thaw.
 2. Transfer the sample to an RNase-free tube (not provided) and add **400µl de RNA Lysis Solution**. Homogenize by passing the lysate 5-10 times through a 20-G (0.90 mm) gauge needle attached to a syringe or an electrical homogenizer. **Incubate at room temperature for 5 minutes.**
 3. **Centrifuge for 3 minutes a maximum speed.**
 4. Transfer de supernatant to another RNase-free microcentrifuge tube (not provided).
- IMPORTANT:** *Note the volume of the supernatant/lysate.*
5. Add a volume of **100% ethanol** that is equivalent to (55 µl of 100 % ethanol is added to every 100µl of supernatant/lysate). Mix by vortex for 10 seconds.
 6. Assemble a **Columns RNA Large Removal** with one the provided **Spin RNA Columns**.
 7. Apply the clarified lysate with ethanol into the columns. **Centrifuge at 10.000 rpm for 1 minute.**
 8. **Retain the liquid that has passed through the column and that is in the Spin RNA Columns**, which contains the small RNA species.
 9. Add **700µl de ethanol 100% to the lysate**. Mix well.
 10. Pass half of the lysate to a **Spin RNA Columns** with its **Spin RNA Columns**. **Centrifuge at 10.000 rpm for 30 seconds.**
 11. Place the column in a new **Spin RNA Columns**. **Pass the other half and centrifuge at 10,000 rpm for 30 seconds.**
 12. Add **100µl of RNA Washing Solution 1**. Centrifuge at maximum speed for 1 minute.

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13. Add **700µl of RNA Washing Solution 2**. Centrifuge at maximum speed for 1 minute.
14. **Centrifuge for 3 minutes at maximum speed (11.000 rpm)** to completely dry the membrane.
15. Place the column in a 1.5 ml (not supplied) microtube.
16. Remove the **Spin RNA Columns** and insert the **Spin RNA Columns** in a 1.5 ml microtube. Add **25-30µl Nuclease-free water** into reservoir of Spin column. Dispense buffer directly onto the silica membrane. Incubate at room temperature for 2 min.
It is very important to add the elution buffer in the centre of the membrane to be completely wet.
17. **Centrifuge at 10.000 rpm for 60 seconds**. Collect 25-30µl and redeposit in the centre of the membrane. This increases yield.
18. **Incubate 2 minutes and centrifuge at maximum speed** for 1 minute

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STORAGE, STABILITY AND PRECAUTION

RNA Lysis Solution and **Washing Solution 1** contain isothiazacyanate guanidinium which is a potent irritant, wear gloves and goggles. Both buffers can form hazardous reactive compounds when combined with bleach.
All components are stable for 12 months from date of purchase when stored as directed

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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