# REAL SPIN DTR

Ref. RBMS01: (para 50 EXTRACCIONES)

(for 50 EXTRATION)

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## 1.INTRODUCTION

**Real Spin DTRs** are ready-to-use chromatographic microcolumns for rapid purification and clean-up of nucleic acids and proteins using a microcentrifuge.

The gel filtration method is used for purification, in which molecules in solution are separated according to their size differences by passing through a microcolumn packed with a chromatographic medium, in the form of a gel. The gel is a heterogeneous phase in which an aqueous liquid phase is contained within the pores of a solid phase, these pores have a controlled size range depending on the gel matrix used.

Small molecules can diffuse through the gel and are delayed in their passage down the column. Conversely, large molecules cannot diffuse through the gel and leave the column quickly followed by smaller molecules in order of their size.

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- \* Microcentrifuga.

### Equipment and reagents needed and not provided

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### 3.1 Specifications:

- **Applications:** Removal of dideoxy terminators and unincorporated dye labelled nucleotides from sequencing reactions prior to analysis on automated sequencers.
- **Gel matrix:** Special grade polyacrylamide
- **Dyes removed:** 98% retention
- Big Dye Terminators, Rhodamine and d-Rhodamine, Well RED dye
- **Capacity:** 10-75 µl
- **Time:** 4 minutes
- **DNA recovery:** 95% DNA > 22 bp
- **Buffer:** Distilled water

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$$\text{Rpm} = 1000 \times (750 / 1,12 r)^{1/2}$$

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1. Resuspend the **Columns Spin DTR** resin by inversion several times or by gentle vortexing.
2. Remove the bottom and top cap and place in a **Whasing Tubes**.
3. Centrifuge for 2 minutes in a microcentrifuge at 750 xg to remove equilibration buffer.
4. Place the Real Spin DTR in the 1.5 ml sample **Collection tube**
5. Carefully apply the sample (10-75 µl) directly to the centre of the column without damaging the gel.
6. Centrifuge for 2 minutes at 750 xg. The purified sample is collected in the **Collection tube**. **Molecules smaller than the exclusion limit will be retained on the Columns Spin DTR.**

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- All components are stable for 12 months from date of purchase when stored as directed.

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### 1. Poor results:

- 1.1. **High centrifugation speed.** It is very important to set the appropriate centrifugation speed. If a higher speed is applied, the gel matrix will not form correctly.
- 1.2. **Incorrect sample placement.** Avoid touching the gel with the micropipette tip.
- 1.3. **Incorrect placement of the Columns Spin DTR.** During the first centrifugation, the matrix gel is formed with a slanted top. **The Columns Spin DTR** shall be placed in the rotor with this slanted part touching the outside of the rotor.

### 2. Sample dilution

Due to applying a lower speed than indicated, so that not all of the equilibration buffer will be removed and in the second centrifugation step it will be collected with the sample and sample dilution will occur.

To avoid errors when setting the working speed, it is recommended that the first time you work with the **Columns Spin DTR**, a short centrifugation of 30-60 seconds is performed after the centrifugation in step 3 of the protocol. If more buffer is observed, it will be necessary to increase the speed slightly, or increase the centrifugation time in step 3. As recommended in the protocol, the appropriate speed is 750 xg; under these conditions no problems should be encountered.

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