# REALSAFE NUCLEIC ACID STAINING SOLUTION

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## 1. INTRODUCTION

*REALSAFE Nucleic Acid Staining Solution (20,000x) is a new and safe nucleic acid stain, an alternative to the traditional ethidium bromide (EtBr) stain for detecting nucleic acid in agarose gels.*

*It emits green fluorescence when bound to DNA or RNA. This new stain has two fluorescence excitation maxima when bound to nucleic acid, one at 309 nm and another at 419 nm. In addition, it has one visible excitation at 514 nm. The fluorescence emission of REALSAFE bound to DNA is centered at 537 nm.*

*REALSAFE Nucleic Acid Staining Solution (20,000x) is as sensitive as EtBr.*

*The staining protocol for REALSAFE Nucleic Acid Staining Solution (20,000x) is similar to that for EtBr. Compared to EtBr, known as a strong mutagen, REALSAFE Nucleic Acid Staining Solution (20,000x) causes much fewer mutations in the Ames test. In addition, REALSAFE Nucleic Acid Staining Solution (20,000x) has a negative result in mouse marrow chromophilous erythrocyte micronucleus test and mouse spermary spermatocyte chromosomal aberration test.*

## 2. COMPONENTES DEL KIT KIT COMPONENTS

1 ml de REALSAFE Nucleic Acid Staining Solution (20,000x)

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*1 ml of REALSAFE Nucleic Acid Staining Solution (20,000x)*

*Store at 4°C and in the dark.*

### 3. PROTOCOLO GENERAL GENERAL PROCEDURE

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#### 3.1 Features and Applications

- *Used for detecting DNA and RNA.*
- *Alternative to the ethidium bromide staining.*
- *As sensitive as EtBr or more sensitive than that.*
- *NON--toxic, NON-mutagenic and NON-carcinogenic.*
- *NO hazard waste .*
- *Visualization of DNA and RNA bands as they separate during agarose gel electrophoresis*
- *Isolation of DNA fragments for subcloning without introducing mutations normally caused by EtBr*

#### 3.2 Preliminary Preparations

- *Resuspend the REALSAFE Nucleic Acid Staining Solution with vortex mixing briefly before use. It is possible to precipitate that the components of this product is high concentrate.*
- *REALSAFE Nucleic Acid Staining Solution (20,000x) is non-carcinogenic but may cause skin and eye irritations. Please wear gloves when working with the product.*
- *REALSAFE does not create toxic waste.*

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### **3.3 Working protocol**

1. Prepare a 100 ml of agarose gel solution (concentration from 0.8 to 3%) in a 250 ml flask and mix thoroughly. Place the flask in the microwave, heat until the solution is completely clear.

**NOTE:** The thickness of gel should be less than 0.5 cm since thick gels may decrease sensitivity.

2. Cool it by about 60°C. Add 5 µl of REALSAFE Nucleic Acid Staining Solution (20,000x) to the agarose solution. Swirl the flask gently to mix the solution and avoid forming bubbles.

3. While the agarose solution cools, pour into the gel tray until the comb teeth are immersed about 1/4-1/2 into the agarose.

**NOTE:** Repeated melting of gels containing REALSAFE may result in low sensitivity.

4. Allow the agarose gel to cool until solidified. Load the samples on the gel and perform electrophoresis.

5. Detect the bands under UV illumination.

**NOTE:** REALSAFE Nucleic Acid Staining Solution (20,000x) allows the visualization of DNA (>50 ng) in the agarose gel under visible light. This eliminates the need for exposure to UV light, which may nick and damage DNA. The intact DNA fragments purified from agarose gel can increase the efficiency of subsequent molecular biology manipulations such as cloning, transformation and transcription.

**NOTE:** REALSAFE can be used as a post-staining. If the band is not bright, you may re-stain the gel about 10-20 minutes more. You can add 5 µl REALSAFE in 100 ml of buffer solution..

6. REALSAFE can perform with low concentration of DNA. However, the smaller fragments of less than 300bp are not as bright as the bigger one.

7. It is not recommended to use for pulse field electrophoresis.

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We recommend contacting REAL's laboratory technical service in Durviz s.l. for any additional question about work protocols or any problem that may appear while you are working.