REAL TOTAL RNA FROM BACTERIA AND YEAST KIT

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1. INTRODUCTION

*This kit allows the extraction of **total RNA** from different bacterial and yeast cultures avoiding the use of toxic reagents.*

The process includes a cell lysis followed by the precipitation of the proteins and part of the genomic DNA. Later, by a precipitation with isopropanol, the total RNA is obtained, which is finally rehydrated.

Gram (+) strains and yeast are more difficult to lysate than Gram (-) strains and for this reason it's necessary the preincubation with lytic enzymes (not supplied with this kit). In case of Gram (+) strains, lysozyme is used, and zymolase or lyticase is used for yeast.

REAL TOTAL RNA from Bacteria and Yeast Kit allows the extraction of total RNA, although contaminant genomic DNA can also be extracted.

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Production system certified under ISO

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3.1 Preliminary conditions

- RNases exist in procariotes as well as in eucariotes, and almost in every cell type. The main RNases sources are microorganisms (bacteria, fungi and their spores) and their products (restriction enzymes, polymerases, etc.). For such reason, they can be present in a Molecular Biology laboratory in many places and it is convenient to be careful with the following recommendations:*
 - Use gloves to prevent contamination by RNases present in hands.*
 - Change frequently your gloves.*
 - Have pipettes for working exclusively with RNA.*
 - Use RNase-free tubes and tips.*
 - Use RNase-free reagents.*
 - Have a determined area in the lab to work with RNA.*

3.2 Preliminary Preparations

- If the lysis solution contains a precipitate due to the low temperatures, incubate at 37°C and mix to dissolve the precipitate*

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3.3 Protocol for total RNA extraction from 1 ml of Gram (+) or (-) bacterial cultures

Gram (+) strains are more difficult to lyse than Gram (-) strains, for this reason it is recommended to do incubation with lytic enzymes. For certain species of Staphylococcus it is necessary a lysis with a mixture of lysozyme (10 mg/ml) and lysostaphine (10 mg/ml)

1. Add 1ml of an overnight culture to a 1.5 ml. microtube.
2. Centrifuge at 14.000 xg for 30 seconds. Remove the supernatant. **For Gram (+) strains proceed with step 3. For Gram (-) strains go directly to step 6.**

3. Resuspend the cells in 540 µl of nuclease-free water or EDTA 50 mM.
4. Add 60 µl of Lysozyme (10 mg/ml). The purpose of this treatment is to weak the cell wall to make the cell lysis more efficient.
5. Incubate the sample at 37°C for 60 minutes. Invert the periodically the sample during the incubation. Centrifuge at 14.000 xg for 2 minutes. Remove the supernatant.
6. Add 600 µl of **Lysis Solution (ER01B)** to the cell pellet. Pipette to resuspend and to lyse the cells.
7. Incubate the sample at 65°C for 5 minutes. Allow to air dry.

8. Add 300 μ l of **Protein Precipitation Solution. (ER02B)** Invert the tube about 8-10 times and incubate in ice for 5 minutes.
9. Centrifuge at 14.000 xg for 5 minutes. The protein precipitate will form a tight pellet.
10. Pour the supernatant containing the RNA in a new microtube with 600 μ l of **Isopropanol**. Mix by inversion for several times.
11. Centrifuge at 14.000 xg for 3 minutes.
12. Remove the supernatant. Add 600 μ l of **Ethanol 70 %** and invert the tube several times to wash the RNA pellet.
13. Centrifuge at 14.000 xg for 2 minutes. Carefully, avoiding losing the pellet. Remove all the ethanol.
14. Invert and drain the tube on absorbent paper and allow to air dry for 15 minutes.
15. Add 100 μ l of **Rehydration Solution. (ER03)**. Use the pipettor for resuspending the pellet.
16. Allow the RNA to rehydrate on ice for at

least 30 minutes. Use the pipettor for resuspending the pellet.

17. Store the samples at -70°C .

- ✓ If the obtained RNA is going to be used in RT-PCR assays the use of the REAL Total RNA from Bacteria and Yeast STAR kit (ref. RBMER04) or the REALSTAR kit (Ref. RBMER10) it is recommended to remove the contaminant DNA.
- ✓ False positives due to contaminant genomic DNA in a RT-PCR can be identified with a "no reverse transcription" control during the RT-PCR assay. The problem can also be minimized with the design of primers which are located at least inside an intron of the genomic sequence. This way, the PCR product obtained from contaminant DNA will be bigger than the one generated from the cDNA.

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3.4 Protocol for total RNA extraction from 1 ml of yeast culture

1. Add 1ml of an overnight culture to a 1.5 ml. microtube.
 2. Centrifuge at 14.000 xg for 2 minutes. Remove the supernatant.
 3. Resuspend the cells in 592µl of nuclease-free water or EDTA 50 mM.
 4. Add 8 µl of Lyticase or Zymolase (20 mg/ml). The purpose of this treatment is to weak the cell wall to make the lysis more efficient.
 5. Incubate the samples at 37°C for 60 minutes. Invert periodically the sample during the incubation. Centrifuge at 14.000 xg for 2 minutes. Remove the supernatant.
 6. Add 600 µl of **Lysis Solution (ER01B)** to the cell pellet. Use a pipettor to resuspend and to lyse the cells.
 7. Incubate the sample at 65°C for 5 minutes. Allow to air dry.
 8. Add 300 µl of **Protein Precipitation Solution (ER02B)**. Invert the tube for 8-10 times and incubate in ice for 5 minutes.
 9. Centrifuge at 14.000 xg for 5 minutes. The protein precipitate will form a tight pellet.
 10. Pour the supernatant containing the RNA in a new microtube containing 600 µl of **Isopropanol**. Mix by inversion several times.
 11. Centrifuge at 14.000 xg during 3 minutes.
 12. Remove the supernatant. Add 600 µl of **Ethanol 70 %** and invert the tube several times to wash the RNA pellet.
 13. Centrifuge at 14.000 xg for 2 minutes. Carefully, avoiding losing the pellet, removing all the ethanol.
 14. Invert and drain the tube on absorbent paper to allow to air dry for 15 minutes.
 15. Add 100 µl of **Rehydration Solution (ER03)**. Use a pipettor for resuspending the pellet.
 16. Allow the RNA to rehydrated on ice for at least 30 minutes. Use a pipettor for resuspending the pellet.
 17. Store the samples at -70°C.
- ✓ If the obtained RNA is going to be used in RT-PCR the use of the REAL Total RNA from Bacteria and Yeast STAR kit (ref. RBMER04) or the REALSTAR kit (Ref. RBMER10) it is recommended to remove the contaminant DNA.
 - ✓ False positives due to contaminant DNA in a RT-PCR can be identified with a “no reverse transcription” control during the RT-PCR assay. The problem can also be minimized with the design of primers which are located at least inside an intron of the genomic sequence. This way the PCR product obtained from contaminant DNA will be higher than the generated from the cDNA

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Due to the wide number of samples that can be treated to extract total RNA using this kit, it's difficult to generalize the possible problems. Please contact REAL Technical Service (in Durviz, s.l.) for any comment or question regarding the protocol.