REAL SPIN GENOMIC DNA KIT

Ref. **RBMEGS01:** (para **50 EXTRACCIONES**)
 Ref. **RBMEGS02:** (para **250 EXTRACCIONES**)
 Ref. **RBMEGS15:** (para **1000 EXTRACCIONES**)

(for **50 EXTRACTION**)
 (for **250 EXTRACTION**).
 (for **1000 EXTRACTION**).

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1. SAMPLE TAKING

This kit is designed for rapid extraction and purification of high quality genomic DNA from a wide variety of samples including whole blood, cultured cells, animal tissues, mouse tails, bacteria, yeast, clinical samples (serum, plasma, stool, urine), forensic samples and paraffin-embedded tissues.

The kit contains sufficient reagents and columns to perform genomic DNA extractions according to the kit from different samples:

- 200 µl whole blood.
- 200 µl buffy coat.
- 104-106 cells in culture.
- 25-50 mg tissue.
- 0.2-0.5 cm mouse tail.
- 108 bacteria.
- 109 yeast.
- Paraffin-embedded tissue sections.
- **For any other type of sample, please request the protocol to DURVIZ S.L. technical service.**

The procedure includes a lysis in SDS and proteinase K, then a chaotropic solution is added to create the necessary conditions for the DNA to bind to the glass fiber membrane and finally the DNA is eluted with an elution buffer.

The DNA obtained is of high quality and can be used directly in PCR, Southern blotting, cloning, any enzymatic reaction, etc...

2. COMPONENTES DEL KIT / KIT COMPONENTS

Componente REAL	Ref.	Envase. (50 Preps)	Envase. (250 Preps)	Envase. (1000 Preps)	T ^a
Lysis Solution	E21	10 ml	50 ml	200 ml	RT
Lysis/Binding Solution	REA01	10 ml	50 ml	200 ml	RT
Proteinase K	REA02	22 mg	110 mg	440 mg	-20 °C
Disinhibition Solution (GHCL 8M)	REA03	16.5 ml	82.5 ml	330 ml	RT
Washing Solution	EP08	10 ml	50 ml	200 ml	RT
Hydration Solution	REA05	10 ml	50 ml	200 ml	RT
Columnas Spin DNA	RSC03	50 unid	250 unid	1000 unid	RT
Collection Tubes	R30	50 unid	250 unid	1000 unid	RT

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Equipment and reagents needed and not provided

- * Isopropanol.
- * Nuclease-free water.
- * 50 mM EDTA and 7.5 µl of Lyticase (75 units/ µl) SIGMA ref. L2524
- * PBS with 20 mg/ml lysozyme and/or lysostaphin
- * 100% ethanol.
- * Microtubes of 1.5 ml.
- * Microcentrifuge or clinical centrifuge.

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The protocol involves the following steps:

- Samples are lysed with the appropriate **Lysis Solution** and **Proteinase K**.
- Nucleic acids are bound to the glass fiber matrix packed in the **Columnas Spin DNA**.
- The nucleic acids are first washed with **Disinhibition Solution (GHCL8M)** to remove PCR inhibitors.
- Nucleic acids are washed to remove salts, proteins and other impurities.
- Nucleic acids are eluted.

3.1 Preliminary preparations

- Dissolve **Proteinase K** in **1.1 ml** of nuclease-free water (kit 50 extractions) and in **5.2 ml** (kit 250 extractions) and store at -20°C . It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- Verify that the **Lysis Solution, Lysis/Binding Solution** do not have precipitates due to low temperatures. If necessary, dissolve by heating to 37°C .
- Add 100% Ethanol to the **Disinhibition Solution (GHCL8M)** indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- Add 100% Ethanol to the **Washing Solution** indicated on the label, about **40 ml** (kit 50 extractions) and about **200 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
- **Pre-warming the Hydration Solution to 70°C may increase the yield of DNA obtained. For some subsequent application it may be necessary for the DNA to be concentrated, elution in volumes smaller than 200 μl will increase the final concentration of DNA in the eluate, but will reduce the overall yield of DNA obtained. For samples containing $< 3 \mu\text{g}$ DNA, elution in 100 μl is recommended. For samples containing $< 1 \mu\text{g}$ DNA, a 50 μl elution is recommended.**

3.2 Protocol for extraction of genomic DNA from whole blood, buffy coat and cells in culture

- 200 μl whole blood.
- 200 μl buffy coat.
- 104- 106 cells in culture

If the material does not reach 200 μl , bring the final volume of the sample to 200 μl with nuclease-free water.

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1. To 200 µl of the indicated material add **200 µl of Lysis/Binding Solution + 20 µl Proteinase K**. Mix well. Incubate at 70°C for 10 minutes.
2. Add **100 µl of Isopropanol**. Mix well.
3. Pipette the lysate into the **Spin Columna DNA** reservoir with its **Collection tube**. Centrifuge at **8,000 rpm for 60 seconds**. Discard the **Collection tube**.
4. Place the **Spin Columna DNA** in a new **Collection tube** and add **500 µl of Disinhibition Solution (GHCL8M)** to the reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Discard the liquid.
5. Add **500 µl of Washing Solution** to the **Spin Columna DNA** reservoir. Centrifuge at 12,000 rpm for 60 seconds. Remove the liquid.
6. 2º Wash. Add 500 µl of **Washing Solution** to the **Spin Columna DNA** reservoir. Centrifuge at **14,000 rpm for 60 seconds**. Eliminate the liquid.
7. Centrifuge at maximum speed for **90 seconds to remove residual ethanol**.
8. Remove the **Collection tube** and insert the **Spin Columna DNA** into a 1.5 ml microtube. Add **50-200 µl of Hydration Solution** (preheated to 70°C) to the **Spin Columna DNA** reservoir. Incubate for 2 minutes.
9. Centrifuge at maximum speed for **60 seconds**. The microtube now contains the genomic DNA.

3.3 Protocol for the extraction of genomic DNA from 25 mg of animal tissues

1. Cut 25 mg of human or animal tissue into small pieces and place into a 1.5 ml microtube. Add **180 µl of Lysis Solution + 20 µl Proteinase K**. Mix well. Incubate at 55°C for 1 hour or until lysis is complete, samples can be incubated "overnight".
Samples that are difficult to lyse can be ground with liquid Ni or can be treated directly with a Polytron type homogenizer.
2. Add **200 µl of Lysis/Binding Solution**. Vortex. Incubate at 70°C for 10 minutes. If insoluble particles are found, centrifuge 5 minutes at maximum speed and transfer the supernatant to a new microtube.
3. Add **100 µl of Isopropanol**. Mix well.
4. Pipette the lysate into the **Spin Columna DNA** reservoir with its **Collection Tube**. Centrifuge at **8,000 rpm for 60 seconds**. Discard the **Collection Tube**.
5. Place the **Spin Columna DNA** in a new **Collection Tube** and add **500 µl of Disinhibition Solution (GHCL8M)** to the reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Discard the liquid.
6. Add **500 µl of Washing Solution** to the **Spin Columna DNA** reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Eliminate the liquid.
7. 2º Wash. Add 500 µl of **Washing Solution** to the **Spin Columna DNA** reservoir. Centrifuge at **14,000 rpm for 60 seconds**. Eliminate the liquid.
8. Centrifuge at maximum speed for **90 seconds to remove residual ethanol**.
9. Remove the **Collection Tube** and insert the **Spin Columna DNA** into a 1.5 ml microtube. Add **50-200 µl of Hydration Solution** (preheated to 70°C) to the spin microcolumn reservoir. Incubate for 2 minutes.

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10. **Centrifuge at maximum speed for 60 seconds.** The microtube now contains the genomic DNA.

3.4 Protocol for extraction of genomic DNA from 25-50 mg mouse tails

1. Cut 0.2-0.5 cm of mouse tail into several pieces and place in a 1.5 ml microtube. Add **180 µl of Lysis Solution + 20 µl Proteinase K**. Mix well. Incubate at 55°C until lysis is complete, samples can be incubated "overnight". To remove bone or hair residues, centrifuge for 5 minutes at maximum speed and transfer the supernatant to a new microtube.
2. Add **200 µl of Lysis/Binding Solution and 100 µl of Isopropanol**. Mix well by vortexing.
3. Pipette the lysate into the **Spin columna DNA** reservoir with its **Collection Tube**. Centrifuge at **8000 rpm for 60 seconds**. Discard the **Collection Tube**.
4. Place the **Spin columna DNA** in a new **Collection Tube** and add **500 µl of Disinhibition Solution (GHCL8M)** to the reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Discard the liquid.
5. Add **500 µl of Washing solution** to the **Spin columna DNA** reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Remove the liquid.
6. 2nd Wash. Add **500 µl of Washing solution** to the **Spin columna DNA** reservoir. Centrifuge at **14,000 rpm for 60 seconds**. Eliminate the liquid.
7. Centrifuge at maximum speed for **90 seconds to remove residual ethanol**.
8. Remove the **Collection tube** and insert the **Spin columna DNA** into a 1.5 ml microtube. Add **50-200µl of Hydration Solution** (preheated to 70°C) to the **Spin columna DNA** reservoir. Incubate for 2 minutes.
9. Centrifuge at maximum speed for **60 seconds**. The microtube now contains the genomic DNA.

3.5 Protocol for the extraction of genomic DNA from 10⁹ bacteria

1. Centrifuge 1-1.5 ml of bacterial culture. Discard the supernatant. Resuspend the pellet in **180 µl of Lysis Solution and then add 20 µl Proteinase K**. Vortex and incubate at 55°C until lysis is complete. For those strains that are difficult to lyse, especially Gram+ strains, prior incubation with lytic enzymes is necessary. Resuspend the pellet with 200 µl of PBS with 20 mg/ml lysozyme and/or lysostaphin (10 mg/ml), incubate 30 minutes at 37°C, then add **20 µl Proteinase K** and incubate at 55°C until lysis is complete.
2. Add **200 µl of Lysis/Binding Solution**. Vortex. Incubate at 70°C for 10 minutes.
3. Add **100 µl of Isopropanol**. Mix well.
4. Pipette the lysate into the **Spin columna DNA** reservoir with its **Collection Tube**. Centrifuge at **8,000 rpm for 60 seconds**. Discard the **Collection Tube**.
5. Place the **Spin columna DNA** in a new **Collection Tube** and add **500 µl of Disinhibition Solution (GHCL8M)** to the reservoir. Centrifuge at **12,000 rpm for 60 seconds**. Discard the liquid.

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6. **Add 500 µl of Washing Solution** to the **Spin columnna DNA** reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Eliminate the liquid.
7. 2º Wash. **Add 500 µl of Washing Solution** to the **Spin columnna DNA** reservoir. **Centrifuge at 14,000 rpm for 60 seconds**. Eliminate the liquid.
8. **Centrifuge at maximum speed for 90 seconds to remove residual ethanol**.
9. Remove the **Collection Tube** and insert the **Spin columnna DNA** into a 1.5 ml microtube. **Add 50-200 µl of Hydration Solution** (preheated to 70°C) to the **Spin columnna DNA** reservoir. Incubate for 2 minutes.
10. **Centrifuge at maximum speed for 60 seconds**. The microtube now contains the genomic DNA.

3.6 Protocol for extraction of genomic DNA from 10⁸ yeast cultures

1. Centrifuge 13000-16000 x g for 2 minutes 1.5-3 ml of yeast culture. Remove the supernatant. **Resuspend the pellet in 293 µl of 50 mM EDTA and 7.5 µl of Lyticase** (75 units/ µl) SIGMA ref. L2524. **Incubate at 37°C for 30-60 minutes**. Centrifuge at maximum speed and remove the supernatant. Resuspend the pellet with **180 µl of Lysis Solution and then add 20 µl Proteinase K**. Vortex. Incubate at 55°C until lysis is complete.
2. **Add 200 µl of Lysis/Binding Solution**. Vortex. Incubate at 70°C for 10 minutes.
3. **Add 100 µl of Isopropanol**. Mix well.
4. Pipette the lysate into the **Spin columnna DNA** reservoir with its **Collection Tube**. **Centrifuge at 8000 rpm for 60 seconds**. Discard the **Collection Tube**.
5. Place the **Spin columnna DNA** in a new **Collection Tube** and add **500 µl of Disinhibition Solution (GHCL8M)** to the reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Discard the liquid.
6. **Add 500 µl of Washing Solution** to the **Spin columnna DNA** reservoir. **Centrifuge at 12,000 rpm for 60 seconds**. Remove the liquid.
7. 2º Wash. **Add 500 µl of Washing Solution** to the **Spin columnna DNA** reservoir. **Centrifuge at 14,000 rpm for 60 seconds**. Eliminate the liquid.
8. **Centrifuge at maximum speed for 90 seconds to remove residual ethanol**.
9. Remove the **Collection Tube** and insert the **Spin columnna DNA** into a 1.5 ml microtube. **Add 50-200µl of Hydration Solution** (preheated to 70°C) to the **Spin columnna DNA** reservoir. Incubate for 2 minutes.
10. **Centrifuge at maximum speed for 60 seconds**. The microtube now contains the genomic DNA.

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- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: *Lysis/Binding Solution* and *Disinhibition Solution (GHCL8M)* contain guanidine hydrochloride which is irritant, wear gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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