REAL SPIN FOOD-STOOL BACTERIA KIT

Ref. **RBMEGS03:** (para **50 EXTRACCIONES**)
 Ref. **RBMEGS04:** (para **250 EXTRACCIONES**)

(for **50 EXTRACTION**)
 (for **250 EXTRACTION**).

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SAMPLE TAKING

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This kit is optimized for efficient and rapid extraction of PCR-ready bacterial DNA (Listeria, Salmonella, E.coli, etc.) from pre-enrichment (peptone water) or enrichment media as stated in the different ISO standards, from different food and stool samples, using MicroSpin columns with glass fiber membranes that selectively bind DNA.

The procedure includes a centrifugation of 1 ml of pre-enrichment (peptonized water) or enrichment medium to concentrate the cells, the cells are lysed during an incubation period with Lysozyme (supplied with the kit). This is followed by digestion with Proteinase K (supplied with the kit) and cleaning of the lysis mixture by centrifugation and binding of the DNA to glass fiber membranes packed in MicroSpin columns. The bound DNA is purified by several washes to remove potential PCR inhibitors and finally eluted. This kit is designed for rapid extraction and purification of high quality genomic DNA from a wide variety of samples including whole blood, cultured cells, animal tissues, mouse tails, bacteria, yeast, clinical samples (serum, plasma, stool, urine), forensic samples and paraffin-embedded tissues.

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Equipment and reagents needed and not provided

- * Isopropanol.
- * 100% ethanol.
- * Microtubes of 1.5 ml.
- * Microcentrifuge or clinical centrifuge

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3. GENERAL PROCEDURE

11. **Centrifugar a máxima velocidad durante 60 segundos.** El microtubo contiene ahora el ADN bacteriano.
 - Dissolve Proteinase K in **1.10 ml** of nuclease-free water (kit 50 extractions) and in **5.25 ml** (kit 250 extractions) and store at -20°C . It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
 - Dissolve Lysozyme in **1.10 ml** of nuclease-free water (kit 50 extractions) and in **5.25 ml** (kit 250 extractions) and store at -20°C . It is recommended to make several aliquots to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
 - Verify that the Lysis/Binding Buffer does not have precipitates due to the low temperatures. If necessary, dissolve and heat to 37°C .
 - Add 100 % Ethanol to the Disinhibition Buffer indicated on the label, about **10 ml** (kit 50 extractions) and about **50 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
 - Add 100% Ethanol to the Wash Buffer indicated on the label, about **40 ml** (kit 50 extractions) and about **200 ml** (kit 250 extractions). Keep the container tightly closed to avoid evaporation of the ethanol.
 - Pre-heat the Elution Buffer to 70°C .

2.2 Protocol for extraction of bacterial DNA from 1-1.5 ml of pre-enrichment or enrichment media

1. Transfer 1-1.5 ml of pre-enrichment or enrichment medium to a 1.5 ml microtube and **centrifuge for 5 minutes at 12,000-14,000 rpm. Discard the supernatant.**
2. Add **280 μl of Lysozyme Reaction Buffer + 20 μl Lysozyme.** Mix well with micropipette to resuspend the bacterial pellet. Incubate at 37°C for 30 minutes.
3. Add **300 μl Lysis/ Binding Buffer + 20 μl Proteinase K.** Mix well. Incubate at 70°C for 10 minutes.
4. Add **150 μl Isopropanol.** Mix well and centrifuge 60 seconds at 14,000 rpm.
5. Add the supernatant to the Spin microcolumn reservoir with the collection tube. **Centrifuge at 8,000 rpm for 60 seconds.** Discard the collection tube.
6. Place the MicroSpin column in a new collection tube and add **500 μl of Deinhibition Buffer to the reservoir. Centrifuge at 8,000 rpm for 60 seconds.** Discard the liquid.
7. Add **500 μl of Wash Buffer** to the Spin microcolumn reservoir. **Centrifuge at 12,000 rpm for 60 seconds.** Remove the liquid.
8. 2º Wash. Add **500 μl of Wash Buffer** to the Spin microcolumn reservoir. **Centrifuge at 14,000 rpm for 60 seconds.** Remove the liquid.
9. **Centrifuge at maximum speed for 2 minutes to remove residual ethanol.**
10. Remove the collection tube and insert the spin microcolumn into a 1.5 ml microtube. **Add 100 μl of Elution Buffer** (pre-warmed to 70°C) to the spin microcolumn reservoir. Incubate for 1 minute.

11. **Centrifuge at maximum speed for 60 seconds.** The microtube now contains the bacterial DNA. 104- 106 cells in culture

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STORAGE AND STABILITY

- All components are stable for 12 months from date of purchase when stored as directed.

CAUTION: Lysis/Union Buffer and Disinhibition Buffer contain guanidine hydrochloride which is irritating, use gloves and goggles.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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