REAL MICROBIAL DNA KIT

Ref. **RBMEGS17:** (para **50 EXTRACCIONES**)

(for **50 EXTRATION**)

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1.INTRODUCTION

REAL Microbial DNA kit is designed for rapid purification of high quality genomic DNA from microorganisms (gram negative bacteria, large positive bacteria, yeasts and fungi).

Microbial samples such as gram positive bacteria, yeasts and spores can be difficult to lyse due to the complex structure of the cell wall. The REAL Microbial DNA kit replaces enzymatic lysis by using a mechanical disruption system that weakens cell wall structures based on beads. In addition, the use of beads can be combined with other cell disruption devices.

The protocol starts with a bead-beating process, the cells are lysed by a combination of mechanical forces, heat and detergents, a vortex with a horizontal adapter for the Genie 2 Vortex. Other common disruption devices can also be used.

The appropriate conditions for DNA to bind to the Microbial DNA columns are achieved by the addition of large amounts of chaotropic salts (**Binding Solution**) to the lysate. Contaminants are removed with two efficient washing steps. The resulting DNA is then immersed in a DNA-free tris buffer in which the following reactions will take place.

Features:

- Designed for **rapid purification of high-quality genomic DNA from microorganisms** (gram-negative bacteria, gram-positive bacteria, yeasts and fungi).
- Suitable for different types of samples. **Cell cultures and agar plates.**
- **Efficient lysing** due to the use of *beads* in combination with liquid *proteinase k*

Applications:

- a) Total DNA for cell culture.
- b) Typical applications: PCR, Real-Time PCR, southern blotting, enzymatic reactions.

2.COMONENTES DEL KIT / KIT COMPONENTS

Componentes REAL	Ref.	Envase	T ^a
Washing Solution*	EP08	10 ml	RT
Bead Microtubes	RBD5	50 unid	RT
Proteinase K*	REA02	30 mg	-20°C
Disinhibition Solution (GHCL 8M)*	REA03	18 ml	RT
Hydration Solution	REA05	10 ml	RT
CTAB Extraction Solution	RES01	45 ml	RT
Binding Solution	RES04	15 ml	RT
Columns Spin DNA	RSC03	50 unid	RT
Collection Tubes	R30	100 unid	RT

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* See the section on preliminary preparations on how to prepare these solutions.

Equipment and reagents needed and not provided

- * Ethanol 100%
- * Microcentrifuge
- * 1.5 ml and 2.0 ml microtubes
- * Water bath, dry bath or heated block (70°C).
- * Vortex for homogenization of the bead microtubes, we recommend the use of the Vortex Genie 2 or similar with a horizontal rack.
- * Bead mill homogenizer (optional).

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3. GENERAL PROCEDURE

3.1 Preliminary preparations

- Dissolve **Proteinase K** in 1.3 ml in nuclease-free water and store at -20°C. **It is recommended to make several aliquots** to avoid too many thaw-freeze cycles. At this temperature it is stable for 1 year.
- **Add 10 ml of 100% Ethanol to the Disinhibition Solution (GHCL 8M)** indicated on the label. Keep the container tightly closed to prevent evaporation of the ethanol.
- **Add 40 ml of 100 % Ethanol to the Washing Solution** indicated on the label. Keep the container tightly closed to avoid evaporation of the ethanol.
- Preheat the **Hydration Solution** to 70°C

3.2 General Considerations

- Cells should be collected from fresh microbial cultures by sedimentation via centrifugation. The supernatant should be removed by aspiration. Microbial cell pellets can be used fresh or stored at -20°C or -80°C before starting DNA extraction.

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- The procedure is optimized by using **beads** in combination with a horizontally shaking vortex (Vortex Genie 2 or similar). Horizontal agitation is essential, so an adapter that allows this orientation is required.
- Adapters with a vertical orientation will not allow proper agitation.
- It is possible to use a **beads** mill for homogenization such as Fast Prep, Precellys and others, but always follow the manufacturer's instructions to optimize sample lysis.

IMPORTANT: many modern disruption devices can apply too much energy. Depending on the tube type and content (**beads**, liquid volume, sample type), especially with high shaking frequencies and/or too long shaking times can cause tube breakage. It is the responsibility of the user to initially test the tubes to be used to ensure their stability and resistance to shaking conditions. This initial test can be performed with water instead of **Lysis Solution** and with low equipment parameters (low frequency, short time) in order to avoid reagent spillage in case of breakage.

- Apart from the standard method, it is possible to modify the protocol to extend the amount, concentration and convenience.

Standard elution: For convenience, the elution can take place in one step by adding 100 µL of elution buffer to the column.

High quantity: two serial elutions of 100 µL each for an eluate volume equal to 200 µL.

High concentration: use the 100 µL of the initial eluate for a second elution. 100 µL total eluate volume, two elutions.

3.3 Protocol for isolation of microbial DNA from cell culture or agar plates

- a) Cell cultures: Collect culture cells by centrifugation in a microcentrifuge tube (not provided). Discard the supernatant. Add 800 µl of **CTAB Extraction Solution**. Resuspend the pellet using a micropipette. No Vortexing. Transfer to the 2ml microtube containing the particles or **beads**. Up to approximately 50 mg of wet cell culture sample can be used.
 - b) Agar plates: Using a seeding loop transfer approximately 50 mg of colonies directly into the microtube with **beads** to which **800 µl of CTAB Extraction Solution** has been added. Be careful not to pick up any agar. Resuspend the pellet with the aid of a micropipette. Do not vortex.
1. **Add 25 µl of Proteinase K. Incubate at 70°C for 10 minutes.**
 2. **Homogenize** by **Beads** for 10 minutes at maximum speed with the Vortex Genie 2 or similar, always with a **horizontal adapter**.
 3. **Centrifuge at 14,000 rpm for 5 minutes.**
 4. Transfer 500 µL of the supernatant to a clean microcentrifuge tube.
IMPORTANT: A layer of debris may appear on the surface of the pellet. Avoid dragging this debris with the supernatant.
 5. Add **250 µl of Binding Solution** and vortex briefly.

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For any further questions or doubts about the protocol, please contact DURVIZ s.l. technical service at durviz@durviz.com.

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